

Product datasheet for **SC327922**

ZNF529 (NM_020951) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF529 (NM_020951) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF529
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC327922 representing NM_020951.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGCCAACTCAAGTTTCATTGGGGATCATGCCATGGAGCTCCTCATGCTGTCATCGCCAGAAGTGGAA
TTCCCTGACCAATTCTTACAGTTCTAACCATGGACCATCGAACTGGTGACACTCAGGGATGTGGTCATC
AACTTTCTCTCAGGAGGAATGGGAATATCTGGATTCGCTCAGAGGAACTGTACTGGGATGTGATGATG
GAGAACTACAGCAACTTACTCTCACTGGATTTGGAGTCCAGGAATGAGACTAAGCATTTATCTGTAGGA
AAGATATTTTCAAACACTGGTTCTCAGTGGGAGTAATGGAAAGTAGCAAGTTATGTGGCCTTGAA
GGTTCCATTTTCAGAATGACTGGCAAAGCAAAGCAAGATTGACTTACAAGGACCTGAAGTGGGATAT
TTCAGTCAAATGAAAATCATCTCTGAAAATGTGCCAGTTACAAAACTCATGAATCTCTTACTTTACCT
CGGAGAACTCATGACAGTGAGAGCCTATGAATACAAGGAATATGAGAAGTCTCAGTGTGACTTA
GAGTTTGATGAATATCAGAAATACACTGGTGGAAAAAACTATGAATGTAATCAATGTTGGAAAACC
TTTGGGATAGATAACTCCAGTATGTTACAACTGAATATCACTGGTGTGAAACCTTGTAATATATATG
GAATATGGGAATACATGTAGTTTTTATAAGACTTTAATGTATACCAGAAATTCATAATGAGAAGTTCT
TATAAATGTAAGGAATACAGAAGACCTTTGAAAGAGTTGGAAAAGTTACTCCACTTTCAAGAGTTCAT
GATGGTGAGAACACTTTGAATGCTATTCTGGGAAATCCTTTAGAGTGCATGCACAACTTACTCGA
CATCAGAAATCCACTGATGAGAAAACTTACAAATGTATGGAATGTGGCAAGGACTTCAGATTTCAT
TCAGCTTACCGAACATCAGAGAATCTACTGGTGAGAACCCTACAAATGTATGCACTGTGAGAG
GTTTTTAGAATAGTTCACAGCTCATTGACATCAGAGAATCACTGGTGAGAACCTTATGCATGT
AAGGAATGTGGGAAGCTTTGGAGTATGAGAAACTTGCTCGCATCAGAGAATCACTGGTAAG
AAACCCTATGAATGCAAAAGCATGTGGAAAGCTTTTAGAAAATAGTTCATCCCTGACTAGCATCAGAG
ATTCACTGGTGAAAAACCCTATAAATGTAAAGAAATGAGAAAGCATTTGGAGTAGGTAGGAACTT
ACTCGACATGAAAGAATCACAGTGGTCAAAAACCTTATGAATGTAAGAGTGTGGAAAGTTCTTTAG
CTTACGCAGCCTTATTTCAACATCAAAGAATTCATAGTGGTGAGAACCTTATGAATGTAAGATGT
GGGAAGGCCTTAGCATAGTTCAGCCTTACAGAACATCAGAGAATCTACTGGGAAAAAACCCTAT
GAATGCAAGGCATGTGGGAAGGCCTTAGCATAGTTCATCCTTACCAAACATCAGCGCATCTACT
GATGATAAACCCTATGAATGTAAGAGTGTGGGAATCCTTAGTGTTGGGCATCTACTGGCA
CCGAAATTTACTCTGGTGAGAAATCATTTGACTGA
ACGCGTACGCGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGCGCCGC
  
```

Restriction Sites:

SgfI-MluI

Plasmid Map:



ACCN:	NM_020951
Insert Size:	1692 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_020951.4</u>
RefSeq Size:	5160 bp
RefSeq ORF:	1692 bp
Locus ID:	57711
UniProt ID:	<u>Q6P280</u>
Cytogenetics:	19q13.12
Domains:	zf-C2H2
Protein Families:	Transcription Factors
MW:	65.9 kDa
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (2) differs in the 5' UTR compared to variant 1. Both variants 1 and 2 encode the same isoform (a). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.