

## Product datasheet for **SC327931**

### EGR4 (NM\_001965) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	EGR4 (NM_001965) Human Untagged Clone
Tag:	Tag Free
Symbol:	EGR4
Synonyms:	AT133; NGFI-C; NGFIC; PAT133
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene sequence for NM_001965 edited

```

ACCACCATGCTCCACCTTAGCGAGTTTTCCGAACCCGACGCGCTCCTCGTCAAGTCCACT
GAAGGCTGTTGCGCCGAACCCAGCGCTGAATTGCCCGGCTGCCTGCCAGGGACGTCCTCC
GCGGCCACCGGCTACCCTGGAGCAGGCGACTTCTTGAGCTGGGCTTTGAACAGCTGCGGC
GCAAGTGGGGAAGTACCGGACTCCTGCTTCTGGAGGGGCTGCGCCACACCCCTCC
GGCCTCAGCTACAGCGGTAGCTTCTTATTAGGCAAGTGCCTGAACACCCGACGACCCG
GAGGCACTCTTCAACCTCATGTGCGGCATCTTAGGCCTGGCACCTTCCCGGTCCAGAG
GCAGCAGCGTCCAGATCCCGCTGGATGCCCCCTTTCTGCGGGTCCGATGCCTTGCTG
CCGGGTCGCGGACCTTACTCCCGGATCTGGGCGCTGCCCTTTCCAGAGGCGTTC
TGGGAGGCTCGCTTGGCGGGTGCCCCCTCGCAGTGCCTGTATGAGCCTCAGCTCTCC
CCGCCGACGTCAAGCCCGGCTCCGGGCGCTCCCGCTCGCCAGCGCTGGACGCTGTC
TCTGCCTTCAAGGTCCCTACGCGCCCTGGGAGCTGCTTCTGTGGGGGCCCCAGGGAAC
TGTGGGTACAGGGAGACTACCAGGCCGCCCGGAGGCTCGTTTTCCCGTAATAGGGACC
AAGATTGAGGACTTGCTGTCCATCAGCTGCCCTGCGGAAGTCCCGGCCGTCCTCCAGCCAAC
AGACTCTATCCAGCGGGGCTATGACGCTTCCCGCTGGCCCCGGTGACTTAGGGGAG
GGGGCTGAGGGCTCCCTGGGCTCCTGACCCCTCCTAGTGGGGAGGGAGGAGTAGCGGC
GACGGCGGAGAGTTTCTGGCCAGTACGACGCTCAGCTTCCCGCTGGGCTTTCGACG
GCCGCCGCGGCGACTTCCCTAAACCTCTGGTGGCGGACATCCCTGGAAGCAGTGGCGTG
GCTGCACACCCGTGCCGCGCGCCGCCACCCCTTTCCCGAGGCCAAGGCGCGACGC
AAGGGGCGCGCGCGGCAAAATGACGACGCGCTGCTTCTGCCCGCGGCGCACGCAAG
GCCTTCGCTTGCCCGGTGAGAGTTGTGTGCGGAGCTTTGCGCGCTCCGACGAGTCAAT
CGCCACCTGCGCATCCACACGGGCCACAAACCTTCCAGTGCCGCATCTGCCTCCGCAAC
TTCAGCCGACGACACCTCACCACGACGTGCGCACCCACACCGGCGAGAAGCCTTTT
GCTTGCGACGTGTGCGGCCGCGCTTCCGCGCAGCGATGAGAAGAAACGGCACAGCAAG
GTGCACCTCAAGCAGAAGGCGCGCGCGAGGAGCGGCTCAAGGGCTCGGCTTTTACTCG
CTGGGCTCTCCTTCGCTTCTCTGAGCAAGAGATGGGTTTATGGGTTAAGGGC
  
```



[View online »](#)

<b>Restriction Sites:</b>	Please inquire
<b>ACCN:</b>	NM_001965
<b>OTI Disclaimer:</b>	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at <a href="mailto:custsupport@origene.com">custsupport@origene.com</a> or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a></p>
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<a href="#">NM_001965.3</a> , <a href="#">NP_001956.3</a>
<b>RefSeq Size:</b>	2379 bp
<b>RefSeq ORF:</b>	1770 bp
<b>Locus ID:</b>	1961
<b>UniProt ID:</b>	<a href="#">Q05215</a>
<b>Cytogenetics:</b>	2p13.2
<b>Gene Summary:</b>	Transcriptional regulator. Recognizes and binds to the DNA sequence 5'-GCGGGGCG-3' (GSG). Activates the transcription of target genes whose products are required for mitogenesis and differentiation (By similarity).[UniProtKB/Swiss-Prot Function]