

Product datasheet for **SC327932**

ALS2CR13 (FAM117B) (NM_173511) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ALS2CR13 (FAM117B) (NM_173511) Human Untagged Clone
Tag:	Tag Free
Symbol:	ALS2CR13
Synonyms:	ALS2CR13
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC327932 representing NM_173511.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTCCAGCGGGTGAGGCGCAATGGGTCCCCACGCCGGCTCCCTTGGGGTGGTGCGGTGGCC
ACGGCCGGGGACCCGGGAGCCGCTTGACGCCATGAGGGCGACGGTTCGTTCCAGCTGAAGCAGCAG
CAGCAGCAGCAACATGGCAGCCCCACGCCGAGCGCGCGCGCGGCAACAACACGGTGGCTGC
TGTGGTGGCGCCTCAGGCCCCGAGGCGCGCGCGCGCGGTGGCCCGCACCGCTCGCGCAGCACC
AGCCCCACGCGCGCGCGGGAACGCGCGCGCGCACCGCCCCACGGTGGCCACGACGAGCGCGCG
TCCGCGACGTCCACGCGAGGCACCGCCCCACGCCGAGCGCGCGCTGGAGCTCGCGGGAGCCCCCA
CGGCGCGCGCGCGCGCGCGCTGCTGGGACCGTGTCTGCGCCAGCTCGTCCGCCACCCACCTGTGG
ACCGGCGAGGTGAGCGCGGCCCCACCCACGCCGCGTCCGGCATCGGAGGAGTCTCCGGAGCAGAGC
CGAAGCTCGCCGAGAAGAGGAGCCCCAGCGCCCCGTTTCAAAAGCAGGTGACAAAACACGACGCT
TCTTCAAGCCCCCTCCAGTATTATCCGACGCACTTCTCCTGGTACTCTTGTGCACCGTATCTTGCT
GGACACTGGCCTCGGGATAGCCATGGGCAAGCTGCACCTTGCATGAGGGACAAAGCTACACAGACAGAG
AGTGATGGGCTGAAGAATACTCTGAAAAGAAGAAAGGTCTCACAAGCGCTCAGCATCTTGGGGCAGT
ACAGATCAACTTAAGGAGATTGCAAAATTACGCCAGCAGTTGCAGAGAAGTAAACACAGCAGTCGGCAT
CATCGAGATAAAGAAAGACAGTCTCCATTTTCATGGCAACCATGCAGCTATTAACAGTGTGAGGCTCCT
GTTCCAAAGAGTGCACTTATTCTGTAATTCCTATCACCAGCTCCCGGTTCCGGAATAGC
GTGGAAGGATTGAATCAGGAGATTGAAATAATAATTAAAGAGACTGGGAAAAGGAAGCAACTTATA
CCGCAAGATATTCAGATGGCCATCGTGCTCACCCCCCTTGTACAGAGAAGTAGCAGCAGCGCAGC
ATTGACACACAGACGCCCTGGTGGGCGAGACAGGGAAGCAACAACAGCAGCCGTTCCAGTCCGTGTCC
CCAACATCGTTTCTCACCATTTCGAATGAAGGTAGCGAGGAGAGTCTTGTCTCAGCGGATGACCTGCTT
GTTGATCCAGAGATAAAGAGAATGGGAACAATTCTCTTTGCCAAATATGCAACCTCACCAAAACCT
AACACAGTTATATGTTCAAAAGGGAACCTCCTGAGGGCTGTGAAAGGGTCAAAGTCTTTGAGGAATGC
TCGCCAAAACAGCTTCATGAAATCCAGCCTTTTATTGCTGACAAAAACAAGTGAAATTCATTCTT
AAAAGTGGATCTGCTTTCTGCTCGTCAGCATCTCAAGCCACTCTTCTACCCCGGATCTTACACTC
AAGGGCTCTGGCCACAGCCTGACAGTCACTGCGCATGACAACCACTCTGCTGCAGCCTATTGCTGTG
GCCTCCCTGTCTACAAACACAGAGCAAGACCGAGTCTCTCGAGGAACAAGTACAGTCATGCCATCAGCT
TCTCTACTCCACCACCAGAACCAATTGAGGAAGCTGAAGGATAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_173511

Insert Size: 1770 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM_173511.3

RefSeq Size: 5795 bp

RefSeq ORF: 1770 bp

Locus ID: 150864

UniProt ID: Q6P1L5

Cytogenetics: 2q33.2

MW: 62 kDa