

Product datasheet for **SC327935**

KLHDC7B (NM_138433) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	KLHDC7B (NM_138433) Human Untagged Clone
Tag:	Tag Free
Symbol:	KLHDC7B
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC327935 representing NM_138433.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGTTCTTAGAAGCCACCCCTTCCCCAGGCAAGACAGGCCCAAGGGAGTGTCCCAGGGCGGTTCCC
GGGAGCCCCGTGGGTCCAGCACTTCCACACACTCTGAGGACAGACACGGCCCCCTTCTTCAGTGGGG
ACAGTCATAGGGACAGGTACAGGGGGCCTGGTTGAGGCTGGAGGTGAGCCACAGCCAAGAAAGCTCCGAG
ACCAACGGATCGCCAGCCAGACCCTCCCCAGGCCTAAGAGGAGAGGGAACAGGGAGAAAAGTCTA
GACCCGCTGCCCCAAGCCGCGATGCCAGGGGGCCCCGACAGCCCCCGCGCAGAGGCCGCTGCCCCC
GCGGCCTCCTCTGCGAGGCGCTCACAGCCGGTACCCAGCTACGAAACGCAGCAGGTGCGAAATC
GCCCCGAGCTCGGAGCAGGAGGTGAGCCGGCCGCTCGGGGACCTCAAGGGGAGGCGCCGGGGAG
GGGGGAGCCCTGCCGGCCGAGCGGGGCGCTACGAAAAGCAGGAGGAGGCCCGAAGCTCATGGTG
TTTCTGAGAGGCCCGGGGTTGGGGGTTGGTGGAGGGGCCCCGAAGCCAGCTCCCGGCCCTGGAG
CCCCCACGGCGCAGCCCTGCCGGCGGCGCTGGACCTGGGAGTTGCCTGGACGTGCTGGCCTTTGCC
CAGCAGACGGAGAGCCCGCCCTGGCGCAGGAGACCTACGCGCTGATGAGCGACAACCTGCTGCGAGTG
CTGGGAGACCCGTGCCTTACCGCCGGCTGAGCGCGGCCGACCGCGAGCGCATCCTCAGCCTGCGGACC
GGCCGGGGCCGGGCGGTGCTGGGCGTCTCGTACTGCCAGCCTTACCAGGGGGCCGCTCAGGGCTC
CCCAGGGGCCCTCGTGGCGAGGAGCCTCCTGCCGGGCCCTGTGTCCCTGCCTTACCTGCGCACCTG
CATGTGTTCAACCCCGGGAGAACACCTGGCGGCCCTGACCCAGGTGCCGAGGAGGCCCGCTTCGG
GGCTGCGGTCTCTGCACCATGCACAACTACCTGTTTCTGGCGGGGGCATCCGTGGCTCCGGTGCCAAG
GCCGTCTGCTCCAACGAGGTCTTCTGCTACAACCTCTGACCAACATCTGGAGCCAGGTTCGGCCCATG
CAGCAGGCCGAGCCAGCTCAAGCTGGTGGCCCTGGACGGGTGCTCTATGCCATCGGTGGCGAATGC
CTGTACAGCATGGAGTGCTACGACCCGCGAACAGACGCCTGGACCCACGCGCGCCACTCCCCGACGGC
ACCTTCCCTGTGGCCACGAGGCTGTGGCCTGCCGTGGGACATCTACGTACCGGGGTACCTTTTC
TACCGCCTGCTCAGGTACAGCCCGTGAAGGATGCTTGGGACGAGTGCCATACAGTGCCAGCCACCGG
CGTTCCAGCGACATCGTGGCACTGGGGGCTTCTGTACCGCTTCGACCTGCTGCGGGGCGTGGGCGCC
GCCGTGATGCGCTACAACACAGTGACCGGCTCCTGGAGCAGGGCTGCCTCCCTGCCCTGCCCGCCCCC
GCCCCACTGCACTGCACACCTGGGCAACACCATTTACTGCCTCAACCCCAAGTCACTGCCACCTTC
ACGGTCTCTGGGGGACTGCCAGTTCAGGCCAAGGAGCTGCAGCCCTTCCCTTGGGGAGCACCAGG
GTCTCTCAGTCCATTATCTGACTCTGCCCCCTGAGGACCGGTGCGACCTCACTCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_138433

Insert Size: 1785 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_138433.3](#)

RefSeq Size: 3008 bp

RefSeq ORF: 1785 bp

Locus ID: 113730

UniProt ID: [Q96G42](#)

Cytogenetics: 22q13.33

MW: 63.3 kDa