

Product datasheet for **SC327939**

CCDC48 (EFCC1) (NM_024768) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CCDC48 (EFCC1) (NM_024768) Human Untagged Clone
Tag:	Tag Free
Symbol:	CCDC48
Synonyms:	C3orf73; CCDC48
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC327939 representing NM_024768.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGAGCCGGTCAGCACGGGCGCGGAGGCCGGCATGGAGGGCGCGGGAGGTGACCCGTACCGGCGACCT
GCGCGGCGCACGCAGTGGCTGCTGAGCGCCCTGGCGCACTACGGGCTGGACCGCGCGTGGAGAAC
GAGATCGTGGTGTGGCCACCGGCTGGACCACTACCTGCAGGAGGTCTTCCACCACTGGACTGCCGC
GGCGCGGCGCGTCTGCCCGCGCGCACTTCCGAGCGCTCTGCGCTGTGCTGGGGTGCAGCGGAGGGG
GCCACCACGGCGGGCAGGCAGCAGGTGACGGGAATCCAGAGATGTGACCCCGGGGATGCGGCCGCT
GAGTTGGCCACGGACGGGACTCAGATACCGATGAAGAGCGCGCCTGGCGCTGCGCGCCGAGCCCGG
GAGCTCACCTCCGCCAGTTCCACGCGCGCTCTGTGGCTACTTCGGCACCCGTGCGGGGCCCCGGCTG
CCCCGCGCGCTCTCAGCGAGCACATCGAGACGAGATCCGCCTGCGCCGTCCGCGCCGCCGCCCGC
CCGCCCTGCGCGCTGGCCCCGACAGCGGTCTGACTGTGAGCGGTTGCGCGGTGGAGGAGGAGAAT
AGCAGCTTGCGCGAGTTGGTGGAGGACCTGCGCGCGCGCTGCAGAGCAGTGATGCGCGCTGCCTAGCA
CTGCAGGTGCGACTCTGGAAGAGCCAGGCGAGCACCCACGAGATGGGGCACGGCGGGCGGAGGCTGCG
GTGCGGGAGCTGCGTCAGGCGCAGGCGCGCCTGGCTGCGCGGAGGCCCGCGCTGGGCGGCTGCGCCGT
GGCCAGGCCGAGGTGCGCGCGCGCGGAGGAGGCCCGGCGAGGTGGTGTGCGCAGCCTGCACCGCGTG
CGAGAGCTGGAGGCGCTGGCGCAACAGGTGCCCGGCTTGACGCGTGGGTGCGGCGGCTGGAGGCGGAG
CTGCAACGCTACAGGTACAGGATTCCAGCTCCCAACCCCGCAGCTAGCCAACCCAGAGCCAGGAGAC
AAGAGTAATGAACCTGAAGATGCTGGGACCAGAGACCCAGACCCCACTCCAGAGGGAGCCTGGCAGTCA
GACAGCAGCTCTGAAGCAGAGCCCTGGATGAAGTGGACGAGCAGCTGTTCCGCTCCGTGGAGGGCCAG
GCCGCTCTGACGAGGAGGAGTGGAGGAGGAGGTGGCAGGAGGAGAAGAAGACGCCGCGCAGCCGAG
GCCAAGACACTGCTGGCCCGGCTCTCAGCTGCAGAGGAGGTGTGATGACCAAGAGCGGCGGAGAAGCTC
ATGACTTACTTTGGTCACTTCGGCGGTGCCAACCATGCCCATACCCTGGGGAGCTGGAGGCCTGCATT
GCCATGCTGGTGGAGCAGCTGAGGACTCAGGGCTGCGGTGGGAGGACCCTGGGGACCTCTGAGGAGGAG
GCAGAGTTGCAGCAGAAGGTGGAAGAGAATGAGCACCTGAGGCTGGAGCTGCAGATGGTAGAGACCGAG
AGGGTGGCGTGTCCCTGCTGGAGGAGAAGCTGGTGGACGTGCTGCAGCTCCTGCAGAGGCTCCGGGAC
CTGAACATATCGAAAAGAGCCCTGGGGAAGATTTTGTGAGCACGCTGGACGCTTCAGGACCCCAACC
CACGAGGGAAGGCCAGCCCTGCAGCCATCCTGGATGCCCTGCACCAAGCCTTGCTGCCTGCCAGCTG
TTGCGGAGACAGCCCTCGGCACCAGCCTCTGCAGCAGTGCCTACCAACCCCTCCTGCTCCTCTGC
TGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_024768

Insert Size: 1797 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_024768.2](#)

RefSeq Size: 2691 bp

RefSeq ORF: 1797 bp

Locus ID: 79825

UniProt ID: [Q9HA90](#)

Cytogenetics: 3q21.3

MW: 65.7 kDa