

Product datasheet for **SC327940**

Apc7 (ANAPC7) (NM_016238) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Apc7 (ANAPC7) (NM_016238) Human Untagged Clone
Tag:	Tag Free
Symbol:	ANAPC7
Synonyms:	APC7
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC327940 representing NM_016238.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGACCCCGCGACGCCGCATTTTGGAGTCTTCCCTAAGGATCCTCTACCGGCTTTTCGAGTCAGTG
CTGCCGCCGCTGCCCGCGCTTTCAGAGCAGGATGAATGTGATAGACCACGTGCCGGACATGGCGGCC
CGGGGGCTGCACTCCAAGTGCGGCTCCTCAGCAGCTTGTACTTACAATGAGTAATAACAACCCTGAG
TTATTCTCCACCTCAGAAGTACCAGCTTTTGGTGTATCATGCAGATTCTCTTTTCATGATAAGGAA
TATCGGAATGCTGTGAGTAAAGTATACCATGGCTTTACAGCAGAAGAAAGCGCTAAGTAAACTTCAAAA
GTGAGACCTTCAACTGGAATTCGCATCTACTCCACAAAGTCAGTGTCTTCCATCTGAAATTGAAGTG
AAATACAAAATGGCTGAATGTTATACAATGCTAAAACAAGATAAAGATGCCATTGCTATACTTGTGGG
ATCCCTTCAAGACAAAGAACTCCAAAAATAACATGATGCTGGCAAACCTGTACAAGAAGGCTGGTCAG
GAGCGCCCTTCAGTCACCAGCTATAAGGAGGTGCTGAGGCAGTGCCATTAGCCCTTGATGCCATTCTA
GGCTTGTGTCCCTTTCTGTAAGGGGGCAGAGGTGGCATCCATGACAATGAATGTGATCCAACCGTG
CCTAACTTGGACTGGCTCTCTGTGTGGATCAAAGCGTATGCTTTTGTGCACACTGGTGACAACTCAAGA
GCAATCAGTACCATCTGTTCAGTACAGAAAAAATCCTTATTGCGAGATAACGTGGACCTATTGGGAAGC
TTGGCAGATCTGACTTCAGAGCTGGAGACAATAAAAATCTGTCTCAAGTTTGAACAGGCACAGATG
TTGGATCCTTATCTGATAAAAGGAATGGATGTATATGGCTACCTACTGGCAGGAGAGGGCGGCTAGAG
GATGTTGAGAACCTTGGATGCCGCCTTTTCAATATCTCTGATCAGCATGCAGAACCGTGGGTGGTTTCT
GGCTGTCACAGCTTCTATAGCAAACGCTACTCCGGGCCCTCTATTTAGGAGCCAAGGCCATTGAGCTG
AACAGTAATAGTGTTCAAGCTCTGCTACTTAAGGGAGCAGCACTTAGGAACATGGGCAGAGTCCAAGAA
GCAATAATCCACTTTCGGGAGGCCATACGGCTCGCACCTTGTGCTTAGATTGTTATGAAGTCTTATC
GAATGTTACTTAGCCTCCAACAGTATTCGAGAAGCAATGGTAATGGCTAACACGTTTACAAAATCTG
GGAGCAAATGCACAGACCCTTACCCTTTTAGCCACCGTTTGTCTTGAAGACCCAGTGACACAGGAGAAA
GCCAAAACATTATTAGATAAAGCCCTGACCCAAAGGCCAGATTACATTAAGGCTGTGGTAAAAAAGCA
GAACTACTTAGCAGAGAACAGAAATATGAAGATGGAATTGCTTTGCTGAGGAACGCACTGGCTAATCAG
AGTGACTGTGCTCCTGCATCGGATCCTAGGAGATTTCTTGTAGCTGTCAATGAGTATCAGGAGGCAATG
GACCAGTATAGTATAGCACTAAGTTTGGACCCCAATGACCAGAAGTCTCTAGAGGGGATGCAGAAGATG
GAGAAGGAGGAGAGTCCCACGGATGCCACTCAGGAGGAGGATGTGGACGACATGGAAGGGAGTGGGGAA
GAAGGGGACCTGGAGGGCAGCAGAGTGGGGCGCCAGTGGGCTGACCAGGAGCAGTGGTTCGGCATG
CAGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_016238

Insert Size: 1800 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_016238.2](#)

RefSeq Size: 3070 bp

RefSeq ORF: 1800 bp

Locus ID: 51434

UniProt ID: [Q9UJX3](#)

Cytogenetics: 12q24.11

Domains:	TPR
Protein Families:	Druggable Genome
Protein Pathways:	Cell cycle, Oocyte meiosis, Progesterone-mediated oocyte maturation, Ubiquitin mediated proteolysis
MW:	66.9 kDa
Gene Summary:	<p>This gene encodes a tetratricopeptide repeat containing component of the anaphase promoting complex/cyclosome (APC/C), a large E3 ubiquitin ligase that controls cell cycle progression by targeting a number of cell cycle regulators such as B-type cyclins for 26S proteasome-mediated degradation through ubiquitination. The encoded protein is required for proper protein ubiquitination function of APC/C and for the interaction of APC/C with certain transcription coactivators. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Nov 2008]</p> <p>Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (a).</p>