

Product datasheet for **SC327943**

GOLGA8B (NM_001023567) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GOLGA8B (NM_001023567) Human Untagged Clone
Tag:	Tag Free
Symbol:	GOLGA8B
Synonyms:	GOLGA5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC327943 representing NM_001023567.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGCAGAAGAACTGGACAGAGTAAATTAGCTGCAGCCAAGAAAAAGTTCAAAGAATATTGGCAGAGA
AACCGCCCTGGTGTTCAGCAGCAGCGAAGAGGAACACGAAAGCAAATGGCAGTAGCCCTGAGACGGCC
GCTTCTGGTGGTTGCCACTCATCTGAGGCTTCCTCCTCCGCTCCTCCTCTGCTGCGCGTCAGAGC
CCGTGCCAAGAACAAGCAGCAGTCTGAACTCGAGGTCCATAAAAAATCAGTCGACTGAATGACACCATC
AAATCTTTGAAACAACAGAAGAAACAAGTGAACATCAGCTGGAAGAAGAAAAGAAAGCAAAATGAG
AAACAGAAAGCTGAAAGGGAGCTAGAGGGTCAAATCCAGAGATTGAACACAGAGAAAAAGAACTAAAT
ACGGACCTGTATCACATGAAACATTCTCTCAGATACTTTGAAGAAGAGTCCAAGGATCTGGCCGGCCGC
CTGCAACGTTTCATCGCAGCGTATAGGAGAGTTAGAGTGGTCTCTCTGTGTGTGCGCCGCACACAGAAG
AAGAAGCCGGATGGGTTCTCGAGCCGAGTAAAGCACTTCTCAAGCGGCAGTTAGAGCAGTCCATACGG
GAGCAGATACTGTGAAAGGACAGTGACACAGTTGAAGGAGTCGCTTAAAGAAGTCCAGCTGGAGAGA
GATCAATATGCTGAACAAATAAAAGGAGAGAGGGCCAGTGGCAGCAGAGGATGAGGAAAAATGTCCGAG
GAGGTTTGACATTGAAGGAGGAGAAGAAGCATGATACGCATCGGTAGAGGAGCTGGAGAGGAGCTTG
TCCAGACTCAAAAACAGATGGCTGAGCCACTGCCCCGGATGCCCCAGCAGTGTCTCTGAGGTGGAG
CTGCAAGACCTGAGGAAGGAGCTGGAGAGAGTGGCAGGAGAGCTCCAGGCTCAGGTGGAACAAATCAG
TGCATCAGTCTCTGAACCGTGGGCAAAAGGAGAGGCTTCGCGAGCAGGAGGAGAGGCTTCAGGAGCAG
CAGGAGAGGCTTCGGGAACGGGAGAAGAGGCTTCAGCAGCTGGCCGAGCCACAGAGCGACTTGAGGAG
CTGAAGCACGAGAACAAGAGCGCACTGCAGTTGGAGCAGCAAGTAAAGGAGCTGCAGGAGAAGCTGGGC
CAGGTGATGGAGACGCTCACCTCGGCCGAGAAGGAGCCAGAGGCAGCAGTCCCAGCCTCAGGGACTGGG
GGCGAGTCTAGCGGCCTTATGGACCTCCTGGAGGAGAAGGCGGACCTGAGGGAGCATGTGGAGAACTG
GAACCTGGATTATCCAGTACCGGAGAGAGAGATTGCCATCAGAAAGTACATCGCCTTCTAACAGAGCCA
GGGACAGTGCCAAAGATGCGTCACCGGGAGGAGGCCATCATCAGGCTGGCCAGGACAAGGAGGAGAG
GAAGGTGAAGCTGCTGGAGCTGCAGGAGATGGTGTTCGCGCTTGTGGCAGCTACAGCGAGGGGCACGGC
AAATTCCTGGCGCTGCCCGGAACCTGCTGCTGAACCCAGTCCAGGAGCCCCAGCCCCCAGGAGCTC
GGGGCTGCCGACAAGCATGGTGATCTTTGTGAGGCGAGCCTACCAACAGCGTGGAGCCTGCACAAGGA
GAAGCCAGGGAGGTTCTTCCAGGACAACCTACTGCACAGCCAGTCTGCAGTCTCTTGGTGAGATG
CAGGACCACCAGGAGACCCAGGCTTGGGCAGCAACTGCTGTGTGCCATGCTTTTGTGGGCTTGGCTG
CCGAGAAGAAGGAGATAA
ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: Sgfl-MluI

ACCN: NM_001023567

Insert Size: 1812 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001023567.4](#)

RefSeq Size: 4271 bp

RefSeq ORF: 1812 bp

Locus ID: 440270

UniProt ID: [A8MQT2](#)

Cytogenetics: 15q14

MW: 67.3 kDa

Gene Summary: May be involved in maintaining Golgi structure.[UniProtKB/Swiss-Prot Function]
 Transcript Variant: This variant (1) represents the shorter transcript and encodes the functional protein. Sequence Note: The RefSeq transcript and protein were derived from genomic sequence to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on alignments.