

Product datasheet for **SC327948**

LRRC70 (NM_181506) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LRRC70 (NM_181506) Human Untagged Clone
Tag:	Tag Free
Symbol:	LRRC70
Synonyms:	SLRN
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC327948 representing NM_181506.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTGTGGATTACAGTTTTCTCTGCCTTGCCTACGACTGTTTCTGGTTGTTACCTGTTATCTTTTATTA
TTACTCCACAAAGAAATCTTGGATGTTCTGTCTGTTTGTACAGCTCTGCACTGGGAGACAAATTAAGTGC
CGTAACCTTAGGCCTTTTCGAGTATTCCTAAGAATTTTCTGAAAGTACAGTTTTTCTGTATCTGACTGGG
AATAATATATCTTATATAAATGAAAGTGAATTAACAGGACTTCATTCTCTGTAGCATTGTATTTGGAT
AATTCTAACATTCTGTATGTATATCCAAAAGCCTTTGTTCAATTGAGGCATCTATATTTCTATTCTA
AATAATAATTTTCATCAACGCTTAGATCCTGGAATATTTAAGGGACTTTTAAATCTTCGAATTTATAT
TTACAGTATAATCAGGTATCTTTTGTCCGAGAGGAGTATTTAATGATCTAGTTTCAGTTCAGTACTTA
AATCTACAAAGGAATCGCCTCACTGTCTTGGGAGTGGTACCTTTGTTGGTATGGTTGCTCTTCGGATA
CTTGATTTATCAACAATAACATTTTGGAGATATCAGAATCAGGCTTTCAACATCTTGAAAACCTTGCT
TGTTTGTATTTAGGAAGTAATAATTTAACAAAAGTACCATCAATGCCTTTGAAGTACTTAAAGTCTT
AGAAGACTTTCTTGTCTCATAATCCTATTGAAGCAATACAGCCCTTTGCATTTAAAGGACTTGCCAAT
CTGGAATACCTCCTCTGAAAAATTCAAGAATTAGGAATGTTACTAGGATGGGTTTAGTGAATTAAT
AATCTTAAACATTTGATCTTAAGTCATAATGATTTAGAGAATTTAAATCTGACACATTCAGTTTGTTA
AAGAATTTAATTTACCTTAAGTTAGATAGAAACAGAATAATTAGCATTGATAATGATACATTTGAAAAT
ATGGGAGCATCTTTGAAGATCCTTAATCTGTCAATTAATAATCTTACAGCCTTGCATCCAAGGTCCTT
AAGCCGTTGTCTTCATTGATTCATCTTCAGGCAAATCTAATCCTTGGGAATGTAAGTCAAACTTTTG
GGCCTTCGAGACTGGCTAGCATCTTCAGCCATTACTCTAAACATCTATTGTCAGAATCCCCATCCATG
CGTGGCAGAGCATTACGTTATATTAACATTACAAATTGTGTTACATCTTCAATAAATGTATCCAGAGCT
TGGGCTGTTGTAATAATCTCCTCATATTCATCACAAAGACTACTGCGCTAATGATGGCCTGGCATAAAGTA
ACCACAAATGGCAGTCCTCTGAAAAACTGAGACTGAGAACATTACTTTCTGGGAACGAATTCCTACT
TCACCTGCTGGTAGATTTTTTCAAGAGAATGCCTTTGGTAATCCATTAGAGACTACAGCAGTGTTACCT
GTGCAAAACAACTTACTACTTCTGTTACCTTGAAGTGGAAAAAACAGTGCTCTACCGAATGATGCT
GCTTCAATGTCAGGGAACATCTCTAATTTGTACACAAGAGTTGAGAAGTTGAATGAGGCTTTTGAC
ATTTTGCTAGCTTTTTTCATCTTAGCTTGTTTAAATCATTTTTTGTCTACAAAGTTGTTCAAGTTT
AAACAAAACTAAAGGCATCAGAAAACCAAGGAAAAATAGACTTGAATACTACAGCTTTTATCAGTCA
GCAAGGTATAATGTAAGTGCCTCAATTTGTAACACTTCCCAAAATCTCTAGAAAGTCTGGCTTGGAG
CAGATTGCACTTATAACAAATTGTTCTGAAATGAGGCACAGGTATTCTTTTGAACATTCTGCT
TTATAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: Sgfl-MluI

ACCN: NM_181506

Insert Size: 1869 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_181506.4</u>
RefSeq Size:	2290 bp
RefSeq ORF:	1869 bp
Locus ID:	100130733
UniProt ID:	<u>Q7Z2Q7</u>
Cytogenetics:	5q12.1
MW:	70.3 kDa
Gene Summary:	Renders cells highly sensitive to the activation by cytokines and lipopolysaccharide (LPS). [UniProtKB/Swiss-Prot Function]