

Product datasheet for **SC327956**

ARHGAP25 (NM_001007231) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ARHGAP25 (NM_001007231) Human Untagged Clone
Tag:	Tag Free
Symbol:	ARHGAP25
Synonyms:	HEL-S-308; KAIA0053
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC327956 representing NM_001007231.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGTCCTAAAAATTGCCAAGGAAGTGGGATTTCAACCTGAAAGTGGAGGCTGCGAAAATAGCTCGGTCA
AGGAGTGTGACTGGCGAGCAGATGGCTGCCTTCCATCCATCGTCCACCCCAACCCGCTGGAGAGG
CCCATCAAGATGGGCTGGCTGAAGAAGCAGAGGTCATCGTGAAGAACTGGCAGCAGAGGTAATTTGTG
CTGAGGGCGCAGCAGCTCTACTACTACAAGGATGAAGAGGACACGAAGCCCAAGGCTGCATGTATCTA
CCAGGATGTACAATCAAGGAGATCGCCACAACCAGAAGAAGCTGGGAAGTTTGTCTTTGAAATCATT
CCAGCCTCATGGGACCAGAATCGCATGGGACAGGACTCCTATGTCCTCATGGCCAGCTCTCAGGCGGAG
ATGGAGGAGTGGGTTAAATTCCTCAGGAGAGTTGCTGGCACACCCCTGTGGAGCAGTGTGGCCAGCGC
TTGGATGAGACTGTGGCCTATGAACAGAAATTCGGCCCCATCTGGTGCCATCTGGTGGAGAAATGT
GCAGAGTTCATCCTGGAGCACGGCCGAATGAAGAGGGCATCTCCGTCTGCCTGGGCAGGACAACCTG
GTGAAGCAGCTGAGAGACGCTTTTGTGCTGGGAGCGGCCCTCCTTTGACAGAGACACAGATGTGCAC
ACTGTGGCTTCCCTGTTAAAGCTCTACCTCCGAGACCTCCAGAGCCCGTGGTCCCTGGAGCCAGTAC
GAAGGGTTCCTGCTCTGTGGGCAGCTCACGAATGCGGATGAGGCAAGGCTCAGCAGGAGTTGATGAAG
CAGCTCTCCATCCTTCTCGTGACAATAAGTCTCCTGAGCTACATCTGCAGGTTCTACATGAAATA
CAGCTGAACTGTGCTGTTAAACAAGATGAGTGTGGACAACCTGGCTACTGTGATTGGTGTGAATCTCATC
AGGTCGAAGGTCGAAGACCCTGCCGTGATCATGAGAGGGACTCCTCAGATCCAAGAGTGTGACTATG
ATGATCAGAGACCATGAAGTCTTCCCAAGTCCAAGGATATACCCCTGTCACCCCTGCCAGAAA
AATGACCCCAAGAAAGCTCCAGTGGCCGAAGCTCTGTAGGCTGGGATGCCACTGAAGACCTCCGAAT
TCTAGGACAGACAGCTTCAGTAGCATGACAAGCGACTCTGATACAACCAGCCCAACCCGACAGCGCG
AGCGATGCGTTTTCCGGAGGACAGCAGCAAAGTACCCAGGAAAAGCCAGGAGACTGGAAAATGCAATCT
CGTAAAAGGACTCAAACACTCCCTAACCAGAAATGTTTCTTGACATCAGTTTTTCAGGGTCCAACAGC
AGCAAAAATGGAGATCTTTAAAAATGAATTCGTGTCGCTTCTCAGAGGCTAAGGCAGGGGAAGGGCAC
AGGAGAACGATGTCTCAAGACTTGCGCCAACCTTCTGACTCCCAACGGACTTCCACCTACGATAACGTC
CCTTCCCTGCCAGGGTCCCCTGGGGAGGAAGCCAGTGCCTCTTCCCAAGCCTGTGACTCCAAGGGA
GATACTCTTGCCAGTCCAAACTCTGAACTGGGCCTGGAAAAAAGAACTCTGGAGAAGAGGAAATGAT
TCTTTGCAGAGGATGGTCCAAGAGCTACGAAAGGAAATAGAAACACAGAAGCAAATGTATGAGGAACAG
ATTA AAAACCTTGAGAAGGAAAATATGACGTTTGGGCTAAAGTGGTGGGCTCAATGAAGAACTGGAG
AAGGAAAAGAAGAAGTCTGCAGCCCTAGAGATCAGCCTCCGCAACATGGAGCGCTCCCGGAGGATGTT
GAGAAGAGGAACAAGGCCTTGAAGAAGAAGTCAAGGAATTTGTCAAATCCATGAAGGAACCCAAGACC
GAGGCTTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001007231

Insert Size: 1941 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001007231.2</u>
RefSeq Size:	3020 bp
RefSeq ORF:	1941 bp
Locus ID:	9938
UniProt ID:	<u>P42331</u>
Cytogenetics:	2p13.3
MW:	73.5 kDa
Gene Summary:	<p>ARHGAPs, such as ARHGAP25, encode negative regulators of Rho GTPases (see ARHA; MIM 165390), which are implicated in actin remodeling, cell polarity, and cell migration (Katoh and Katoh, 2004 [PubMed 15254788]).[supplied by OMIM, Mar 2008]</p> <p>Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (a).</p>