

Product datasheet for **SC327959**

MON1A (NM_032355) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MON1A (NM_032355) Human Untagged Clone
Tag:	Tag Free
Symbol:	MON1A
Synonyms:	SAND1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC327959 representing NM_032355.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCACCCCGGGGCGGGCCGAGTCGGGCTGAGCGGCTGGAGCTGGCCCTTGGCCGCGAGAGCCGGCC
AAAGCTATCTTCTGCACCGCGCGGGTGGGGCGGTGGGAGGAGCGGTGCCTGCGCTGTGGCCAT
GTCTGTGTACGCAGGGGGCCGGGGCCCGGAAGCTGTGCCCTCAGGCCGGCCCGCGGACACCCCTT
ACCCCTCCATGGGTCCGGCAGAGGGCAGTGACCGGCACCTTCTGTGCGTCTGGACTCCTCTGAGAAAC
AGGAGAGCTCAAAGGATGGTACTGACATGCAGAGGAAGAGAAGCAGCAATGCCTTGATGGCACATTG
ACTCCTTCTGATGGACAGAGTATGGAGAGAGCTGAGAGCCCCACACCAGGAATGGCCAGGGAATGGAG
CCAGGTGCGGGCCAGGAGGGTGCATGTTCTGTCATGCCGTTCTACGAGGACCTGACTGAGTCAGAG
GATGGGGCAGCTTCTGGGACAGCCACAAGGAGGGTACCAGGGTCCCCCGCGTGCCTACAGACATG
CGCCAGATCAGCCAGGACTTTAGCGAGCTAAGCACCCAGCTGACGGGTGTGGCCCGGACCTGCAGGAG
GAGATGCTGCCAGGAAGCTCTGAGGATTGGCTGGAACCCAGGGGCAGTTGGGCGACCAGCCACAGAG
CCCCCAGGGAGGGCACAAACGAGGGGGATGAGGAGGATGCCACGGAGGCATGGCCCTGCACCAGAAG
CATGTCTTTGTGCTGAGTGAGGCAGGGAAGCCTGTGTAATCCCGTATGGGTCTGAGGAGGCACTTTCC
AGCACTATGGGTGTTATGGTGGCCCTGGTGCCTTCTGGAGGCAGACAAGAAGCCATCCGCTCCATC
CATGCAGATGGCTACAAGGTAGTATTCGTGCGCCGGAGCCCGTGGTGCTAGTGGCGGTGGCTCGTACG
CGGCAGTCGGCACAAGAGCTGGCGCAGGAGCTGCTCTACATCTACTACCAGATCCTAAGCCTTCTTACC
GGTGCGCAGCTGAGCCACATCTTCCAGCAGAAGCAGAATATGATTTGCGGCGCCTACTCTCGGGCTCA
GAGCGCATCACCACAACCTGCTGCAGCTCATGGCAGAGACCCAGCTTCTGATGGGGCGGCACGG
TGCTGCCCTGGCGGGCGCGTGCAGCAGACTGTGAGCGCCAGCCTGCAGCAGGCGCGTGCAGCAGC
CTGGTCTTCTCCATCCTGCTGGCCCGCAACCAGCTCGTGGCACTCGTGGCCGAAAGGACCAATTTCTG
CACCCCATCGACCTGCACCTGCTCTTCAACCTCATTAGTTCCTCCTCGTCTTTTCGCGAGGGCGAGGCC
TGGACGCCGTGTGCCTGCCAAATTCACGCAGCCGGCTTCTTCCAGCACACATCTCTTACCTAGAG
CCTGACACTGACCTCTGCCTGCTGCTTGTCTCCACTGACCGTGAGGACTTCTTTGAGTCTCTGACTGC
CGCCGCCGCTTCCAGGAGCGCTTCGAAGCGGGAGCCACCTGGCCCTGCGAGAGGCACTGCGCACA
CCCTACTACAGCGTTGCCAAGTGGCATCCCTGACCTGCGTCACTTCTCTATAAGTCAAAGAGCTCG
GGACTCTTACCAGCCCTGAGATTGAGGCCCATACACCAGTGAAGAGGAGCAGGAGCGGCTGCTGGGC
CTCTACCAGTACTTGACAGTCGTGCCACAATGCCTCTCGCCACTCAAGACCATTTACTACACGGGC
CCCAACGAGAACCTCCTGGCCTGGGTGACAGGCGCTTTGAGCTTACATGTGTTACAGCCCCCTGGGG
ACCAAGGCGTCAGCCGTAGTGCATCCATAAGCTGATGCGCTGGATCCGCAAGAGGAAGACCCGCTC
TTCATTCTCACGCCCTCACCTATGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_032355

Insert Size: 1959 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_032355.3](#)

RefSeq Size: 2244 bp

RefSeq ORF: 1959 bp

Locus ID: 84315

UniProt ID: [Q86VX9](#)

Cytogenetics: 3p21.31

MW: 72.9 kDa

Gene Summary: Plays an important role in membrane trafficking through the secretory apparatus. Not involved in endocytic trafficking to lysosomes (By similarity). Acts in concert with CCZ1, as a guanine exchange factor (GEF) for RAB7, promotes the exchange of GDP to GTP, converting it from an inactive GDP-bound form into an active GTP-bound form (PubMed:23084991).
[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (a).