

Product datasheet for **SC327960**

ZNF630 (NM_001037735) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF630 (NM_001037735) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF630
Synonyms:	dj54B20.2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC327960 representing NM_001037735.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGATTGAGTCCCAGGAACCACTGACATTTGAGGATGTGGCTGTGGACTTCACGCAGGAAGAGTGGCAG
CAGTTGAATCCTGCTCAGAAAGCCCTGCATAGGGATGTGATGCTGGAGACCTATAATCACCTGGTCTCC
GTGGGGTGTTCAGGTATAAAACCAAGATGTAATCTTTAAGTTGGAACATGGAAGGACCCATGGATCATA
GAGAGTGAGTTGTCAAGGTGGATCTACCCAGACAGAGTGAAAGGCCTTGAATCTTCCAGCAGATCATT
TCTGGAGAACTTTTATTTCAAAGGGAGATACTAGAAAGAGCCCCAAAGGATAATTCATTGTACTCTGTT
TAAAAATCTGGCATATTGATAATCAGATGGATAGATATCAAGGAAATCAAGACAGAGTTTGGAGCAG
GTCACAGTCATCAGTCGTGAAACATTGACTGATGAGATGGGTTCCAAGTACAGTGCATTTGGGAAAATG
TTCAATCGGTGCACAGACCTTGCTCTTTAAGTCAAAAATTCATAAGTTTGATTCATGTGAAAATAGC
TTGAAGTCTAATTCAGACTTACTAAATTATAACAGGAGCTATGCAAGAAAGAACCCCACTAAGAGATTT
AGATGTGGGAGACCACCTAAGTATAATGCTTCTGTTCTGTGCCTGAGAAGGAAGGCTTCATTCACT
GGAATGGAGCCCTATGGAGATAGTCAATGTGAAAAAGTTCTCAGTCATAAGCAAGCCCATGTTCACTAT
AAGAAATTTCAAGCCAGAGAGAAACCCATGTTTGTAGTATGTGTGGAAAGCCTTTATCAAGAAGTCA
CAGCTCATTATACATCAAGAATTCATACTGGAGAGAAACCATATGTATGTGGAGATTGTAGGAAAGCC
TTCAGTGAGAAATCACACCTCATTGTGCATCAGAGGATTCATACTGGGGAGAAACCCATGAATGTACT
AAGTATGGAAGAGCATTCTCCCGAAGTCACCTTTCACTGTTTCATCAGAGAGTCCATACTGGAGAGAAA
CCCTATGAGTGTTTTGAGTGTCAAAAGCTTTCTCCAGAAGTCACATCTAATTATACATCAGAGAGTT
CATACCAGAGAGAAGCCCTTTGAATGCAGTGAATGCAGGAAAGCCTTCTGTGAGATGTCTCACCTTTTT
ATACACCAGATAACTCATACTGGGAAGAAGCCCTATGAATGTACTGAATGTGGGAAGACCTTCCCTCGG
AAAACACAGCTCATTATACATCAGAGAACGCATCTGGAGAGAAGCCCTATAAGTGTGGTGAATGTGGG
AAAACCTTTCTGCCAACAGTCCCACCTCATAGGACATCAAGAATTCATACAGGAGAAAAACCTTATGTG
TGTAAGTACTGTGGGAAGGCCTTTTCCAGAAGTCACACCTCACTGGCCATCAAAGACTTCATACTGGA
GAGAAACCTTATATGTGTAAGTGTGGAAAATCCTTCTCTCAGAAATCACCTCTTATCATACACCAG
AGAATTCATACAGGGGAGAAACCTTATCAGTGTGGTGAATGTGGCAAACCTTCTCCAGAAATCACTC
CTCATTATTCATCTGAGAGTTCACACAGGGGAGAAACCTTATGAGTGTACTGAGTGTGGGAGGCGCTTT
TCCCTGAAGTCACATCTCATTCTACATCAGAGAGGTCATACTGGAGAGAAACCTATGAATGTAGTGA
TGTGGAAAGGCCTTCTGTGGAAAGTCTCACTCATTATACATCAGAAACTCATCTAGGAGAGAAACC
CCTGAATGTGCTGAGTCTGGAATGACTTTTTCTGGAAATCACAGATGATTACATATCAGAGAAGACAC
ACTGGGGAGAAACCTCCAGATGCAGTGACTGTGGGAAGGCATTCTGCCAGCATGTATACTTTACTGGG
CATCAGAATCCATATAGGAAAGACACCTTGATATATGCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001037735

Insert Size: 1974 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001037735.3</u>
RefSeq Size:	2529 bp
RefSeq ORF:	1974 bp
Locus ID:	57232
UniProt ID:	<u>Q2M218</u>
Cytogenetics:	Xp11.23
MW:	76.1 kDa
Gene Summary:	This gene encodes a protein containing an N-terminal Kruppel-associated box-containing (KRAB) domain and 13 Kruppel-type C2H2 zinc finger domains. This gene resides on an area of chromosome X that has been implicated in nonsyndromic X-linked cognitive disability. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2013] Transcript Variant: This variant (1) represents the shortest transcript and encodes the longest isoform (1). Variants 1 and 3 encode the same isoform (1).