

Product datasheet for **SC327965**

GOLGA6L1 (NM_001001413) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GOLGA6L1 (NM_001001413) Human Untagged Clone
Tag:	Tag Free
Symbol:	GOLGA6L1
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

Fully Sequenced ORF:	>NCBI ORF sequence for NM_001001413, the custom clone sequence may differ by one or more nucleotides <pre> ATGCTGATGTGGCCCAACCCACCTCCCTACCCATCCCCACCTCCCTACCCATCCCCAC CTCCCTACCCATCCCCACCTCCCTACCCATCCCCACCTCCCTACCCATCCCCACCTCCCT ACCCATCCCCACCTCCCTACCCATCCTATGATGTCCAAAGAAACCAGACAGCAAAATTG GCCGAGGCCAAGGAACAGTTGACAGACCATCATCCCCAGACCAACCCCTAGTGTGGTACA GCAGCAAGCGACACCAAAAAGAAGAAAATAAATAATGGCACTAACCCCTGAGACAACCACT TCTGGTGGTTGCCACTCGCTGAGGATGAACAGAAGGCAAGCCACCAACATCAGGAAGCC CTAAGGAGGGAGCTAGAGGCCAGGTTCTATACCATACGAATCCTTACATGTCAGAAAACT GAGCTTCAGATGGCACTCTACTACAGCCAGCATGCTGTCAAGCAGTTGGAAGGAGAGGCC AGGGATCTGATCAGCCGCTGCATGATTCATGGAAGTTTGCAGGAGAGTTAGAGCAGGCT CTCTCTGCTGTCGCTACACAGAAGAAGAAGGCGGATAGGTACATTGAGGAGTTAACAAAG GAGAGGGACGCCCTGAGTCTGGAAGTGTACAGGAACACCATAACTGATGAGGAGCTGAAG GAGAAAAATGCCAACTACAAGAAAAACTTCAACTTGTAGAATCTGAAAAGTCTGAGATC CAGCTCAACGTAAGGAGCTAAAAAGGAACTGGAGAGGGCCAAGCTCCTGCTGCCACAG CAGCAGCTGCAGGCGGAGGCTGACCACCTGGGTAAAGAGCTGCAGAGTGTGTGAGCAAG CTCCAAGCCCAGGTGGAAGAGAACGAGTTGTGGAACCGCCTGAACCAGCAACAGGAGGAG AAGATGTGGAGGCAGGAGGAGAAGATACAGGAGTGGGAGGAGAAGATACAGGAGCAGGAG GAGAAGATACGGGAGCAGGAGGAGAAGATACGGGAGCAGGAGGAGAAGATGCGGAGGCAG GAGGAGATGATGTGGGAGAAGGAGGAGAAGATGCGGAGGCAGGAGGAGATGATGTGGGAG AAGGAGGAGAAGATGCGGAGGCTGGAGGAGATGATGTGGGAGAAGGAGGAGAAGATACGG GAGCTGGAAGAGAAGATGCACGAGCAGGAGAAGATACGGGAGCAGGAAGAGAAGAGGCAG GAGGAGGAGAAGATACGGGAGCAGGAGAAGAGGCAGGAGCAGGAGGCGAAGATGTGGAGG CAGGAGGAGAAGATACGGGAGCAGGAAGAGAAGATACGGGAGCAGGAGAAAAAGATGTGG AGGCAGGAGGAGAAGATTACGAGCAGGAGAAGATACGGGAGGAGGAGAAGAGGCAGGAG CAGGAGGAGATGTGGAGGCAGGAGGAGAAGATAAGGGAGCAGGAGGAGATATGGAGGCAA AAGGAGAAGATGCACGAGCAGGAGAAGATACGGAAGCAGGAGGAGAAGGTGTGGAGGCAG GAGGAGAAGATGCACGACCAGGAGGAGAAGATACGGGAGCAGGAGGAGAAGATGTGGAGG CAGGAGGAGAAGATAAGGGAGCAGGAGGAGAAGATACGGGAGCAGGAGGAGAAGATACGA GAGCAGGAGGAGATGATGCAGGAACAGGAAGAGAAGATGGGGAGCAGGAAGAGAAGATG CAAGAACAGGAGAAGATGCGGAGGCAGGAGGAGAAGATAAGGGAGCAGGAGGAGAAGATA CGGGAGCAGAAGGAGAAGATACGAGAGCAGGAGGAGAAGATATGGGAGCAGGAGGAGAAG ATACGAGAGCAGGAGGAGATGATGCAGGAACAGGAAGAGAAGATGTGGGAGCAGGAGGAG AAGATGTGTGAGCAGGAAGAGAAGATGCAAGAACAGGAGGAGAAGATGCGGAGGCAGGAG GAGAAGATGTGGGAGCAGGAAGTGAGGCTGCGGCAGCAGGAGGAGAAGATGCAGGAACAC CAGGAGCACCTGGAAGCTGCCATC </pre>
Restriction Sites:	Please inquire
ACCN:	NM_001001413
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM_001001413.3, NP_001001413.3

RefSeq Size: 3822 bp

RefSeq ORF: 2007 bp

Locus ID: 283767

UniProt ID: Q8N7Z2

Cytogenetics: 15q11.2