

Product datasheet for **SC327971**

C8orf45 (MCMDC2) (NM_173518) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	C8orf45 (MCMDC2) (NM_173518) Human Untagged Clone
Tag:	Tag Free
Symbol:	C8orf45
Synonyms:	C8orf45
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC327971 representing NM_173518.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGTCAAATCTAAAAATGAAAGAGGCGGCCCTCATCTATCTTGACAGAAGTGGAGGCCTCCAAAAGTTT
ATAGATGATTGCAAGTACTACAATGATTCAAAACAAAGCTATGCTGTCTATCGATTCAAAATTTAATA
AATCCCTCTGATGTTGTTGAATTAGATGCTGAGCTTGAAATCACATTTTACACCAACCTTTAAAAGCT
GCTGAAGTCTTTCAATCAGTCTGTTTTATTGCTGTTAAGACTCTCTCATTAAATTGGACAATTGCAGACT
GAAACGCAAAATTAATATAGTGCTGAAATTAACACATTTACCTCCCCTGCCAAGTTATGGTCTTGATCTT
TGTGAGTTTCCACTTGATTATACATCTCAGAGATTTTATATGATGCAAGGAATTGTGATTGCAATGACA
ACTATAACCAAGTATACACAAGGGGCAAGATTTCTTTGTCAGATGAAGCATGCCCTCTTCAAAAGGA
TTTCAGTATATAAGAGTGCATGTGCTGGTGCTACAGAATCTGCAACGATAAGAAATGACTTTTTGTGT
AATCTATGTGCATCTTCACTTCAAGAAGACAGAAAATTTAGAGTACTTGGTGATAACAAATAGTTGAA
ATAATTGCCCAAAGGCCTTCGTGCTTTCAAGGATATTCTAACACCAGCCATTTAGGTTTCAATCA
CTTACAATTTTCCCTAAGAGATGAATCAGTGAATAAAATGAATATAGGAAATGAATATAAAATTATTGGA
ATTCCAACCTGTGTAACCACTCACAACTGCTGTCTGTATAGAAGCAAAATAGCATAACTTTTTGTAAAT
TCAAAAGTTCTTCAGGAATTAGTGACAACTTCAGGTGTCTCCTCTCTTAACCTCCAGCTCATGCTGG
AAGTTTACAGCAATACTTGCCAATATCTTTGCATCACAAATTACCCCTCTGGGACTTACAATTTGCTC
AAGCTCTGTTTGTGATGAGTCTAGTACAGACAACTGACCGTAACAAGGAAGTGGAGATTGCCTGGAT
ATTTTAATTATAACAAGTGATACTCTACTCATAGACAGGCTTCTGAATTTTAGCATAAACCTTGCCCC
CGTGGTATACGTCATCTAGTCTCTACTGAAATTTTCCCACTCTATCCAGGAATAAGTATGGAAGTGA
GCAGTTAGCATTACAGGCTGGCAGTGCTTTGCTAGCTAAAGGTGGTATCTGCTTTATAGGAGACTTGGCT
TCACACAAAAAGATAAACTTGAACAGCTTCAACAGTTCTGGAGAGCAGAAGCATCACAGTGTACATC
CCAGGAAGAAGTTTGGGAGGATATTGATCAACAGATGACTTTTCCAGTTCAGTGCAGTTTTTGGTCT
TTTGTGATGTGGATTCTCTCAAGGAGAAATGCACAGAAAATCAACACTCTAATTGGTCAGATGGAT
TGCAGTTTGATTCCAGCTAACCTTGTTGGAGGCATTTGGTTTATTGATTAAGTGAACGAGTCATCTCCC
TGCCACCCATTTCTCTACTGTGCAACACACTTTGAACAAAGCCATTAATCCTGAAGGGCTGTTTTAT
GCGGCTTCTAGACAGTTCACAACTGAAGATTTTGAAAAGCTATTGGCTTTTGCAAGAATTTGAATGTA
GAATTCAGCTTGAAGCAGAAAGAATGACCCATGGCTATTATCTAGCAAGTCGAGAATCAGAACAGGC
TCTGTATGTGGATCAAAGCTGTGAGCATCTGCATTAATAATCTCGTTTTCTATCTGAAGCCCATGCA
CGACTGAACCTTAAGGAACAAAGTCTTAAAGAAGATGTGCTGATTGCAGCCTTACTATTTGAAACATCC
CTCACATTGAAATATGGGGCCACTGTATTTGTGTTGCTCCGAATGCAGTATTTCCATTTGAACTGTAT
AATGAAGAATACTTAGAGCAAAGGGATCTCTATCTGACTCAGTGCCAACAACAGCTGGAACAATTTATT
GCCACATATGGACCTGGGACAATTTTTCTCTAGTGATGAATAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_173518

Insert Size: 2046 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_173518.4</u>
RefSeq Size:	5114 bp
RefSeq ORF:	2046 bp
Locus ID:	157777
UniProt ID:	<u>Q4G0Z9</u>
Cytogenetics:	8q13.1
MW:	76.2 kDa
Gene Summary:	Plays an important role in meiotic recombination and associated DNA double-strand break repair.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (a). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The extent of this transcript is supported by transcript alignments.