

## Product datasheet for **SC327978**

### MID2 (NM\_052817) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	MID2 (NM_052817) Human Untagged Clone
Tag:	Tag Free
Symbol:	MID2
Synonyms:	FXY2; MRX101; RNF60; TRIM1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

**Fully Sequenced ORF:** >SC327978 representing NM\_052817.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGGTGAAAGCCAGCCTCCGTGGTTCTTAATGCCTCAGGAGGACTATTTTCACTAAAGATGAAACA
CTGGAGCTGAATTGACCTGTCCAATCTGCCTAGAGTTGTTTGAAGACCCCTTCTGCTCCCTTGCT
CACAGCCTCTGCTTCAGCTGTGCCATCGCATTTTGGTATCAAGCTGCAGCTCTGGTGAATCCATTGAA
CCCATTACTGCTTCCAGTGTCTACCTGCAGGTATGTTATCTCGCTGAACCACCGGGCCCTGGATGGC
CTCAAGAGGAATGTGACTCTGCAGAACATTATTGATCGCTTCCAGAAGGCTTCACTCAGTGGGCCAAT
TCCCTAGTGAGAGCCGCCGGAAAGGACTTACAGGCCACCACTGCCATGTCTAGCGAGCGAATTGCT
TGCCAATTCTGTGAGCAGGACCCGCCAAGGGATGCAGTAAAAACATGCATCACCTGTGAGGTCTCTAC
TGTGACCGTTGCCTGCGGGCCACGCACCCCAACAAGAACTTTACCAGCCACCGCCTGGTGAACCA
GTGCCAGACACATCTCGAGGGATCACCTGCCTGGACCATGAGAATGAGAAAGTGAACATGTACTGT
GTATCTGATGACCAATTGATCTGTGCCATTATGCAAACTGGTGGTTCGTCACCGAGACCATCAGGTCGA
TCCTGAAATGATCGATTTGAGAACTCAAGCAACTCTGGAGATGAACCTCACCACCTGGTTAAGCGC
AACAGCGAACTAGAAAATCAAATGGCCAACTAATACAGATCTGCCAGCAGGTTGAGGTGAATACTGCT
ATGCATGAGGCCAAAATTTATGGAAGATGTGACGAGTTGGTAGAGATCATCCAGCAGAGGAAGCAAATG
ATCGCTGTCAAAATCAAAGAGACAAAGTTATGAACTGAGAAAGTTGGCACAGCAGGTTGCTAATTGC
CGCCAGTGTCTTGAACGGTCAACAGTCTCATCAACCAAGCTGAGCATATCCTGAAAGAAAATGACCAG
GCACGGTTTCTACAGTCTGCAAAAAATATTGCTGAGAGGGTCGCTATGGCAACTGCATCTTCTCAAGTT
CTGATTCAGACATCAATTTTAAATGATGCCTTTGAAAACCTTTCCTTTAGATTTTTCCAGAGAAAAGAAA
CTCTAGAGGGGTTAGATTATTTAACAGCCCCAAACCCACCATCTATCCGAGAAGAACTCTGACTGCC
TCCCATGACACCATTACAGTCCACTGGATCTCGGATGATGAGTTCAGCATCAGCTCCTATGAGCTTCAG
TACACCATATTCAGTGGCCAGGCTAACTTCATCAGCCTGTATAATTCAGTAGACAGCTGGATGATTGTT
CCCAACATTAACAGAACCATTACACAGTGCATGGACTCCAGAGCGGGACTCGCTACATCTTCATCGTT
AAAGCCATAAACCAAGCCGGCAGCCGGAACAGTGAACCTACCCGACTAAAAACAACAGCCAACCCCTT
AAATTGGATCCAAAATGACTCACAGAAGTTGAAGATCTCCAATGATGGATTGCAGATGGAGAAGGAT
GAAAGCTCTCTAAAGAAGAGCCACACCCAGAGAGGTTTAGTGGCACAGGGTGCATGGGGCAGCAGGA
AATATATTCATTGACAGTGGCTGCCACTATTGGGAGGTGGTCATGGGTTCTCAACATGGTATGCAATT
GGCATTGCCTACAAATCAGCTCCAAAGAATGAATGGATTGGCAAGAATGCCTCCTCATGGGTCTTCTCT
CGCTGCAATAGTAACTTCGTGGTGGAGACACAACAAGGAAATGCTGGTGGATGTGCCCCACACCTG
AAGCGTCTGGGTGTCTCTGGATTATGACAACAATATGCTGTCTTTCTATGACCCAGCTAACTCTCTC
CATCTTCACTTTTGTGATGTGACCTTATTCTTCCAGTTTGTCCAACATTTACAATCTGGAACAAATCC
CTAATGATCCTGTCTGGCTTGCCTGCCAGATTTTATTGATTACCCTGAGCGGCAGGAATGCAACTGC
AGGCCTCAAGAATCCCTTATGTTTCTGGGATGAAAACCTGTCATTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

**Plasmid Map:**


**ACCN:** NM\_052817

**Insert Size:** 2118 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_052817.2](#)

**RefSeq Size:** 7255 bp

**RefSeq ORF:** 2118 bp

**Locus ID:** 11043

**UniProt ID:** [Q9UJV3](#)

**Cytogenetics:** Xq22.3

**Domains:** zf-B\_box, RING, BBC, SPRY, FN3

**Protein Families:** Druggable Genome

**MW:** 79.9 kDa

**Gene Summary:** The protein encoded by this gene is a member of the tripartite motif (TRIM) family. The TRIM motif includes three zinc-binding domains, a RING, a B-box type 1 and a B-box type 2, and a coiled-coil region. The protein localizes to microtubular structures in the cytoplasm. Alternate splicing of this gene results in two transcript variants encoding different isoforms. [provided by RefSeq, Feb 2009]

Transcript Variant: This variant (2) uses an alternate in-frame splice site in the coding region, compared to variant 1. It encodes isoform 2 which is shorter than isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.