

Product datasheet for **SC327979**

LMF2 (NM_033200) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LMF2 (NM_033200) Human Untagged Clone
Tag:	Tag Free
Symbol:	LMF2
Synonyms:	TMEM112B; TMEM153
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF:

>SC327979 representing NM_033200.

Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGCGGGCTCCCGGCTCCCGCGGCAGCTCTTCTCCAGGGCGTGGCGGCGCTTTCATGTTTGCTTTC
GCTTCCCTCTACACGCAATCCAGGCCGTGTATGGCCCCGAGGGCATCCTACCTGCAAGGAGGACGCTG
CGGCCTCAGGGCAAGGGGCGCTGGCAGCAGCTGTGGGAGACCCCGACGCTGCTGTGGGAAGCGCCGAGA
CTGGGGGTGGACACGGCCAGGGCCTGGAGCTGCTGAGCCTGCTGGGTGCACTAGTGCCCTGGGAGCC
CTGCTGCTGAGCCCACTGCGCCACCCTGTCTACTTGTGCTTTGGGCGCCTACCTGTGAGCCTGC
CAGGTGGGCCAGGTGTTCTTTATTTCCAGTGGGACTCCCTGCTGCTAGAGACTGGCTTCTGGCCGTG
CTGGTGGCCCCGCTGAGGCCAGCCTCCACCGCAAGGAGGCCCGCCAGGGCAGGCAGGCGCCCTG
CCCCACGAAGACCTCCCTTCTGGCTGGTGGATGGCTGCTGTTCCGCTCATGTTGCGCTCAGGCGTG
GTCAAGCTGACCAGCCGCTGCCCTGCGTGGTGGGGCTCACTGCCCTCACCTACCACTACGAGACCCAG
TGCTTGGCCACGCGCCGCGCTGGTTGCGACACCCTGCCGCTGGCTGCACAAGCTCAGCGTGGTG
GCCACCTTCTAATTGAGATCGCTGTGCCGCCCTGTTCTTCGCCCCATTTCAGCGCTGCGCTTGGCT
GCTTTCTACTCGCAGGTGCTGCTGCAGGTCTGATTATCATCACCGGCAACTACAACCTTCTCAACCTG
ATGACGCTGGTGCTTACCACTGCGCTGCTGGACGACCAGCACCTGGCTGCTGAGCCTGGCCACGGCAGC
CGCAAGAAGACGGCCACCTCCTGGCCCAAGGCCCTGCTGGCCACCCTGTCGCTGCTGCTGGAAGTAGCC
GTCTACGGGCTTCTGGCCTATGGCACTGTGCACTACTTTGGCCTGGAGGTTGACTGGCAGCAGCGCACC
ATCCACTCCAGAACCTTTACCTTCCACAGTTTTCTCAGTGGCTGAAGACTGACGCTGCCACT
GTGTGGCTGGGTGTGGCCTCCCTGGTCTGGGAGCTGCTGAGTGCCTGTGGAGGTGGACCCAGGTGCGG
GGCTGGCTACGGAAGCTCAGTGTGTAGTCCAAGTGTCCCTGTGGGCACTGCGACCGTGGCCTTGTTC
CTGATTAGCCTGGTGCCGTACTCCTACGTGGAGCCCGGGACCCACGGGCGCCTCTGGACCGGGGCCAC
CGCCTGTTTGGTGCCGTGGAGCACCTACAGCTGGCCAAGTCTACGGCCTCTCCGCCGATGACTGGG
CTTGGTGGACGGCCTGAGGTGGTGTGAGGGGAGTTACGACGGCCACCCTGGACGGAGATCGAGTTC
ATGTACAAGCCTGGGAACCTGAGCCGGCCGCCCGGTTGTGGTGCCCCACCAGCCACGCTGGACTGG
CAGATGTGGTTTGACGCCCTGGGCCCACACAGCAGCCCGTGGTTACAAGCCTGGTCTTGCGCCTG
CTGCAGGGCAAGGAGCCAGTGATCCGCTTGTCCAGAGCCAAGTGGCCAGGTATCCCTTCCACAAGCAG
CCGCCCCACCTACGTCCGAGCCAGCGCTACAAGTACTGGTTCTCCAGCCTGGGGAGCAGGGCCAGTGG
TGGCGGCGCCAGTGGGTGGAGGATTCTCCCATCCGTGTCCCTGGGGGACCCACGCTGGAGACGCTG
CTCAGGCAAGTTTGACTACAGGAGAAAAGCCACCTCGCACCCGACGCGCAACAGCACCTTGGCCAG
GCCCTCCACTGGACTCGCTCTCAGCTGTCTCCCTGGAGGCCCGCCCTGCTCTGGGGGCTCCTCATG
GCCGTGGGGGCTGTGAGATTTGTGCAAGCCCTGTAGCACCTGTTCTCTCCGGTCTCCCGCTGGCA
CCAGTCAGCGGGGAGAAGCGCAGGCCAGCCTCCGAAAGACTCCGGAGCTGCCTCCGAACAGGCCACC
GCAGCCCCCAACCCCTGCTCCAGTAGTTCGAGGACCAACCGGCGAAAGAAGTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites:

Sgfl-MluI

ACCN:

NM_033200

Insert Size:

2124 bp

OTI Disclaimer:

Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_033200.2
RefSeq Size:	2617 bp
RefSeq ORF:	2124 bp
Locus ID:	91289
UniProt ID:	Q9BU23
Cytogenetics:	22q13.33
Protein Families:	Transmembrane
MW:	79.7 kDa
Gene Summary:	Involved in the maturation of specific proteins in the endoplasmic reticulum. May be required for maturation and transport of active lipoprotein lipase (LPL) through the secretory pathway (By similarity).[UniProtKB/Swiss-Prot Function]