

Product datasheet for **SC327981**

RTF1 (NM_015138) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RTF1 (NM_015138) Human Untagged Clone
Tag:	Tag Free
Symbol:	RTF1
Synonyms:	GTL7; KIAA0252
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC327981 representing NM_015138.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGCGCGGTCCGCTTTGTGTGGTTCGAGCAGCGCGCGCGCGGCAGTGGCGGTCCACTGGCAGGC
GGCAAGAGGGGAGTCCGGCGCGCGCGCGCGTGGGAGCCGGGGACCACCATGGTAAAGAAGCGGAAA
GGCCGCGTCGTGATCGACTCGGACACAGAGGACAGCGGCAGCAGGAGAACCTGGATCAGGAGCTTTG
TCCTGGCAAAGCGAAAGCGCAGTACTCTGAGGAGAAGGAGCCGCTGTGAGTACGCTGCAGCCTCG
TCAGACTCGGAGACGTCTGACAGTGACGATGAGTGGACATTTGGGAGCAATAAAAAAAGAAGAAAGGA
AAAGCCAGAAAAATAGAGAAGAAAGGAACCATGAAGAAACAGGCCAACAAAACCTGCCTCCTCAGGCAGT
TCAGACAAAGACAGTTCAGCTGAGAGCTCAGCCCTGAGGAAGGTGAAGTGTGAGACTCTGACAGCAAC
AGCTCCTCTTCCAGTTCAGATTCAGACTCTTCTCAGAAGATGAAGAGTCCATGATGGCTATGGAGAA
GACCTCATGGGAGATGAGGAAGACAGGGCCGCTGGAACAGATGACAGAGAAAGAGAGAGCAAGAA
CTGTTCAATCGCATAGAGAAGAGGGAGGTGTTGAAAAGAAGATTTGAAATCAAGAAAAAACTAAAAACA
GCCAAAAAGAAAAGAAAAGAAAAGAAAAGAAAAGCAAGAAAGAGGAGCAAGAAAAGAAAAAACTGACA
CAGATTCAGAATCTCAGGTAACATCCCAACAACAGGAACGGCGTTCCAAGCGGGATGAGAACTAGAC
AAGAAATCTCAAGCCATGGAGGAGCTAAAAGCAGAGCGAGAAAAACGAAAGAACGAAACAGCTGAGCTC
TTGCCAAAAACAGCCATTA AAAACCAGTGAGGTCTACTCTGATGATGAAGAGGAGGAAGAGGATGAC
AAATCCAGTGAAGAGTACAGCCGCTCATCACGAACATCATCGTCTGATGAAGAAGAGGAGAAAGAAGAG
ATCCCTCCCAATCCCAACCAGTTTCTTACCTGAAGAATTGAATCGGGTTCGATTATCACGGCATAAG
CTAGAAGCTGGTGTACATGCCCTTCTTTGCTAAAACGTGACAGGATGTTTTGTGCGGATTGGCATC
GGAAACCAACAGCAAACCAGTTTACCGGGTCTGAGATTACGGGTGTTGTGGAACCTGCCAAAGTT
TACCAACTAGGTGGCACCAGAACAACAAGGGCTGCAACTACGGCATGGCAATGACCAACGCGTGTTC
CGTTTAGAGTTTGTCTCAAACCAAGAATTCACCGAAAAGTGAAGTTATGAAGTGGAAAAGAAGCGATGTTT
TCTGCTGGCATGCAGTTGCCACTCTAGATGAAATCAATAAAAAGGAATTATCTATTAAGAAGCTCTT
AATTATAAATTAATGATCAGGACATTGAAGAGATTGTAAGAGAGAAAGAAAGGTTTCAAGAAAGCTCCA
CCCAACTACGCTATGAAGAAGACTCAGCTACTGAAGGAAAAGGCCATGGCTGAGGACCTGGGGGATCAG
GACAAGGCCAAACAAATCCAAGATCAACTGAATGAGCTGGAGGAACGGGCAGAGGCCCTGGACCGCCAG
CGGACCAAGAACATATCCGCTATCAGTTACATCAACCAGCGGAACCGGGAGTGAACATTGTAGAGTCT
GAGAAGGCCCTTGTGGCTGAAAGTCACAACATGAAAAACCAACAGATGGATCCCTTTACTCGCGGCAG
TGCAAGCCTACCATCGTTTCTAATTCAGAGACCCAGCTGTTCAAGCTGCCATCTTGCCCAAGCTGAAT
GCAAAATACGGTTCTGGAGTGTACCAGATGCTCAAAGGAAATGAGCAAGGGTCAAGGCAAGATAAA
GATTTGAATCTAAGTCAGCCAGTGACCTCTCAGAAGATCTGTTCAAAGTACACGATTTTGTGTAAG
ATTGACTTACAAGTTCCAGCTCAGAGTCAAAGGCTTTAGCCATCACCTCCAAGGCTCCGCCAGCCAAG
GATGGGCTCCAAGGAGATCTCTGAACTTGAAGACTACAAAAACGACGAGGGCTTATTTGA
AGCGGACCGACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGAT
ATCCTGGATTACAAGGATGACGACGATAAGGTTTAA
  
```

Restriction Sites: Sgfl-RsrII

Plasmid Map:


ACCN: NM_015138

Insert Size: 2133 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_015138.4](#)

RefSeq Size: 5038 bp

RefSeq ORF: 2133 bp

Locus ID: 23168

UniProt ID: [Q92541](#)

Cytogenetics: 15q15.1

MW: 80.3 kDa

Gene Summary: This locus may represent a gene involved in regulation of transcription elongation and chromatin remodeling, based on studies of similar proteins in other organisms. The encoded protein may bind single-stranded DNA. [provided by RefSeq, Sep 2010]