

Product datasheet for **SC327984**

TMPRSS7 (NM_001042575) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TMPRSS7 (NM_001042575) Human Untagged Clone
Tag:	Tag Free
Symbol:	TMPRSS7
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC327984 representing NM_001042575.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTTTCGCATCACCACATTGAGTTTCTCCGAATACCGACAAAAGGAGTCCAGGGAATTTCTTTCA
GTGTCACGGACTGTGCAGCAAGTGATAAACCTGGTTTATACAACATCTGCCTTCTCCAAATTTTATGAG
CAGTCTGTTGTTGCAGATGTGCAGCAACAACAAAGCGGCCCTCTTGTCACATTTTGGATTGTTTTGTC
ATGCCACGTGCCAAAGGCCACATCTTCTGTGAAGACTGTGTTGCCGCATCTGAAGGACTCCATCCAG
ACAAGCATCATAAACCGGACCTCTGTGGGAGCTTGACGGGACTGGCTGTGGACATGGACTCTGTGGTA
CTAAATGACAAAGGCTGCTCTCAGTACTTCTATGCAGAGCATCTGTCTCTCCACTACCGCTGGAGATT
TCTGCAGCCTCAGGGAGGCTGATGTGCACTTCAAGCTGGTGGCCATAGTGGGCTACCTGATTCTGCTC
TCAATCAAGTCCATCCAAATCGAAGCCGACAACGTGTCACTGACTCCCTGACCATTACGACTCCCTT
TTGCCCATCCGGAGCAGCATCTGTACAGAATTTGTGAACCCACAAGAACATTAATGTCATTGTTTCT
ACAAATAATCTCATGTTGGTGACATTTAAGTCTCTCATATACGGAGGCTCTCAGGAATCCGGGCATAT
TTTGAGGTCATTCCAGAACAAAAGTGTGAAAACACAGTGTGGTCAAAGACATCACTGGCTTTGAAGGG
AAAATTTCAAGCCCATATTACCCGAGCTACTATCTCCAAATGCAAGTGTACCTGGAAATTTGAGACT
TCTCTATCAACTCTTGGCATAGCACTGAAATCTATAACTATTCAATAACCAAGAAGAGTATGAAAGGC
TGTGAGCATGGATGGTGGGAAATTAATGAGCACATGTACTGTGGCTCCTACATGGATCATCAGACAATT
TTTCGAGTGCCCGACCCCTCTGGTTCACATTGAGTCCAGTGCAGTTCAAGGCTTTGAGACAAGCCACTT
TTGGCAGAATATGGCAGTTACAACATCAGTCAACCCTGCCCTGTTGGATCTTTAGATGCTCCTCCGGT
TTATGTGTCCCTCAGGCCAGCGTTGTGATGGAGTAAATGACTGCTTTGATGAAAGTGTGAACTGTTT
TGCGTGAGCCCTCAACCTGCCTGCAATACCAGCTCCTTCAGGCAGCATGGCCCTCTCATCTGTGATGGC
TTCAGGGACTGTGAGAATGGCCGGGATGACAAAACGCACTCAAAGTATTCCATGCAACAACAGAAGT
TTTAAGTGTGGCAATGATATTTGCTTTAGGAAACAAAATGCAAAATGTGATGGGACAGTGGATTGTCCA
GATGGAAGTGTGAAGAAGGCTGCACCTGCAGCAGGAGTTCCTCCGCCCTTCACCGCATCATCGGAGGC
ACAGACACCCCTGGAGGGGGGTTGGCCGTGGCAGGTGAGCCTCCACTTTGTTGGATCTGCCTACTGTGGT
GCCTCAGTCATCTCCAGGGAGTGGCTTCTTTCTGCAGCCCACTGTTTTCATGGAACAGGCTGTCAGAT
CCCACACCATGGACTGCACACCTCGGGATGTATGTTGAGGGGAATGCCAAGTTTGCTCCCCGGTGAGA
AGAATTGTGGTCCACGAGTACTATAACAGTCAGACTTTTGATTATGATATTGCTTTGCTACAGCTCAGT
ATTGCCTGGCCTGAGACCCTGAAACAGCTCATTAGCCAATATGCATTCTCCCACTGGTCAGAGAGTT
CGCAGTGGGGAGAAGTGTGTTAACTGGCTGGGGGCGAAGACACGAAGCAGATAATAAAGGCTCCCTC
GTTCTGCAGCAAGCGGAGGTAGAGCTCATTGATCAAACGCTCTGTGTTTCCACCTACGGGATCATCACT
TCTCGGATGCTCTGTGCAGGCATAATGTGAGGCAAGAGAGATGCCTGCAAAGGAGATTCCGGTGGACCT
TTATCTTGTGCAAGAAAAAGTGTGAAAAATGGATTTTACTGGCATTGTTAGCTGGGGACATGGAAGT
GGACGACCAAACTTTCTGGTGTACACAAGGTGTCAAACCTTTGTTCCCTGGATTATAAATATGTC
CCTTCTCTTTGTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001042575

Insert Size: 2154 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001042575.2</u>
RefSeq Size:	2509 bp
RefSeq ORF:	2154 bp
Locus ID:	344805
UniProt ID:	<u>Q7RTY8</u>
Cytogenetics:	3q13.2
Protein Families:	Druggable Genome
MW:	80.4 kDa
Gene Summary:	Serine protease which preferentially hydrolyzes peptides with Arg at the P1 position. [UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) encodes a 717aa protein.