

Product datasheet for **SC327992**

ATHL1 (PGGHG) (NM_025092) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ATHL1 (PGGHG) (NM_025092) Human Untagged Clone
Tag:	Tag Free
Symbol:	PGGHG
Synonyms:	ATHL1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC327992 representing NM_025092.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGGAGGACGCCGGCGAGGACCCACCAGTTTGCTGCCACTCTCTGCCAGTGACCCCGTCTCTTG
GCCACTGTGACCAACGCATACCTGGGCACACGAGTGTTCACGACACGCTGCACGTGAGCGGCGTGTAC
AATGGGGCTGGCGGGACACGCACCGGGCATGCTGCCAGCCCTCAACGTCCGGCTGGAGGCCCT
GCAGGGATGGGGAGCAGCTGACCGAGACCTTCGCCCTGGACACCAACACAGGCTCCTTTCTCACACC
CTGGAGGGCCCCCGCTTCGGGGCTCCAGTGCATCTATGCGCATCGCACGCTGCCCCACGTGCTGGCT
TTCCGAGTGTCCATCGCCCGCTGGCCCCGGGAGCGGGCCATCACGCTGCTCTGCGGTGAGCTTC
TCCCCAGAAAGCCAGACCTGGACCTGCATCAGGGTCTGACTTCAGGGAGCCCGGTACCTGTATGGC
CACACCCTCACCCCTGAGCAGCCCGGGGGCCACAGCAAGAGGTACACATGCTGTGGACACCAGCACCC
CCAGACCTGACCCCTGGGAAGGTGAGGAGCTAGGACGTGGACTTCTGACAGCAGTGGCGGCAGC
CAGGCTGAGGCTCAGGCCTGCCTACTGAGGCCCTGCAGCTGCAGGCCAGGGAGCTCTGTATACGGCT
CAGGCACAGGCCCTGGGCCAGCTCTGGGTAGAATGTGGCTTGGACGTGGTGGGGCCCTGCAGCTGCCG
CAGGCCCTGCGTGGCTCCCTCTACTACCTGCTCAGTGCCCTGCCCCAGCCAAAGGCCCCAGGATACATC
TGCCATGGCCTCAGTCTGGGGCTCTCCAATGGGAGCCGTGAGGAATGCTACTGGGGCCACGTCTTC
TGGGACCAGGACCTCTGGATGTTCCCGAGTATCCTGATGTTCCACCCAGAAGCCGCCAGGGCCATCCTG
GAGTACCGCATCCGCACGCTGGACGGGGCCCTGGAGAAGCCCAAGCTGGGCTACCAGGGAGCCAAAG
TTTGCCTGGGAGAGTGCAGACTCCGGCTTAGAGGTTGCCCTGAGGACATTTACGGAGTCCAGGAGGTG
CAGCTCAACGGGGCCGTGGTGTGGCTTCGAGCTGACTACCATACCACCCAGGACCTGCAGCTATTT
CGAGAGGCTGGTGGCTGGGACGTGGTCAAGGCTGTGGCCGAGTTTGGTGCAGTGTGTTGAGTGGAC
CCAGGGAGGAAAAGTACCACCTGAGGGGAGTCAATGTCACCCGACGAGTACCATTAGGGGTCAACAA
TCTGTGTACACCAACGTCTGGTCCAGAACAGCCTGCGCTTTGCTGCTGCCCTGGCCAGGACCTGGGT
CTTCCCATCCCCAGCCAGTGGCTGGCGGTGGTGACAAGATCAAGGTACCCTTTGACGTGGAGCAGAAC
TTCCACCCGGAGTTCGATGGGTATGAGCCTGGAGAGGTGGTGAAGCAGGCAGACGTCGTGCTCCTGGGA
TACCCAGTCCCCTTCTCCTGAGTCTGATGTTGCGAGGAAAAATCTGGAGATTTACGAGGCTGTGACG
TCCCCCAGGGCCCCGCCATGACCTGGAGCATGTTTGTGTGGGCTGGATGGAGCTGAAGGACGCAGTG
CGGGCCCCGGGCCTCCTGGACAGGAGCTTTGCCAATGGCTGAACCCTCAAGGTGTGGACGGAGAAT
GCAGACGGGTGAGCGCTGTGAACCTCCTGACAGGCATGGGGGGCTTCTGCAGGCGGTGGTCTTCGGG
TGCACGGGGTTCAGGGTACCCGAGCGGTGTGACCTTTGACCCTGTGTCTGTGCGGGATCTCCAGA
GTGAGCGTCTCCGCATCTTCTACCAGGGGAACAAGCTCAACTTCTTTTTCCGAGGACTCCGTGACC
GTGGAGGTACAGCTCGAGCAGGGCCCTGGGCTCCTCACCTGGAGGCTGAGCTGTGGCCATCCCAGTCC
CGGCTCTCCCTGTTGCCAGGACACAAGTCTCCTTTCCCGCTCGGCTGGCCGGATACAAATGTCACCC
CCGAAGCTGCCTGGAAGTTCAGCTCCGAGTCCCTGGGAGGACTTTTTCAGATGTTAGGGACCCGCTC
CAGAGCCCTCTGGGTACCCCTGGGTCTCCAGCCCCACCGAGTCACTCACTGTGGACCTGCCTCT
GAATAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_025092

Insert Size: 2214 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_025092.4](#)

RefSeq Size: 3270 bp

RefSeq ORF: 2214 bp

Locus ID: 80162

UniProt ID: [Q32M88](#)

Cytogenetics:	11p15.5
Domains:	Glyco_hydro_65m
MW:	80.7 kDa
Gene Summary:	Catalyzes the hydrolysis of glucose from the disaccharide unit linked to hydroxylysine residues of collagen and collagen-like proteins.[UniProtKB/Swiss-Prot Function]