

Product datasheet for **SC327998**

CCDC157 (NM_001017437) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CCDC157 (NM_001017437) Human Untagged Clone
Tag:	Tag Free
Symbol:	CCDC157
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC327998 representing NM_001017437.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGATCGCC
ATGGCGCACCTGCTGGGCAGCCAGGCCTGCATGGAGAGCCTGCGCACAGACCTCACCGACCTGCAGGGT
GCCATCGTAGACGTCTTCTCCCGCGCCGGGCCTGTGCGCTTCCCCTCCTGGAAGTTCCCTGACCGCATG
GCCTGTGACCTCGACATGGTGGCCCTGCTGGAGCACTATGACCATGTTCCGGGTGACCCCGAGTTCACG
CAGCTGTCCCATGCCGTGCTGCTGGAGCTGGTCACTCGACAGGCTCCTGCTGCTGCTTCAAAGCTGTATG
AGCTACTTGGAGAACCTTGGCTCAGAGCAGATGATGCCCCCTGCACAGGCTGCGGGGCCCTGCATGTCC
GTGGGGCTCACGGTGGCGCTTCTGGGACAGCCTGCTGAGGCTGGGCACGCTCCACCAGCAGCCACTC
CCCCAGAAAGGGCAAACCAAAGGGAGACTCCACCTCCAAGCCACCACCAAGGGCGAGCCAGCCAGG
AGCCCTGAATATCTGACTACCAAGTTAATCAAGCCCTCCTCCCACTGCTAGGCTTGCCCCAGACCTGC
CAAGAGCCAGAGAGCATCCCTGTGAGAGCCTCCCTGCAGTTCACAGCCACGACCTTCAAGAACACCAGG
AGTGTCCACTCCAGACATTGAGACGGCCCTGGTGGCCTGTGACGCATGCGCCAGCGTCCAGGGAAGC
CTGCAGAAAGTGGCAAGGTGGTCATCAGCCTGTGTGAGAGCCAGAACCTGCCCTCGTCTTAGGCCAG
TTCCAGCAACTGGTGCAGGACAGCATGGGGCTCAGGCCACTGCCGGCTGCCACCGTGGGCGCTGGGCA
GCAGAGCAGAGGAAAGACCTGACGCGCCTCAGTAAGCATGTGGAGGCCCTCAGGGCCAGCTGGAGGAG
GCTGAAGGGCAGAAGGATGGCCTGAGGAAGCAGGCGGGCAAGCTGGAGCAGGCGCTGAAACAGGAGCAG
GGGGCACGGCGGGCAGAGGCGGAGGAGGATGAGCAGTGCCTGTCTGAGTGGGAGCACGACAAACAGCAG
CTGCTCACAGAAACAAGTGACCTAAAGACAAAGATGGCCACCCTGGAGAGAGAACTGAAACAGCAGCGG
GAGTCCACACAGGCTGTGGAGGCAAAGGCCAGCAGCTGCAGGAGGAAGGTGAGCGCAGGGCGGCAGCG
GAGAGGAGGTGCAGCAGCTGGAGGAGCAGGTGCAGCAGTTGGAGGCGCAGGTGCAGCTGTTGGTGGGT
CGGCTGGAGGGCGCTGGCCAGCAGGTCTGCTGGGCCAGCACGGAGCTGGATAAGGAGAAGGCCCGTGTG
GACAGCATGGTCCGCCACCAGGAGTCTCTGAGGCCAAGCAGCAGGCCCTGCTAAAGCAGCTGGACAGC
CTGGACCAGGAACGTGAGGAGCTGCGGGGAGCCTGGACGAGGCTGAGGCCAGCGGGCCCGCTGGAG
GAGCAGCTGCAGAGCGAGCGGGAGCAGGGCAATGCCAGCTCAGGGCCAGCAGGAGCTGCTGCAGAGC
CTGCAGAGGGAGAAGCAAGGCCTGGAGCAGGCGACTACGGACCTGCGGCTAACCATCCTGGAGCTAGAA
CGGGAGCTGGAGGAGCTGAAGGAGCGGGAGCGGCTGCTGGTGGCCTTCCAGACCTGCACAGGCCACC
GAGACCCAGATCCATGGAGGCAGATCCAGCAGCGTGAATCCAGATAACATGTCCACAGACTCTGGC
AACGTCACAGATCACATGGAGAGCAAGTGCAGGCCAACGACATCCGCATCCGGGTCTACAGGAGGAG
AACGGGCGGCTCCAATCAATGCTGTCCAAAATCCGGGAAGTGGCCAGCAGGGTGGCCTCAAGCTGATC
CCGAGGACCGGCTCTGGTCCCCTTCCAGCAAGGGAACCCAGGGAGCAACACCACAGTCCAGGCCAAG
AGCACATCCCCAGGGCCTCTGGGAGACAGCACCTTCTAGCAGCAGGACGGGTAGGACCTGCTGGGC
CAGCCCTGCACATCCCCACCTCGGCAGCCCTGCACATCCCCACCTCGGCAGCCCTGCACATCCCCACCT
CGGCAGCCCTGCACATCCCCATCTCGGCAGCCCTGCAGCCAGCCAGCAAGTCTTGCTGGAGGGTGTG
ACCCACTTGGACACCTGCACCCAGAACCCCATCAAGGTCTTGGTCAGGCTGAGAAAGAGACTGTACCT
GGCCGGGACAGGCCAGCTCTGCACACCAGCCCCAGGAGCGGCCATGTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI
ACCN: NM_001017437
Insert Size: 2259 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001017437.3](#)

RefSeq Size: 4832 bp

RefSeq ORF: 2259 bp

Locus ID: 550631

UniProt ID: [Q569K6](#)

Cytogenetics: 22q12.2

MW: 83.9 kDa