

## Product datasheet for **SC328001**

### Zinc finger protein 287 (ZNF287) (NM\_020653) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Zinc finger protein 287 (ZNF287) (NM_020653) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF287
Synonyms:	ZKSCAN13; ZSCAN45
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

**Fully Sequenced ORF:** >SC328001 representing NM\_020653.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGTTAGCCTCAAGCAAGAGGATGAACAGTTCTCACGTTCTCAAATCCTTCTAAGTGGAAGTCAGAC
AAGGCTCAGAGTGGACCCTACAATGTTGAGAAGGAAATCCTTACTTCAAGATTCTTGCCTGACACTGAG
ACCTGTGCACAGAATTTTAGGAATTTTCCATACCCAGACCTGGCTGGTCTCGAAAGGCATTGAGTCAA
CTCCGAGAGCTCTGCCTTAAGTGCTGAGACCTGAGATTCACTCAAAGGAACAAATTTTGGAGCTGCTG
GTGCTGGAGCAATTCCTGACCATCTGCCTGGTGGAGTTAGGACTTGGGTAAGTCCCAGTATCCAGAG
AGCAGCGAGGAAGCAGTACTCTGGTGGAGATTGACTCAGATTCTAGAAGAGGAAGCTCCTCAAAAC
TCTACCTTTCCCAAGATACCCAGAGGAAGACCCAGAGGAAAACATGCTTTCCAGACAGGATGGCTA
AATGACTTGGTGACCAAGAATCGATGACATTCAAAGATGTGGCTGTAGACATCACCCAGGAGGACTGG
GAGTTAATGCGTCTGTGCAGAAGGAATTATACAAGACTGTGACGTTACAGAACTATTGGAACATGGTT
TCTCTGGGACTTACAGTGTACAGACCAACTGTGATTCCCATATTGGAAGAACCATGGATGGTGATAAAA
GAAATTTTAGAAGGCCCTAGTCCAGAATGGGAAACTAAAGCCCAAGCATGTACTCCAGTGGAGGATATG
TCTAAACTCACAAGGAAGAAACCCATACCATCAAATTAGAAGACTCATATGACTACGATGATAGACTA
GAGAGGCGAGGAAAAGGTGGCTTCTGGAAAATTCACACTGATGAAAGAGGTTTCAGTTTGAAGTCAGTC
CTTTCACAAGAATATGATCCTACAGAAGAATGTCTTAGTAAATATGATATATAGAAAATATTTTGAA
AAGCATCAAACCTAATGTACAGTTTGATACCAATTAGATAATAAACTTCTGTGTATAATGAAGGC
AGGGCAACCTTCAATCATGTCTCATATGGTATTGTACATAGGAAAATACTTCTGGAGAGAAGCCTTAC
AAGTGTAATGTGTGTGGGAAAAAATTTAGGAAATACCCATCCCTCCTGAAACACCAAAGTACCCATGCC
AAAGAAATCGTATGAATGTGAAGAATGTGGGAAAGAGTTTAGGCATATCTCATCCCTTATTGACAT
CAGAGAATGCACACTGGAGAAAAACCATATGAATGCCACCAGTGTGGTAAAGCCTTCAGCCAGCGTGCA
CACCTTACTATACATCAGAGAATTCATACTGGAGAGAAACCTATAAGTGTGATGACTGTGGGAAAGAC
TTCAGTCAGCGTGACACCTTACCATCCATCAAAGGACACATACTGGAGAGAAACCATATAAATGCTTG
GAATGTGGTAAAACCTTCAGTCATAGTTCATCACTGATTAATCATCAGAGAGTTCATACTGGAGAAAAA
CCTTATATATGCAATGAATGTGGGAAAGACTTTCAGTCAGAGTACACACCTTCTTCAGCATCAAAAAATA
CATACTGGGAGAAACCATATAAATGCAATGAATGTTGGAAAGTGTTAGTCAGAGTACTTACCTTATT
CGACATCAGAGAATTCATTCTGGAGAGAAGTGTATAAATGTAATGAATGTGGAAAAGCCTTTGCTCAT
TCCTCAACCTTATTCAACATCAAACCACTCACACTGGAGAGAAATCCTATATATGTAATATATGTTGG
AAAGCCTTCAGCCAGAGTGCAAACTTACTCAACATCATAGAACACATACTGGAGAGAAACCATATAAA
TGCAAGTGTGTGGGAAAGCATTAGCCAGAGTGTGCACCTTACTCAACATCAGAGGATTCATAATGGA
GAAAAACCTTTAAATGCAATATATGTGGGAAAGCATATAGACAAGGCGCAAATCTTACTCAGCATCAA
AGGATTCATACTGGAGAAAAACCTATAAATGTAATGAATGTGGGAAAGCTTTCATTTATTCCTCATCA
CTTAATCAACATCAGAGAATCATACTGGAGAGAGACCCTATAAATGTAATGAATGTGATAAGGATTTT
AGCCAGAGAACATGCCTTATTCAACACCAGAGAATTCACACAGGAGAGAAACCTATGCATGTCGTATA
TGTGGTAAAACCTTCACCCAGAGTACAAACCTTATTGAGCATCAACGTGTTTATACAGGTGCCAAACAT
CGTAATTA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

**Plasmid Map:**


**ACCN:** NM\_020653

**Insert Size:** 2286 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_020653.2](#)

**RefSeq Size:** 4278 bp

**RefSeq ORF:** 2286 bp

**Locus ID:** 57336

**UniProt ID:** [Q9HBT7](#)

**Cytogenetics:** 17p11.2

**Domains:** KRAB, LER, zf-C2H2  
**Protein Families:** Transcription Factors  
**MW:** 88.3 kDa

**Gene Summary:** This gene encodes a member of the krueppel family of zinc finger proteins, suggesting a role as a transcription factor. Its specific function has not been determined. This gene is located near the Smith-Magenis syndrome region on chromosome 17. [provided by RefSeq, Jul 2008]  
Transcript Variant: This variant (2) differs in the 5' UTR compared to variant 1. Variants 1-4 all encode the same isoform (a). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.