

Product datasheet for **SC328008**

CDHR4 (NM_001007540) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CDHR4 (NM_001007540) Human Untagged Clone
Tag:	Tag Free
Symbol:	CDHR4
Synonyms:	CDH29; PRO34300
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC328008 representing NM_001007540.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGTGCTGCTCAGGCTCCTCGTGTCTCTTTGCTCCGGTGGTCTCTGACCTCTGCAGCCTGCCCTGC
TTTATAAATGTCTCTGAGAGCCAGGGCCCTGGCAGTCCTTCAGTTTTATCCTTCAACTGCTCCTCC
TACACGCCACACCCACCCTGGAGTTGCTCAATGTCCAGCCACCACCACCTTCTTCAACCACCCAGC
TTGGCCAGGTGGCAAGGGACCTATGTGGCAAGTTGACCTTGAGCAGCTCTGCTCAGTTGGATGCCCTG
ATGGTGAACCACTACAAGGTGCAGCTGAAGTTCACATGTGGCAACCATGTGATGGAGGGCTCACTCTCT
GTGGATGTGCAGCGGGACCTTAGCCATATCCAGTGTGCTGGTCAATTTGCCAGCCAGCTGGGAAATG
ATTCAGGTGCCAGAGACAGTCACACCTGGGGCTCGGCTGTACACTCTGCTCCTCCAGGCCAGAACTC
CACGGAGCGCAGATGAGCATTATCAGTGGCCAGGACCTGCCACACTTCCCTGGACCTTTCTCCATCAAT
GAGCAAGGTTGGCTGCAGGCACCATCCAGGGCCTCTAGGCCAGGCTCAAAAGGTCTTCCAGCTGCAA
ATCTCAGTGTCTTTGGACAAAGGCAAGCTGCCAAGGGATGGTGATAGTGAAGTTTTGCTGTCTCC
TCCAGCCAGGTCTCCTTCTCGAGCAGGCTCAGAATATACCATCCCTGAGAATCTGGCCCCCGGTAGT
GAGGTGGTTGAGTCCAGGCCCGGGGTGTGACCTGCGCTATGAAATCCTGTCTCCGGTGGCCAGCCCA
CTCTTCTCCATTGGTCTGTCAGACGGTGTGGTCCGGACCACCACGCCCTAGAGTTAGCTCGCACCTCA
GGCACCAGCGGTCTCCAGGCTGCAGGTGAAGGCCTTTGAGCAGGGCCAGCTGTGGGCCAGTGCCAAGCTC
AATCTCACCATGAATGTGCAGCTGGTCAACCTCTGGCCTCCACGCTGCCTCCAGCGCTTCTGGTGTCC
CAAATCCCGGAGACTGCACCCGTGGGCACCGTGTGAATACTCTCACTTGCGAAGATCCGGACTCTGTT
GGTGCCACCCTGGACTACAAGCTGTGGTCCGCAGCTCTTCAACCTGCCAGCCTCTGCCTTTATGAC
AGAGTCTTGAGGTGAATGCCACACTGGACTGTGACACTCCTGGAGCCTGCTTCCAGCATGCAGCCTCC
ATCCTGGTGTCTGATGGTGGCCAGCCCCAGATGACCACTGAGGTGCCGGTACTGGTGATGGTGACACCC
ATCAACGAGTTCTCCCCAGCCTGTGCCCTCGCACGTTCCGGGTTGAGGAGGATGCGGCGCCCCACACT
CTACTGGGCTCCGTGGTGGGCACGGATATGGATTACCCTCATGACAACATTGAGTACTACACCTCTGGT
GGTCTACACCTTTGCTGTGGACCGTCTCAGCGGGGAAGTTACCTCCTGGGACCTTTGGACTATGAG
CAGCAGAGGCTGTACAGGCTCACTGTCTTGTGATTGACCATGGCCAAGACCAGAACCCCAACCATCAC
CTCTCAGGCTCCTGTACCATTACCATCGAGGTTGAGGATGTGAATGACCATGCCCCGAGTGTGAGCCC
CCATTTAGGAACCTACCATCTATGCTCCTCTGGGCCGTAGCGTGGAGGTGACCAAGATGTCATGCCAG
ATCCCTCAGGAGCCACAGCGCCTGATCTACTCCTATAGCATCGTGGGAGGGAATAGCCAGAACCATTTC
ATCCTGCAAGGGGCCATCCTGGTGCACAGTGACCTTGTGTGGGGCCGTTCTGGCCAGAGCAGCCCCGT
ACCTATGAGCTACTGATCTGTGTGGCCGATGCAGGCCCTCCACCCCCACCTCAGCACCACAGCCACC
ATTATTGTGATCTAGTTCCCGGAGGGCCAGCACAGTGGCCACCAGCACCACAGAACCACAGTGCCC
TCAACGATGACACCCATGCTCGTGACAGACACAGAGGCTTTCTGGCAGCCACAGCCCTGGTTTGGTG
GTGTTGACAGCAACTGGTGTCTTCTCCTCTTGGCCCTAGGCTGGCTTCTTGGCAGGCTCCTCCAGGGG
TTGGCCAGCTGCTGCAAGCACCAGCAAAACAGCCAGGCTTTGCTGCTAAACAGCATCCAGGGAAT
GAGGGATCCATCGAGGGTTTCTGGAGGCACCGAAGATGGAGATGTCCAGGCACCCAGCAGTGTATG
AGTCTGCATTTTATGAGGAGCAGCAGGACTCCCGTACCGGAAGAGACTACCTGTTTAAACACACACACA
GGAGCCCGGCGCTGGCTCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI
ACCN: NM_001007540
Insert Size: 2367 bp

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001007540.3</u>
RefSeq Size:	2514 bp
RefSeq ORF:	2367 bp
Locus ID:	389118
UniProt ID:	<u>A6H8M9</u>
Cytogenetics:	3p21.31
MW:	85.8 kDa
Gene Summary:	Cadherins are calcium-dependent cell adhesion proteins. They preferentially interact with themselves in a homophilic manner in connecting cells; cadherins may thus contribute to the sorting of heterogeneous cell types (By similarity).[UniProtKB/Swiss-Prot Function]