

Product datasheet for **SC328026**

USP20 (NM_001008563) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	USP20 (NM_001008563) Human Untagged Clone
Tag:	Tag Free
Symbol:	USP20
Synonyms:	hVDU2; LSFR3A; VDU2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC328026 representing NM_001008563.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGGGGACTCCAGGGACCTTTGCCCTCACCTTGACTCCATAGGAGAGGTGACCAAGAGGACTTGCTG
CTCAAATCTAAGGGAACCTGTCACTCGTGTGGGTCAACGGACCAACCTATGGGCTGTCTGCAGTT
GCCTGCCCTATGTTGGCTGCGGAGAATCCTTTGCTGACCACAGACCATTTCATGCACAGGCAAAAAG
CACAACCTGACCGTGAACCTGACCACGTTCCGACTGTGGTGTACGCCTGTGAGAAGGAGGTATTCTG
GAGCAGCGGTGGCAGCCCTCTGCTGGGCTCCTCTTCCAAGTTCTCTGAACAGGACTCCCCGCCACCC
TCCCACCTCTGAAAGCTGTTCTATTGCTGTGGCTGATGAAGGAGAGTCTGAGTCAGAGGACGATGAC
CTGAAACCTCGAGGCTCACGGGCATGAAGAACCTCGGGAACCTCTGCTACATGAACGCTGCCCTGCAG
GCCCTGTCCAATTGCCCGCCGCTGACTCAGTTCTTCTGGAGTGTGGCGGCTGGTGCACAGATAAG
AAGCCAGCCCTGTGCAAGAGCTACCAGAAGCTGGTCTCTGAGGTCTGGCATAAGAAACGGCAAGCTAC
GTGGTCCCACCACTGTCTCATGGGATCAAGTTGGTCAACCAATGTTCCGAGGCTATGCCAGCAG
GACACCCAAGAGTTCTTCGCTGCTGATGGACCAGCTGCACGAGGAGCTCAAGGAGCCGGTGGTGGCC
ACGGTGGCGCTGACGGAGGCTCGGGACTCAGATTCGAGTGACACGGATGAGAAACGGGAGGGTGACCGG
AGCCCATCAGAAGATGAGTTCTTGTCTGTGACTCGAGCAGTGACCGGGGTGAGGGTGACGGGCAGGGG
CGTGGCGGGGACGCTCGCAGGCCGAGACGGAGCTGCTGATCCCAGATGAGGCGGGCCGAGCCATCTCT
GAGAAGGAGCGGATGAAGGACCGCAAGTTCTCTGGGGCCAGCAGCGTACAAACTCGGAGCAAGTGGAC
GAGGACGCTGATGTGGACTGCCATGGCTGCCCTTACGACACAGCCCGGAGGGCCAGCCCCGTCAC
CCACGGTCTCCAGCCCCTGCCGACGCCAGAGCCGACAATGATGCTCACCTACGCAGCTCCTCTCGC
CCCTGCAGCCCCGTCCACCACCAGAGGGCCATGCCAAGCTGTCTAGCAGCCCCCTCGTGAAGCCCC
GTGAGGATGGCACCGTCTGACGTCTCAAGAAAGCCAGGATTGAGTGTGGCAGCCGGAGGCGGAAG
GAGCAGCGCTACCGCAGCGTCTCTCAGACATCTTTGACGGCTCCATTCTCAGCCTGTGCAAGTGTCTC
ACCTGTGACCGGGTATCCACCACAGTGGAAACGTTCCAGGACTTATCACTGCCATTCTGGAAAGGAG
GACCTGGCCAAGCTCCATTCAGCCATCTACCAGAATGTGCCGGCCAAGCCAGGCGCTGTGGGGACAGC
TATGCCGCCAGGGCTGGCTGGCCTTATTGTGGAGTACATCCGACGGTTTGTGGTATCTGTACCCCC
AGCTGGTTTTGGGGCCTGTGTCACCCTGGAAGACTGCCTTGTGCCTTCTTGGCGTGTGAGTTA
AAGGGTGACAACATGTACAGCTGTGAGCGGTGTAAGAAGCTGCGGAACGGAGTGAAGTACTGAAAGTC
CTGCGGTTGCCGAGATCCTGTGCATTCACCTAAAGCGCTTTCGGCAGGAGTGTACTCATTCAAG
ATCAACAGCCACGCTCCTTCCCCTCGAGGGGCTCGACCTGCGCCCTTCTTGGCAAGGAGTGCACA
TCCCAGATCACCACTACGACCTCCTCTCGGTATCTGCCACCACGGCACGGCAGGACAGTGGGCACTAC
ATCGCCTACTGCCAGAACGTGATCAATGGGCAGTGGTACGAGTTTGTGACCAAGTACGTCACAGAAGTC
CACGAGACGGTGGTGCAGAACGCCGAGGGCTACGTAATCTTACAGGAAGAGCAGCAGGAGGCCATG
CGGGAGCGACAGCAGGTGGTGTCCCTGGCCGCATGCGGGAGCCAGCCTGTGCGGTTCTACGTGTCC
CGCGAGTGGTCAACAAGTTCAACACCTTCGCGGAGCCAGGCCCCATCACCAACCAGACCTTCTCTGC
TCCCACGGAGGCATCCC GCCCAAAAATACCACTACATCGACGACCTGGTGGTCTCTGCCCCAGAAC
GTCTGGGAGCACCTGTACAACAGATTTCGGGGTGGCCCCGCCGTGAACCACTGTACGTGTCTCCATC
TGCCAGGTGGAGATCGAGGCACTGGCCAAGCGCAGGAGGATCGAGATCGACACCTTCATCAAGTTGAAC
AAGGCCTTCCAGGCCGAGGAGTCCCGGGCGTCTACTGTCATCAGCATGCAAGTGGTCCGGGAGTGG
GAGGCGTTCTGTAAGGGGAAGGACAACGAGCCCCCGGGCCATTGACAACAGCAGGATTGCACAGGTC
AAAGGAAGCGCCATGTCCAGCTGAAGCAGGGAGCTGACTACGGGAGATTTCCGAGGAGACCTGGACC
TACCTGAACAGCCTGTATGGAGTGGCCCCGAGATTGCCATCCGCCAGAGTGTGGCGCAGCCGCTGGGC
CCAGAGAACCTGCACGGGAGCAGAAGATCGAAGCCGAGACGCGGGCCGTGA
ACGGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001008563

Insert Size: 2745 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001008563.4](#)

RefSeq Size: 4459 bp

RefSeq ORF: 2745 bp

Locus ID: 10868

UniProt ID: [Q9Y2K6](#)

Cytogenetics: 9q34.11

Protein Families: Druggable Genome, Protease

MW: 102 kDa

Gene Summary: This gene encodes a ubiquitin specific processing protease that was first identified as a substrate of the VHL (von Hippel-Lindau disease) protein E3 ubiquitin ligase complex. In addition to being ubiquitinated by the VHL-E3 ligase complex, this enzyme deubiquitinates hypoxia-inducible factor (HIF)-1 alpha and thereby causes increased expression of HIF-1 alpha targeted genes which play a role in angiogenesis, glucose metabolism, cell proliferation and metastasis. The enzyme encoded by this gene also regulates G-protein coupled receptor signaling by mediating the deubiquitination of beta-2 adrenergic receptor (ADRB2). This enzyme is a ubiquitously expressed thiolester hydrolase. Alternative splicing results in multiple transcript variants encoding the same protein. [provided by RefSeq, Jan 2013]
Transcript Variant: This variant (2) has different 5' and 3' non-coding, terminal exons, compared to transcript variant 1. All three transcript variants encode the same protein.
Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The extent of this transcript is supported by transcript alignments and orthologous data.