

## Product datasheet for **SC328155**

### BLCAP (NM\_001167823) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	BLCAP (NM_001167823) Human Untagged Clone
Tag:	Tag Free
Symbol:	BLCAP
Synonyms:	BC10
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	<p>&gt;NCBI ORF sequence for NM_001167823, the custom clone sequence may differ by one or more nucleotides</p> <p>ATGTATTGCCTCCAGTGGCTGCTGCCCCGTCCTCCTCATCCCCAAGCCCTCAACCCCGCC            CTGTGGTTCAGCCACTCCATGTTTCATGGGCTTCTACCTGCTCAGCTTCCTCCTGGAACGG            AAGCCTTGCAACAATTTGTCCTTGGTTTTCTGGCAGCCCTGTTCTTATCTGCTATAGC            TGCTGGGGAAACTGTTTCCTGTACCACTGCTCCGATTCCCGCTTCCAGAATCGGCGCAT            GATCCCGGCGTTGTGGGCACCTAA</p>
Restriction Sites:	Please inquire
ACCN:	NM_001167823
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).


[View online »](#)

<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u>NM_001167823.1, NP_001161295.1</u>
<b>RefSeq Size:</b>	2119 bp
<b>RefSeq ORF:</b>	264 bp
<b>Locus ID:</b>	10904
<b>UniProt ID:</b>	<u>P62952</u>
<b>Cytogenetics:</b>	20q11.23
<b>Protein Families:</b>	Transmembrane
<b>Gene Summary:</b>	<p>This gene encodes a protein that reduces cell growth by stimulating apoptosis. Alternative splicing and the use of alternative promoters result in multiple transcript variants encoding the same protein. This gene is imprinted in brain where different transcript variants are expressed from each parental allele. Transcript variants initiating from the upstream promoter are expressed preferentially from the maternal allele, while transcript variants initiating downstream of the interspersed NNAT gene (GeneID:4826) are expressed from the paternal allele. Transcripts at this locus may also undergo A to I editing, resulting in amino acid changes at three positions in the N-terminus of the protein. [provided by RefSeq, Nov 2015]</p> <p>Transcript Variant: This variant (5) differs in the 5' UTR compared to variant 1. Variants 1 to 7 all encode the same protein.</p>