

## Product datasheet for SC328201

### IFITM10 (NM\_001170820) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	IFITM10 (NM_001170820) Human Untagged Clone
Tag:	Tag Free
Symbol:	IFITM10
Synonyms:	DSPA3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC328201 representing NM_001170820. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG  
 GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC  
 ATGAGGGAGGGAAAACGGGGGCCCATGCATCCTCAGCTTCCGGGGGACTTTGGAGAGGGTCGAGGCT  
 CAGTGGGAGCTGGAGGCCAGGGCCCGGCCAGTCCCAGCCCCGCTGGGAGACCCGGCCAGCACCACG  
 GACGGCGCCAGGAAGCCGAGTCCCCCTGGACGGGGCCTTCTGGATTCCGAGGCCCGGCAGGTTTCG  
 CCCAAGGGCTGCTTCGTTGCGTGTCAAGCCCCCTGCCCTGCAGGCTCCGGCGGCCCTGCCCTGAG  
 CCCTCGGCCTCTCCCCGATGGCGCCCACTGTTCCCATGGAGTCCAAGAGCAGCAAGACCGACAGC  
 GTGCGGGCTGCCGGCGCGCCCCCTGCCTGCAAGCACCTAGCCGAGAAGAAGACGATGACCAACCCACG  
 ACCGTATCGAGGTCTACCCGGACACCACCGAGGTGAACGACTATTACCTGTGGTCCATCTCAACTTC  
 GTCTACCTCAACTTCTGCTGCCTGGGCTTCATCGCCTTGGCTACTCCCTCAAAGTGCGAGACAAGAAG  
 CTTCTCAATGACCTGAATGGAGCCGTGGAGGATGCAAAGACGGCCCGGCTGTTCAACATCACCAGTTCT  
 GCCCTGGCAGCTCCTGCATCATCCTCGTCTTCATCTTCTGCGGTACCCCTCACCAGTACTAA  
 ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT  
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC

Restriction Sites:	SgfI-MluI
ACCN:	NM_001170820
Insert Size:	687 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).


[View online »](#)

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
RefSeq:	<a href="#">NM_001170820.3</a>
RefSeq Size:	3721 bp
RefSeq ORF:	687 bp
Locus ID:	402778
UniProt ID:	<a href="#">A6NMD0</a>
Cytogenetics:	11p15.5
MW:	24.6 kDa