

Product datasheet for **SC328292**

TSPAN3 (NM_001168412) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TSPAN3 (NM_001168412) Human Untagged Clone
Tag:	Tag Free
Symbol:	TSPAN3
Synonyms:	TM4-A; TM4SF8; TSPAN-3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC328292 representing NM_001168412. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGGCCAGTGCGGCATCACCTCCTCAAGACCGTGCTGGTCTTCTCAACCTCATCTTCTGGTTTGTC
ATCATCCTGCTCTTGGTTTTTGTACAGAAGTTGTTGTAGTGGTTTTGGGATATGTTTACAGAGCAAAG
GTGAAAATGAGGTTGATCGCAGCATTAGAAAGTGTATAAGACCTACAATGGAACCAACCTGATGCT
GCTAGCCGGGCTATTGATTATGTACAGAGACAGCTGCATTGTTGTGGAATTCACAACCTCAGACTGG
GAAAATACAGATTGGTTCAAAGAAACCAAAACCAGAGTGTCCCTTAGCTGCTGCAGAGAGACTGCC
AGCAATTGTAATGGCAGCCTGGCCACCCTTCCGACCTCTATGCTGAGGGGTGTGAGGCTCTAGTAGTG
AAGAAGCTACAAGAAATCATGATGCATGTGATCTGGGCCGCACTGGCATTTCAGCTATTCAGCTGCTG
GGCATGCTGTGTGCTTGCATCGTGTGTGCAGAAGGAGTAGAGATCCTGCTTACGAGCTCCTCATCACT
GGCGAACCTATGCATAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-MluI



[View online »](#)

Plasmid Map:



ACCN: NM_001168412

Insert Size: 570 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001168412.1](#)

RefSeq Size: 3679 bp

RefSeq ORF: 570 bp

Locus ID: 10099

UniProt ID: [O60637](#)

Cytogenetics: 15q24.3

Protein Families: Transmembrane

MW: 21.2 kDa

Gene Summary: The protein encoded by this gene is a member of the transmembrane 4 superfamily, also known as the tetraspanin family. Most of these members are cell-surface proteins that are characterized by the presence of four hydrophobic domains. The proteins mediate signal transduction events that play a role in the regulation of cell development, activation, growth and motility. The use of alternate polyadenylation sites has been found for this gene. Multiple alternative transcripts encoding different isoforms have been described. [provided by RefSeq, Dec 2009]

Transcript Variant: This variant (3) lacks an alternate in-frame exon, compared to variant 1, resulting in a shorter protein (isoform 3), compared to isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.