

Product datasheet for **SC328329**

Claudin 7 (CLDN7) (NM_001185022) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Claudin 7 (CLDN7) (NM_001185022) Human Untagged Clone
Tag:	Tag Free
Symbol:	Claudin 7
Synonyms:	CEPTL2; claudin-1; CLDN-7; CPETL2; Hs.84359
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	<p>>NCBI ORF sequence for NM_001185022, the custom clone sequence may differ by one or more nucleotides</p> <pre> ATGGCCAATTCGGGCTGCAGTTGCTGGGCTTCTCCATGGCCCTGCTGGGCTGGGTGGGT CTGGTGGCCTGCACCGCCATCCCGCAGTGGCAGATGAGCTCCTATGCGGGTGACAACATC ATCACGGCCCAGGCCATGTACAAGGGGCTGTGGATGGACTGCGTCACGCAGAGCAGCGGG ATGATGAGCTGCAAAATGTACGACTCGGTGCTCGCCCTGTCCGCGGCTTGACAGGCCACT CGAGCCCTAATGGTGGTCTCCCTGGTGTGGGCTTCTGGCCATGTTTGTGGCCACGATG GGCATGAAGTGCACGCGCTGTGGGGGAGACGACAAAGTGAAGAAGGCCCGTATAGCCATG GGTGGAGGCATAATTTTCATCGTGGCAGGTCTTGCCGCTTGGTAGCTTGCTCCTGGTAT GGCCATCAGATTGTCACAGACTTTTATAACCCTTTGATCCCTACCAACATTAAGTATGAG TTTGGCCCTGCCATCTTTATTGGCTGGGCAGGGTCTGCCCTAGTCATCCTGGGAGGTGCA CTGCTCTCCTGTTCTGTCTGGGAATGAGAGCAAGGCTGGGTACCGTGTACCCCGCTCT TACCCTAAGTCCAACCTCTTCCAAGGAGTATGTGTGA </pre>
Restriction Sites:	Please inquire
ACCN:	NM_001185022
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.


[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001185022.1, NP_001171951.1</u>
RefSeq Size:	1462 bp
RefSeq ORF:	636 bp
Locus ID:	1366
UniProt ID:	<u>O95471</u>
Cytogenetics:	17p13.1
Protein Families:	Transmembrane
Protein Pathways:	Cell adhesion molecules (CAMs), Leukocyte transendothelial migration, Tight junction
Gene Summary:	<p>This gene encodes a member of the claudin family. Claudins are integral membrane proteins and components of tight junction strands. Tight junction strands serve as a physical barrier to prevent solutes and water from passing freely through the paracellular space between epithelial or endothelial cell sheets, and also play critical roles in maintaining cell polarity and signal transductions. Differential expression of this gene has been observed in different types of malignancies, including breast cancer, ovarian cancer, hepatocellular carcinomas, urinary tumors, prostate cancer, lung cancer, head and neck cancers, thyroid carcinomas, etc.. Alternatively spliced transcript variants encoding different isoforms have been found. [provided by RefSeq, May 2010]</p> <p>Transcript Variant: This variant (2) has an alternate 5' UTR sequence, as compared to variant 1. Variants 1 and 2 encode the same isoform (1).</p>