

Product datasheet for **SC328340**

CD84 (NM_001184882) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CD84 (NM_001184882) Human Untagged Clone
Tag:	Tag Free
Symbol:	CD84
Synonyms:	hCD84; LY9B; mCD84; SLAMF5
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	<p>>NCBI ORF sequence for NM_001184882, the custom clone sequence may differ by one or more nucleotides</p> <pre> ATGGCTCAGCACCTATGGATCTTGCTCCTTTGCCTGCAAACCTGTCGGCTTGGGAAA CCAAAAATTACACAGAGTTTAATGGCATCTGTGAACAGCACCTGTAATGTCACACTGACA TGCTCTGTAGAGAAAGAAGAAAAGAAATGTGACATACAATTGGAGTCCCCTGGGAGAAGAG GGTAATGTCCTTCAAATCTTCCAGACTCCTGAGGACCAAGAGCTGACTTACACGTGTACA GCCCAGAACCTGTCTCAGCAACAATTCTGACTCCATCTCTGCCCGGCAGCTCTGTGCAGAC ATCGCAATGGGCTTCCGTAATCACCACACCGGGTTGCTGAGCGTGCTGGCTATGTTCTTT CTGCTTGTCTCATTCTGCTTTCAGTGTTTTTGTCCGTTTGTTCAAGAGAAGACAAGAT GCTGCCTCAAAGAAAACCATATACATATATCATGGCTTCAAGGAACACCCAGCCAGCA GAGTCCAGAATCTATGATGAAATCCTGCAGTCCAAGGTGCTTCCCTCCAAGGAAGAGCCA GTGAACACAGTTTATTCCGAAAGTGCAGTTTGCTGATAAGATGGGGAAAGCCAGCACACAG GACAGTAAACCTCCTGGGACTTCAAGCTATGAAATTGTGATCTAG </pre>
Restriction Sites:	Please inquire
ACCN:	NM_001184882
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.


[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001184882.1, NP_001171811.1</u>
RefSeq Size:	7903 bp
RefSeq ORF:	645 bp
Locus ID:	8832
UniProt ID:	<u>Q9UIB8</u>
Cytogenetics:	1q23.3
Protein Families:	Druggable Genome, Transmembrane
Gene Summary:	<p>This gene encodes a membrane glycoprotein that is a member of the signaling lymphocyte activation molecule (SLAM) family. This family forms a subset of the larger CD2 cell-surface receptor Ig superfamily. The encoded protein is a homophilic adhesion molecule that is expressed in numerous immune cells types and is involved in regulating receptor-mediated signaling in those cells. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Oct 2011]</p> <p>Transcript Variant: This variant (4) lacks two alternate segments, compared to variant 1. The resulting protein (isoform 4) is shorter when it is compared to isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>