

Product datasheet for **SC328347**

PRRG1 (NM_001173490) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PRRG1 (NM_001173490) Human Untagged Clone
Tag:	Tag Free
Symbol:	PRRG1
Synonyms:	PRGP1
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001173490, the custom clone sequence may differ by one or more nucleotides ATGGGGAGGGTTTTCTCACGGGAGAAAAAGCCAATTCCATATTTAAAACGCTACCCAAGA GCTAATGGGTTTTTTGAAGAAAT AAGACAGGGCAACATTGAGCGTGCAAGAAGAA TTCTGTACATTTGAAGAAGCAAGAGAAGCTTTTAAAAAT AATGAAAAACTAAGGAGTTT TGGAGCACCTACACAAAAGCGCAACAAGGGGAGAGTAACCGAGGAAGTACTGGTTTCAG TTTTACCTTACCTTTCCGTTAATCTTTGGCCTCTTCATTATCCTCCTTGTCATTTTCCTA ATCTGGAGATGCTTCTAAGAAACAAAACCTCGTAGACAGACAGTACTGAAGGCCACATT CCTTTCCCTCAGCACCTTAATATTATCACCCACCCCCCACCAGATGAAGTGTGGAC AGCAGTGGATTGTCTCCAGGCTTTCTGGGATATGTAGTTGGGCGCTCAGATTCGTCTCT ACTCGCCTGTCCAATTGTGATCCCCGCCAACCTATGAGGAAGCCACTGGCCAAGTGAAC CTGCAGAGGAGTGAAACAGAACCTCATTTAGACCCACCCCAGAGTATGAGGACATAGTC AACTCCAACCTCAGCCAGTGCCATTCTATGGTGCCTGTGGTCACCACCATCAAATGA
Restriction Sites:	Please inquire
ACCN:	NM_001173490
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001173490.1, NP_001166961.1</u>
RefSeq Size:	4527 bp
RefSeq ORF:	657 bp
Locus ID:	5638
UniProt ID:	<u>O14668</u>
Cytogenetics:	Xp21.1
Protein Families:	Transmembrane
Gene Summary:	<p>This gene encodes a vitamin K-dependent, gamma-carboxyglutamic acid (Gla)-containing, single-pass transmembrane protein. This protein contains a Gla domain at the N-terminus, preceded by a propeptide sequence required for post-translational gamma-carboxylation of specific glutamic acid residues by a vitamin K-dependent gamma-carboxylase. The C-terminus is proline-rich containing PPXY and PXXP motifs found in a variety of signaling and cytoskeletal proteins. This gene is highly expressed in the spinal cord. Several alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Mar 2010]</p> <p>Transcript Variant: This variant (4) uses an alternative donor splice site at the 5' terminal non-coding exon, and is missing another 5' non-coding exon compared to variant 1. Variants 1-4 encode the same isoform (1).</p>