

Product datasheet for SC328370

VEGFA (NM_001171623) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	VEGFA (NM_001171623) Human Untagged Clone
Tag:	Tag Free
Symbol:	VEGFA
Synonyms:	MVCD1; VEGF; VPF
Mammalian Cell Selection:	None
Vector:	pCMV6-XL5
E. coli Selection:	Ampicillin (100 ug/mL)

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

Fully Sequenced ORF:

>OriGene sequence for NM_001171623 edited GCAGCCGGGTAGCTCGGAGGTCGTGGCGCTGGGGGGCTAGCACCAGCGCTCTGTCGGGAGG CGCAGCGGTTAGGTGGACCGGTCAGCGGACTCACCGGCCAGGGCGCTCGGTGCTGGAATT CTGTATTGTTTCTCGTTTTAATTTATTTTTGCTTGCCATTCCCCACTTGAATCGGGCCGA CGGCTTGGGGAGATTGCTCTACTTCCCCAAATCACTGTGGATTTTGGAAACCAGCAGAAA GGGGTGACCGCCGGAGCGCGGCGTGAGCCCTCCCCCTTGGGATCCCGCAGCTGACCAGTC GCGCTGACGGACAGACAGACA: GACACCG: CCCCCAGCCCCAGCTACCACCTCCTCCCCG GCCGGCGGCGGACAGTGGACGCGGCGGCGAGCCGCGGGCAGGGGCCGGAGCCCGCGCCCG GAGGCGGGGTGGAGGGGGTCGGGGGCTCGCGGCGTCGCACTGAAACTTTTCGTCCAACTTC TGGGCTGTTCTCGCTTCGGAGGAGCCGTGGTCCGCGCGGGGGAAGCCGAGCCGAGCGGAG AAGAGAAGGAAGAGGAGAGGGGGCCGCAGTGGCGACTCGGCGCTCGGAAGCCGGGCTCAT TGGCCCGGGCCTCGGGCCGGGGAGGAAGAGTAGCTCGCCGAGGCGCCGAGGAGAGCGGGC CGCCCCACAGCCCGAGCCGGAGAGGGAGCGCGAGCCGCGCCCGGCCCCGGTCGGGCCTCCG AAACCATGAACTTTCTGCTGTCTTGGGTGCATTGGAGCCTTGCCTTGCTGCTCTACCTCC ACCATGCCAAGTGGTCCCAGGCTGCACCCATGGCAGAAGGAGGAGGGCAGAATCATCACG AAGTGGTGAAGTTCATGGATGTCTATCAGCGCAGCTACTGCCATCCAATCGAGACCCTGG TGGACATCTTCCAGGAGTACCCTGATGAGATCGAGTACATCTTCAAGCCATCCTGTGTGC CCCTGATGCGATGCGGGGGGCTGCTGCAATGACGAGGGCCTGGAGTGTGTGCCCACTGAGG AGTCCAACATCACCATGCAGATTATGCGGATCAAACCTCACCAAGGCCAGCACATAGGAG AAGAAAAAAAATCAGTTCGAGGAAAGGGAAAGGGGCAAAAACGAAAGCGCAAGAAATCCC GGTATAAGTCCTGGAGCGTGTACGTTGGTGCCCGCTGCTGTCTAATGCCCTGGAGCCTCC CGCAGACGTGTAAATGTTCCTGCAAAAACACAGACTCGCGTTGCAAGGCGAGGCAGCTTG AGTTAAACGAACGTACTTGCAGATGTGACAAGCCGAGGCGGTGAGCCGGGCAGGAGGAAG GAGCCTCCCTCAGGGTTTCGGGAACCAGATCTCTCACCAGGAAAGACTGATACAGAACGA TCGATACAGAAACCACGCTGCCGCCACCACCACCATCACCATCGACAGAACAGTCCTTAAT CCAGAAACCTGAAATGAAGGAAGAGGAGACTCTGCGCAGAGCACTTTGGGTCCGGAGGGC GAGACTCCGGCGGAAGCATTCCCGGGCGGGTGACCCAGCACGGTCCCTCTTGGAATTGGA TTCGCCATTTTATTTTCTTGCTGCTAAATCACCGAGCCCGGAAGATTAGAGAGTTTTAT TTCTGGGATTCCTGTAG

Restriction Sites:	Please inquire
ACCN:	NM_001171623
Insert Size:	2100 bp

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

Service VEGFA (NM_001171623) Human Untagged Clone – SC328370	
OTI Disclaimer:	Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at <u>custsupport@origene.com</u> or by calling 301.340.3188 option 3 for pricing and delivery.
	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <u>More info</u>
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	 Centrifuge at 5,000xg for 5min. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. Close the tube and incubate for 10 minutes at room temperature. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM 001171623.1, NP 001165094.1</u>
RefSeq Size:	3677 bp
RefSeq ORF:	699 bp
Locus ID:	7422
UniProt ID:	<u>P15692</u>
Cytogenetics:	6p21.1
Protein Families:	Druggable Genome, Secreted Protein
Protein Pathways:	Bladder cancer, Cytokine-cytokine receptor interaction, Focal adhesion, mTOR signaling pathway, Pancreatic cancer, Pathways in cancer, Renal cell carcinoma, VEGF signaling pathway

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

VEGFA (NM_001171623) Human Untagged Clone – SC328370

Gene Summary:

This gene is a member of the PDGF/VEGF growth factor family. It encodes a heparin-binding protein, which exists as a disulfide-linked homodimer. This growth factor induces proliferation and migration of vascular endothelial cells, and is essential for both physiological and pathological angiogenesis. Disruption of this gene in mice resulted in abnormal embryonic blood vessel formation. This gene is upregulated in many known tumors and its expression is correlated with tumor stage and progression. Elevated levels of this protein are found in patients with POEMS syndrome, also known as Crow-Fukase syndrome. Allelic variants of this gene have been associated with microvascular complications of diabetes 1 (MVCD1) and atherosclerosis. Alternatively spliced transcript variants encoding different isoforms have been described. There is also evidence for alternative translation initiation from upstream non-AUG (CUG) codons resulting in additional isoforms. A recent study showed that a C-terminally extended isoform is produced by use of an alternative inframe translation termination codon via a stop codon readthrough mechanism, and that this isoform is antiangiogenic. Expression of some isoforms derived from the AUG start codon is regulated by a small upstream open reading frame, which is located within an internal ribosome entry site. The levels of VEGF are increased during infection with severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2), thus promoting inflammation by facilitating recruitment of inflammatory cells, and by increasing the level of angiopoietin II (Ang II), one of two products of the SARS-CoV-2 binding target, angiotensin-converting enzyme 2 (ACE2). In turn, Ang II facilitates the elevation of VEGF, thus forming a vicious cycle in the release of inflammatory cytokines. [provided by RefSeq, Jun 2020] Transcript Variant: This variant (1) represents the longest transcript. This variant can initiate translation from four non-AUG (CUG) sites, and also from a downstream, in-frame AUG. The isoform (i, also referred to as VEGF206) represented in this RefSeg is derived from the AUG start codon, and is shorter at the N-terminus, compared to isoform a. CCDS Note: This CCDS, which is supported by the mRNA S85192.1, represents a VEGFA isoform derived from an AUG start codon. Alternative translation initiation from multiple upstream non-AUG(CUG) codons is also possible for this gene, as indicated in PMIDs:11352659, 11563986 and 11731620. The longest isoform encoded by this variant, which is derived from the 5'-most CUG start codon, is represented by CCDS 34457.1.

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US