

Product datasheet for **SC328404**

PHETA1 (NM_001177997) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PHETA1 (NM_001177997) Human Untagged Clone
Tag:	Tag Free
Symbol:	PHETA1
Synonyms:	FAM109A; IPIP27A; SES1
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001177997, the custom clone sequence may differ by one or more nucleotides

```

ATGAAGCTGAACGAGCGCAGCCTGGCCTTCTACGCCACCTGTGACGCCCCGGTGGACAAT
GCAGGCTTCTGTACAAGAAGGGTGGGCGGCACGCGCCTACCACCGCGCTGGTTCGTG
CTGCGCGGGAACATGCTCTTCTACTTCGAGGACGCTGCCAGCCGTGAGCCCGTGGGCGTC
ATCATCTGGAGGGCTGCACTGTGGAGCTGGTGGAGGCCGCCGAGGAGTTCGCCTTCGCT
GTGCGCTTTCGCGGGACCCGGGCGCGCACCTACGTGCTGGCCGCTGAGAGTCAGGATGCC
ATGGAGGGCTGGGTCAAGGCCCTGTCGCGTGCCAGCTTCGACTACCTGCGGCTGGTGGTG
CGCGAGCTGGAGCAGCAGCTGGCGGCTGTACGTGGCGGGGGTGGCATGGCCCTGCCCCAG
CCCCAGCCCCAGTCCCTGCCCTTGCCCCCGTCCCTGCCCTCTGCCCTGGCCCCAGTCCA
TCCCTGCCTTCTGCCCCAGCCCCGGTCCCAGCCCTGCCCTGCCCGCCGCCCCAGTGCC
CTCCCGCCCAAGGAGAATGGCTGCGCTGTCTGGAGCACTGAGGCCACCTTCAGGCCTGGA
CCCGAGCCCCCTCCACCACCGCCTCGCCGCCGGGCTCGGCACCCACGGGCCCTTGGAC
ATGGCCCCCTTCGCGCGCTGCACGAGTGCTATGGCCAGGAGATCCGGGCCCTGCGTGCC
CAGTGGCTCAGCAGCCGGGTCCAGCCCTGA
  
```

Restriction Sites:	Please inquire
ACCN:	NM_001177997
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).


[View online »](#)

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001177997.1</u> , <u>NP_001171468.1</u>
RefSeq Size:	3149 bp
RefSeq ORF:	750 bp
Locus ID:	144717
UniProt ID:	<u>Q8N4B1</u>
Cytogenetics:	12q24.12
Gene Summary:	<p>This gene encodes a protein that localizes to the endosome and interacts with the enzyme, inositol polyphosphate 5-phosphatase OCRL-1. Alternate splicing results in multiple transcript variants. [provided by RefSeq, May 2010]</p> <p>Transcript Variant: This variant (3) differs in the 5' UTR, lacks a portion of the 5' coding region and initiates translation at a downstream start codon, compared to variant 1. The encoded isoform (2) has a shorter N-terminus, compared to isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>