

Product datasheet for **SC328417**

MBNL3 (NM_001170704) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MBNL3 (NM_001170704) Human Untagged Clone
Tag:	Tag Free
Symbol:	MBNL3
Synonyms:	CHCR; MBLX; MBLX39; MBXL
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001170704, the custom clone sequence may differ by one or more nucleotides ATGTTCCGCCAGCAGATGCAGCTTATGCTCCAAAACGCTCAAATGTCATCACTTGGTTCT TTTCCTATGACTCCATCAATTCCAGCTAATCCTCCCATGGCTTTCAATCCTTACATACCA CATCCTGGGATGGCCTCGTTCCTGCAGAACTGTACCAATACACCTGTTCTGATTCTCT GGAAACCCACCTCTTGAATGCCAGGAGCTGTTGGCCAAAACCTGATGCGTTCAGATAAA CTGGAGGTTTGCCGAGAAATTCAGCGTGGAAATGTACCCGTGGGAGAAATGATTGCCGC TATGCTCACCTACTGATGCTTCCATGATTGAAGCGAGTGATAACTGTGACAATCTGC ATGGATTACATCAAAGGTCGATGCTCGCGGGAGAAATGCAAGTACTTTCATCCTCCTGCA CACTTGAAGCCAGACTCAAGGCAGCTCATCATCAGATGAACCATTGAGCTGCCTCTGCC ATGGCCCTGCAGCCTGGTACACTGCAACTGATACCAAAGAGATCAGCACTGGAAAAGCCC AATGGTGCCACCCGGTCTTTAATCCCACTGTTTTCCACTGCCAACAGGCTCTGACTAAC CTGCAGCTCCACAGCCGGCATTATCCCTGCAGGGCCAATACTGTGCATGGCACCCGCT TCAAATATTGTGCCCATGATGCACGGTGTACACCTACCACTGTGTCTGCAGCAACAACA CCTGCCACCAGCGTTCGGTTCGCTGCACCAACTACAGGCAATCAGCTGAAATTCTGA
Restriction Sites:	Please inquire
ACCN:	NM_001170704
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001170704.1 , NP_001164175.1
RefSeq Size:	11539 bp
RefSeq ORF:	777 bp
Locus ID:	55796
UniProt ID:	Q9NUK0
Cytogenetics:	Xq26.2
Gene Summary:	<p>This gene encodes a member of the muscleblind-like family of proteins. The encoded protein may function in regulation of alternative splicing and may play a role in the pathophysiology of myotonic dystrophy. Alternatively spliced transcript variants have been described. [provided by RefSeq, Dec 2009]</p> <p>Transcript Variant: This variant (6) differs in the 5' UTR, lacks a portion of the 5' coding region, and initiates translation at a downstream start codon, compared to variant 1. The encoded isoform (5) is shorter than isoform 1. Variants 5 and 6 encode the same isoform (5). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>