

Product datasheet for **SC328427**

PQBP1 (NM_001167989) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PQBP1 (NM_001167989) Human Untagged Clone
Tag:	Tag Free
Symbol:	PQBP1
Synonyms:	MRX2; MRX55; MRXS3; MRXS8; NPW38; RENS1; SHS
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001167989, the custom clone sequence may differ by one or more nucleotides ATGCCGCTGCCCGTTGCGCTGCAGACCCGCTTGGCCAAGAGAGGCATCCTCAAACATCTG GAGCCTGAACCAGAGGAAGAGATCATTGCCGAGGACTATGACGATGATCCTGTGGACTAC GAGGCCACCAGGTTGAGGGCCTACCACCAAGCTGGTACAAGGTGTTCCGACCCCTCCTGC GGGCTCCCTTACTACTGGAATGCAGACACAGACCTTGATCCTGGCTCTCCCCACATGAC CCCAACTCCGTGGTTACCAAATCGGCCAAGAAGCTCAGAAGCAGTAATGCAGATGCTGAA GAAAAGTTGGACCGGAGCCATGACAAGTCGGACAGGGGCCATGACAAGTCGGACCGCAGC CATGAGAACTAGACAGGGGCCACGACAAGTCAGACCGGGGCCACGACAAGTCTGACAGG GATCGAGAGCGTGGCTATGACAAGGTAGACAGAGAGAGAGAGCGAGACAGGAACGGGAT CGGGACCGCGGTATGACAAGGCAGACCGGAAGAGGGCAAGAACGGGCCACCATCGC CGGGAGGAGCTGGCTCCCTATCCAAGAGCAAGAAGGTAAGCCGAAAGGATGAAGAGTTA GACCCCATGGACCTAGCTCATACTCAGACGCCCCCGGGGCACGTGGTCAACAGGACTC CCCAAGCGGAATGAGGCCAAGACTGGCGCTGACACCACAGCAGCTGGGCCCTCTCCAG CAGCGCCGTATCCATCCCCAGGGGCTGTGCTCCGGGCCAATGCAGAGGCCTCCCGAACC AAGCAGCAGGATTGA
Restriction Sites:	Please inquire
ACCN:	NM_001167989
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001167989.1</u> , <u>NP_001161461.1</u>
RefSeq Size:	994 bp
RefSeq ORF:	795 bp
Locus ID:	10084
UniProt ID:	<u>O60828</u>
Cytogenetics:	Xp11.23
Protein Families:	Transcription Factors
Protein Pathways:	Spliceosome
Gene Summary:	<p>This gene encodes a nuclear polyglutamine-binding protein that is involved with transcription activation. The encoded protein contains a WW domain. Mutations in this gene have been found in patients with Renpenning syndrome 1 and other syndromes with X-linked cognitive disability. Multiple alternatively spliced transcript variants that encode different protein isoforms have been described for this gene.[provided by RefSeq, Nov 2009]</p> <p>Transcript Variant: This variant (8) differs in the 5' UTR and uses a different splice site in the 3' coding region, compared to variant 1. The resulting protein (isoform 4) is shorter by 1 aa when it is compared to isoform 1.</p>