

## Product datasheet for SC328441

### C5orf35 (SETD9) (NM\_001171990) Human Untagged Clone

#### Product data:


Product Type:	Expression Plasmids
Product Name:	C5orf35 (SETD9) (NM_001171990) Human Untagged Clone
Tag:	Tag Free
Symbol:	SETD9
Synonyms:	C5orf35
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC328441 representing NM_001171990. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCCTGGCCGTCTGCTGCGGGCCTGTGGCAGCGATGGCGCCGTTACAAGTACCGCTTCGTTCCCTGG
ATCGCACTGAACCTAAGCCACAACCCGAGGACCCTCCGATATGTTCCAGAGGAATCCAAAGACAAAGTT
ATCTCAGATGAAGATGCTCCTAGGAACATTACTGAAAGTTTTCCAGGCTCTATTCTAAATGATTTCAAT
AAACAATCAGAAATCTTGCTATGCTTCCAGAATCTGTTAAATCAAATATCAAGACCTACTGGCAGTT
GAACATCAAGGGGTGAACTGCTTGAAAACAGACATCAACAGCAAAGTACCTTTAAACCAGAAGAAATT
CTTTACAAGACTTTGGGTTTCAGTGTTGCCAAGCAACTAGCTCATTGATTTCTGCTGAAAAGGTGTC
TTCGTTACTAAAGGATTGGTACAAAAGGCGCAGTCGTATCTATGTATCCCTGGTACAGTATATCAGAAG
TATGAGCCGATCTTTTTCCAGTCCATTGGAAATCCGTTTATTTTTAGATGCCTGGATGGGGTACTCATT
GATGGGAATGACAAAGGGATATCAAAGTTGTGTACAGATCTTGCAATGGGAGGGATCGACTCGGCCCT
TTAAAAATGAGTGATAGTACATGGCTAACGTCAGAAATTCATAACCCCTCTGGCTGTGGGACAGTATGTC
ACAATTTGTTCCAATGACAGAGCAGCTAATGTCTGTTATCAGGAATTTGATGTGCCTGCAGTTTTCCCT
ATAGAACTGAAGCAGTATCTTCCAAACATTGCCTACAGCTATGACAAACAAGGTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

**Plasmid Map:**


**ACCN:** NM\_001171990

**Insert Size:** 816 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001171990.2](#)

**RefSeq Size:** 4766 bp

RefSeq ORF: 816 bp  
Locus ID: 133383  
UniProt ID: Q8NE22  
Cytogenetics: 5q11.2  
MW: 31 kDa