

Product datasheet for **SC328446**

PRR7 (NM_001174102) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PRR7 (NM_001174102) Human Untagged Clone
Tag:	Tag Free
Symbol:	PRR7
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001174102, the custom clone sequence may differ by one or more nucleotides ATGGTGATGTCCCAGGGCACCTACACGTTCTCAGTGTCTCGCCGGCTTCTGGCTCATC TGGGGTCTCATCGTCTGCTGCTGCTTCTGCAGTTCTCTGCGCCGCCCTCAAACGG CGCCAGGAGGAGCGACTGCGCGAGCAGAACCTGCGCGCCCTAGAGCTGGAGCCCCCTCGAA CTCGAGGGCAGTCTGGCCGGGAGCCCCCGGGCTGGCGCCGCCGAGCCACCACACAC CGTAGCCGCTGGAGGCGCCGGCTCACGCGCACTCGCATCCGCACGTGCACGTGCACCCG CTGCTGCACCACGGGCCCGCGCAGCCGACGCGCACGCGCACCCACACCCGACCCACCAC GCGCTCCCGCACCCGCGCTACGCACCTGTGGTGCCGCCACGGCCCTGGAGTACCCG CGCCAAGCGGAATCGGACATGTCAAACCACCGTGTACGAAGAGGCGGTGCTGATGGCA GAGCCCGCCGCGCCCTATAGCGAGGTGCTCACGGACACGCGCGCCCTCTACCGCAAGATC GTACGCCCCTTCTGAGTCGCGCGACAGCGCGGAGAAGCAGGAGCAGCCGCTCCAGC TACAAGCCGCTCTTCTGGACCGGGGCTACACCTCGGCGCTGCACCTGCCAGCGCCCT CGGCCCGCGCCGCTGCCAGCCCTCTGCCTGCAGGCCGACCGTGGCCCGGGTCTTC CCCAGCTGGACCGACTCAGAGCTCAGCAGCCGCGAGCCCTGGAGCACGGAGCTTGGCGT CTGCCGCTCCATCCCCTTGTTCGGGAGGACTACAGCCGTATAG
Restriction Sites:	Please inquire
ACCN:	NM_001174102
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.



Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001174102.1, NP_001167573.1</u>
RefSeq Size:	1338 bp
RefSeq ORF:	825 bp
Locus ID:	80758
UniProt ID:	<u>Q8TB68</u>
Cytogenetics:	5q35.3
Protein Families:	Transmembrane
Gene Summary:	<p>Acts as a synapse-to-nucleus messenger to promote NMDA receptor-mediated excitotoxicity in neurons in a JUN-dependent manner (By similarity). Inhibits ubiquitination-mediated degradation and promotes phosphorylation and transcriptional activity of transcription factor JUN (PubMed:27458189). Might play a redundant role in the regulation of T cell receptor signaling (PubMed:21460222). Might promote apoptosis in T cells (PubMed:21460222). [UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (3) differs in the 5' UTR compared to variant 1. All variants (1-3) encode the same protein. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>