

Product datasheet for **SC328470**

MINPP1 (NM_001178118) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	MINPP1 (NM_001178118) Human Untagged Clone
Tag:	Tag Free
Symbol:	MINPP1
Synonyms:	HIPER1; MINPP2; MIPP
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC328470 representing NM_001178118. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTGCTTATTTTCAGCTGTGCGGATTAGTAAGATATATGGAGTTTGGACCTCCAACAGTTAATGATAAA
CTAATGAGATTTTTTGTACTGTGAGAAGTTTTAACTGAAGTAGAAAAAATGCTACAGCTCTTTAT
CACGTGGAAGCCTTCAAACCTGGACCAGAAATGCAGAACATTTTAAAAAAGTTGCAGCTACTTTGCAA
GTGCCAGTAAATGATTTAAATGCAGATTAATTCAAGTAGCCTTTTTCACCTGTTCAATTTGACCTGGCA
ATTAAGGTGTTAAATCTCCTTGGTGTGATGTTTTTACATAGATGATGCAAAGGTATTAGAATATTTA
AATGATCTGAAACAATATTGAAAAGAGGATATGGGTATACTATTAACAGTCGATCCAGCTGCACCTTG
TTTCAGGATATCTTTCAGCACTTGACAAAGCAGTTGAACAGAAACAAAGGTCTCAGCCAATTTCTCT
CCAGTCATCCTCCAGTTTGGTCATGCAGAGACTCTTCTCCACTGCTTTCTCTCATGGGCTACTTCAA
GACAAGGAACCCCTAACAGCGTACAATTACAAAAACAAATGCATCGGAAGTTCCGAAGTGGTCTCATT
GTACCTTATGCCTCGAACCTGATATTTGTGCTTTACCACTGTGAAAAATGCTAAGACTCCTAAAGAACA
TTCCGAGTGCAGATGTTATTAATGAAAAGGTGTTACCTTTGGCTTACTACAAGAACTGTTTCATTT
TATGAAGATCTGAAGAACCCTACAAGGACATCCTTCAGAGTTGTCAAACCAGTGAAGAATGTGAATTA
GCAAGGGCTAACAGTACATCTGATGAACTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_001178118

Insert Size: 861 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM_001178118.1

RefSeq Size: 2225 bp

RefSeq ORF: 861 bp

Locus ID:	9562
UniProt ID:	Q9UNW1
Cytogenetics:	10q23.2
Protein Families:	Druggable Genome
Protein Pathways:	Inositol phosphate metabolism
MW:	33.1 kDa
Gene Summary:	<p>This gene encodes multiple inositol polyphosphate phosphatase; an enzyme that removes 3-phosphate from inositol phosphate substrates. It is the only enzyme known to hydrolyze inositol pentakisphosphate and inositol hexakisphosphate. This enzyme also converts 2,3-bisphosphoglycerate (2,3-BPG) to 2-phosphoglycerate; an activity formerly thought to be exclusive to 2,3-BPG synthase/2-phosphatase (BPGM) in the Rapoport-Luebering shunt of the glycolytic pathway.[provided by RefSeq, Sep 2009]</p> <p>Transcript Variant: This variant (3) differs in the 5' UTR and coding sequence compared to variant 1. The resulting isoform (3) has a shorter and distinct N-terminus compared to isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>