

Product datasheet for **SC328494**

PEX2 (NM_001172086) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PEX2 (NM_001172086) Human Untagged Clone
Tag:	Tag Free
Symbol:	PEX2
Synonyms:	PAF1; PBD5A; PBD5B; PMP3; PMP35; PXMP3; RNF72; ZWS3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC328494 representing NM_001172086. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGCTTCCAGAAAAGAGAATGCGAAGAGTGCAAACAGAGTGCTAAGAATAAGCCAGTTGGATGCACTT
GAATAACAAGGCCCTGGAGCAGCTAGTTTGGTCCCAGTTTACTCAGTGCTTTCATGGATTTAAACCT
GGGCTGTAGCTCGCTTTGAGCCAGAGGTGAAAGCGTGCTTATGGGTTTTCTGTGGAGATTCACCATC
TACTCAAAAATGCCACAGTGGGACAGTCAGTTTGAATATTAAGTACAAAATGATTTTTCCCTAAC
CTGAGATATCAGCCACCCAGTAAAAATCAAAAATCTGGTATGCTGTTTGTACAATTGGTGGCAGGTGG
TTAGAAGAACGATGCTATGATTTGTTTCGAAACCATCATTAGCATCATTGGGAAAGTCAAGCAGTGT
GTGAATTTGTGATTGGACTTTTGAAATTAGGTGGGCTGATTAATTTTTGATTTTCCTTCAGAGGGGA
AAGTTTGCAACTTTGACAGAACGCTCCTAGGTATTCATTCTGTATTTTGCAAGCCTCAAAACATACGT
GAAGTTGGCTTTGAATACATGAATAGGGAACCTCTCTGGCATGGTTTTGTGAATTTCTGATTTTCTC
TTACCATTATCAATGTCCAGAAGTTGAAAGCCAAGCTGTCTTCATGGTGTATTCTCTTACTGGTGCA
CCTAATAGTGACAATACATTAGCCACAGTGGCAAAGAATGCGCTCTATGTGGAGAGTGGCCACCATG
CCTCACACCATAGGATGTGAGCATATTTCTGTATTTCTGTGCTAAGAGTAGTTTCTTATTTGACGTG
TACTTTACTTGTCTAAGTGTGGCACAGAAGTACACAGTCTGCAGCCACTGAAATCAGGAATCGAGATG
TCAGAAGTAAATGCTCTTAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:



ACCN: NM_001172086

Insert Size: 918 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001172086.1](#)

RefSeq Size: 4385 bp

RefSeq ORF: 918 bp

Locus ID:	5828
UniProt ID:	P28328
Cytogenetics:	8q21.13
Protein Families:	Druggable Genome, Transmembrane
MW:	34.9 kDa
Gene Summary:	<p>This gene encodes an integral peroxisomal membrane protein required for peroxisome biogenesis. The protein is thought to be involved in peroxisomal matrix protein import. Mutations in this gene result in one form of Zellweger syndrome and infantile Refsum disease. Alternative splicing results in multiple transcript variants encoding the same protein. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (3) differs in the 5' UTR compared to variant 1. Variants 1, 2, 3 and 4 encode the same protein. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>