

Product datasheet for **SC328581**

SPATA22 (NM_001170695) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SPATA22 (NM_001170695) Human Untagged Clone
Tag:	Tag Free
Symbol:	SPATA22
Synonyms:	NYD-SP20; NYDSP20
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001170695, the custom clone sequence may differ by one or more nucleotides <pre> ATGAAGCGAAGCCTAAATGAAAATTCAGCTCGAAGTACAGCAGGCTGTTGCCTGTTCCG TTGTTCAATCAGAAAAAGAGGAACAGACAGCCATTAACCTCTAATCCACTTAAAGATGAT TCAGGTATCAGTACCCCTTCTGACAATTATGATTTTCCTCCTCTACCTACAGATTGGGCC TGGGAAGCTGTGAATCCAGAGTTGGCTCCTGTAATGAAAACAGTGGACACCGGGCAAATA CCACATTCAGTTTCTCGTCTCTGAGAAGTCAAGATTCTGTCTTTAACTCTATTCAATCA AATACTGGAAGAAGCCAGGGTGGTTGGAGCTACAGAGATGGTAACAAAAATACCAGCTTG AAAACCTTGAATAAAAAATGATTTTAAGCCTCAATGTAACGAACAACTTAGTGGCAAAT GATGGAAAAAATCTTGTCAGTGAGTTCGGGAGCTCAACAACAAAAACAATTAAGAATA CCTGAACCTCCTAACTTATCTCGCAACAAAGAAACCGAGCTACTCAGACAAACACATTCA TCAAAAATATCTGGCTGCACAATGAGAGGGCTAGACAAAAACAGTGCACCTACAGACACTT AAGCCCAATTTTCAACAAATCAATATAAGAAACAAATGTTGGATGATATTCCAGAAGAC AACACCCTGAAGGAAACCTCATTGTATCAGTTACAGTTTAAGGAAAAAGCTAGTTCTTTA AGAATTATTTCTGCAGTTATTGAAAGCATGAAGTATTGGCGTGAACATGCACAGAAAACT GTACTTCTTTTGAAGTATTAGCTGTTCTTGATTACAGCTGTTACACCTGGCCCATATTAT TCGAAGACTTTTCTTATGAGGGATGGGAAAAATACTCTGCCTTGTGTCTTTATGAAATC GATCGTGAACCTCCGAGACTGATTAGAGGCCGAGTTTCATAGATGTGTTGGCAACTATGAC CAGAAAAAGAACATTTTCCAATGTGTTTCTGTCAGACCGCGTCTGTTTCTGAACAAAAA ACTTTCCAGGCATTTGTCAAAATTGCAGATGTTGAGATGCAGTATTATATTAATGTGATG AATGAAACTTAA </pre>
Restriction Sites:	Please inquire
ACCN:	NM_001170695


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001170695.1, NP_001164166.1</u>
RefSeq Size:	1417 bp
RefSeq ORF:	1092 bp
Locus ID:	84690
UniProt ID:	<u>Q8NHS9</u>
Cytogenetics:	17p13.2
Gene Summary:	Meiosis-specific protein required for homologous recombination in meiosis I. [UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (5) differs in the 5' UTR and contains an alternate exon which results in a frameshift, compared to variant 1. The resulting protein (isoform 2) is longer and has a distinct C-terminus, compared to isoform 1.