

## Product datasheet for **SC328626**

### SSH3BP1 (ABI1) (NM\_001178123) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	SSH3BP1 (ABI1) (NM_001178123) Human Untagged Clone
Tag:	Tag Free
Symbol:	ABI1
Synonyms:	ABI-1; ABLBP4; E3B1; NAP1BP; SSH3BP; SSH3BP1
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001178123, the custom clone sequence may differ by one or more nucleotides ATGGCAGAGCTGCAGATGTTACTAGAGGAGGAGATCCCGTCTGGCAAGAGGGCGCTGATC GAGAGTTACCAGAACCTGACTCGGGTGGCAGACTACTGTGAAAACAACACATACAGGCT ACAGACAAGAGAAAAGCTTTAGAGGAGACCAAAGCCTATACAACCAATCTCTAGCTAGT GTTGCTTATCAAATAAATGCATTGGCCAACAATGTACTCCAGTTGCTGGATATCCAAGCC TCTCAGCTTCGGAGAATGGAGTCTCCATCAATCATATCTCACAGACTGTGGATATTCAT AAGGAGAAAAGTGGCACGAAGAGAGATTGGTATTTTGACAACAAATAAGAATACATCAAGA ACTCACAAAATAATAGCACCTGCGAATATGGAGCGCCCTGTAAGGTATATTCGGAAACCT ATCGATTACACAGTTCTGGATGATGTGGCCATGGTGCAAGTGGCTAAAAGCCAAGCAT GGAAATAACCAGCCTGCAAGAAGTGGCACACTGTGAGAACAAATCCTCCTACTCAGAAA CCGCCAAGTCTCCCATGTCAGGCCGGGAACACTGGGACGGAATACTCCTTATAAAACC CTGGAACCTGTTAAACCCCAACAGTTCCTAATGACTATATGACCAAGTCTGCTAGGCTT GGAAGTCAGCATAGTCCAGGCAGGACAGCATCTTTAAATCAGAGACCAAGGACACACAGT GGAAGTAGTGGAGGAAGTGAAGTCGAGAAAACAGTGGTAGCAGTAGTATTGGCATTCCC ATTGCTGTGCCTACACCTTCGCCACCCACTATTGGACCAGTTGCTGATAGTCCAACCTCCA CCGCCACCACCTCCACCAGATGACATTCCTGTTGATGACTCTCCACCTCCCCACCA CCACCACAGTGGATTATGAAGATGAGGAGGCTGCAGTAGTTCAGTATAATGATCCATAT GCAGATGGGGATCCTGCTTGGGCCCAAGAATTATATTGAGAAAGTTGTTGCAATATAT GATTATACAAAAGACAAGGATGATGAGCTGTCTTTATGGAGGGTGAATCATTATGTT ATAAAGAAGAATGATGATGGCTGGTATGAAGGAGTCTGCAATCGAGTGAAGTGGTCTGTT CCTGGAACTATGTTGAATCAATCATGCACTATACTGATTAA
Restriction Sites:	Please inquire
ACCN:	NM_001178123



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001178123.1</a></u> , <u><a href="#">NP_001171594.1</a></u>
<b>RefSeq Size:</b>	3380 bp
<b>RefSeq ORF:</b>	1182 bp
<b>Locus ID:</b>	10006
<b>UniProt ID:</b>	<u><a href="#">Q8IZP0</a></u>
<b>Cytogenetics:</b>	10p12.1
<b>Gene Summary:</b>	<p>This gene encodes a member of the Abelson-interactor family of adaptor proteins. These proteins facilitate signal transduction as components of several multiprotein complexes, and regulate actin polymerization and cytoskeletal remodeling through interactions with Abelson tyrosine kinases. The encoded protein plays a role in macropinocytosis as a component of the WAVE2 complex, and also forms a complex with EPS8 and SOS1 that mediates signal transduction from Ras to Rac. This gene may play a role in the progression of several malignancies including melanoma, colon cancer and breast cancer, and a t(10;11) chromosomal translocation involving this gene and the MLL gene has been associated with acute myeloid leukemia. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene, and a pseudogene of this gene is located on the long arm of chromosome 14. [provided by RefSeq, Sep 2011]</p> <p>Transcript Variant: This variant (10) lacks three alternate in-frame coding exons, compared to variant 1. The resulting protein (isoform j) has the same N- and C-termini but is shorter when it is compared to isoform a. This isoform has also been called 'variant 21'. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>