

Product datasheet for SC328635

FAM53A (NM_001174070) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	FAM53A (NM_001174070) Human Untagged Clone
Tag:	Tag Free
Symbol:	FAM53A
Synonyms:	DNTNP
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC328635 representing NM_001174070. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGTCACACTCATCACTGAGAAGCTGCAGAGCCAGAGCCTGGACGACCTCACCTGCAAGGCGGAGGCT
GGCCCGTTGCAGTATTCTGCGGAAACCCTGAACAAGAGCGGTCGTCTGTTCCCTTTGGAGCTCAACGAC
CAGAGTCCCTGGAAGGTCTTCAGTGGAGGACCGCCGTCAGAAGCCAGGCAGCCACGGGCCCTGATTTTC
TCCTTCTGCCGGCCTGTCTGCTGCCGCTCACACCATGGGTCTTCAGTGGCAGCCACAGTCCCCGCGC
CCAGGCGCAGGCCTGGGTGCAGCCAGCACCGTGGACCCAGTCAAAGCACAGGCTCGTCCACGGCCCCA
CCGACCAAGCGCATTGCCGGTCTTGTGACAACCGAGGAGCTTGTGCGCTGCCGGTCCCCTGGCGC
CCCGGCAGCTCCAAGGTCTGGACTCCAGTCTCAAAGAGCGGTGCGACAGCGGGGAGTCCACGCGG
CAGGGAAGCCCCGGCGCCGCTCTGCCGAGGAGTGTGTGGTGCACCGTCCCACCTCGCCCGCCACG
CCCCGGCCGCTCTCCGCCAGCGCGGCTTCGTGGACAGCAGCAGGGGAGTGCAGGCTCAGGCCCGCTC
TGGTGTCCCGGAGTCTGCTTGCCTCCACGAGGCGCCCGCTCCCTCTCACAGGAGCGACTCGCG
GGTGGGGCACTCCCCTGCCCTGGGCCAGCAGCAGCCCCACGTCCACGCTGCGCTGGCGGGGCGCCGT
GGGCTGCTCCGGTGGCGCTCACAGCCTTGGCTGCTCAGTGGGAAGAGGAGCCGGCGCAAACGGAGGCGT
GAGGAGGACGCCAGGTGGACACGCCATCCTTGGACTTCTGAAAATGACCCAGACTTTAAAAAATCA
AAAAGCCTTTGCTCCCTCAATTACGAAGATGACGATGAGGATGACACCCAGTGAAGCGTTCTGTCC
TCCCCATGTGACTCCCGGGGCTCCCTGGCATCACCATGCCTGGCTGACGCCAGAGGGGCTCAGGACC
AGCCCTGTCCACCCCAACCTGTGGGCTCTAGGGAGTCGGTGACCAAGTGTGGTCCCGCAGGAGCAGC
GGGACCCCGTGTGGGACAGTGTGGGGAGGAGGGGCTTCCCCGGGCGCCGCTGGGAGCTGGAC
CTGGAGCAGATCGAGAACAAC TGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites: Sgfl-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_001174070

Insert Size: 1197 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001174070.2](#)

RefSeq Size: 2854 bp

RefSeq ORF: 1197 bp

Locus ID: 152877

UniProt ID: [Q6NSI3](#)

Cytogenetics: 4p16.3

MW: 42.6 kDa

Gene Summary: May play an important role in neural development; the dorsomedial roof of the third ventricle.[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (2) differs in the 5' UTR, compared to variant 1. Both variants 1 and 2 encode the same isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.