

## Product datasheet for **SC328652**

### SRPX (NM\_001170751) Human Untagged Clone


#### Product data:

Product Type:	Expression Plasmids
Product Name:	SRPX (NM_001170751) Human Untagged Clone
Tag:	Tag Free
Symbol:	SRPX
Synonyms:	DRS; ETX1; HEL-S-83p; SRPX1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC328652 representing NM_001170751. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGGGAGCCCCGCACATCGGCCCGCGTCTGCTGCTGCTGCCGCTCTGCTGCTGCTGCTGCTGCTG
CGCGTCCCAGCCAGCCGAGCTTCCCAGGATCGGGAGACTCACCCTAGAAGACGATGAAGTCGGGTAT
TCACACCCTAGATATAAAGATACCCCGTGGTCTCCCCATCAAGGTGAAGTATGGGGATGTGTACTGC
AGGGCCCTCAAGGAGGATACTACAAAACAGCCCTGGGAACCAGGTGCGACATTCGCTGCCAGAAGGGC
TACGAGCTGCATGGCTTCCCTACTGATCTGCCAGTCAAACAACGATGGTCTGACAAGGTCATCTGC
AAACATATGGAACCTCCTAGAATCAAGTGCCCAAGTGTGAAGGAACGCATTGCAGAACCCAACTG
ACAGTCCGGGTGTCTGGGAGACACCCGAAGGAAGAGACACAGCAGATGGAATCTTACTGATGTCATT
CTAAAAGGCCCTCCCCCAGGCTCCAACCTTCCAGAAGGAGACCACAAGATCCAGTACACAGTCTATGAC
AGAGCTGAGAATAAGGGCACTTGCAAAATTCGAGTTAAAGTAAGAGTCAAACGCTGTGGCAAACCTCAAT
GCCCCAGAGAATGGTTACATGAAGTGCTCCAGCGACGGTGATAATTATGGAGCCACCTGTGAGTTCTCC
TGCATCGGCGGCTATGAGCTCCAGGGTAGCCCTGCCGAGTATGTCAATCCAACCTGGCTGGTCTGGC
ACGGAGCCACCTGTGCAGCCATGAACGTCAATGTGGGTGTCAGAACGGCAGCTGCATCTCTGGATCAG
TTTTATGAGAAAAGGAGACTCCTCATTGTGTCCACACCCACAGCCGAAACCTCCTTACCGGCTCCAG
CTAGGAATGCTGCAGCAAGCACAGTGTGGCCTTGATCTTCGACACATCACCGTGGTGGAGCTGGTGGGT
GTGTTCCCGACTCTCATTGGCAGGATAGGAGCAAAGATTATGCCTCCAGCCCTAGCGCTGCAGCTCAGG
CTGTTGCTGCGAATCCCCTCTACTCCTCAGTATGGTGTAGTGGATAAGCATGGCATGGACAAAGAG
CGCTATGTCTCCCTGGTGTGCTGTGGCCCTGTTCAACCTGATTGACACTTTTCCCTTGAGAAAAGAA
GAGATGGTCTACAAGCCGAAATGAGCCAGACCTGTAACACCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



**Plasmid Map:**


**ACCN:** NM\_001170751

**Insert Size:** 1218 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001170751.1](#)

**RefSeq Size:** 1743 bp

**RefSeq ORF:** 1218 bp

**Locus ID:** 8406

**UniProt ID:** [P78539](#)

**Cytogenetics:** Xp11.4

**MW:** 45.1 kDa

**Gene Summary:** May be involved in phagocytosis during disk shedding, cell adhesion to cells other than the pigment epithelium or signal transduction.[UniProtKB/Swiss-Prot Function]  
Transcript Variant: This variant (3) lacks an alternate in-frame exon in the 5' coding region compared to variant 1. The resulting protein (isoform 3) is shorter compared to isoform 1.  
Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.