

Product datasheet for **SC328659**

PCYT2 (NM_001184917) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	PCYT2 (NM_001184917) Human Untagged Clone
Tag:	Tag Free
Symbol:	PCYT2
Synonyms:	ET; SPG82
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC328659 representing NM_001184917. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGATCCGGAACGGGCGCGGGGCTGCAGGCGCGCAGAGCAGCCGGGCCCGGGGGCAGGCGCGCCGTG
AGGGTGTGGTGCATGGCTGCTATGACATGGTGCATTACGGCCACTCCAACCAGCTGCGCCAGGCACGG
GCCATGGGTGACTACCTCATCGTAGGCGTGCACACCGATGAGGAGATCGCCAAGCACAAGGGGCCCCCG
GTGTTCACTCAGGAGGAGAGATAACAAGTGGTGCAGGCCATCAAATGGGTGGACGAGGTGGTGCCAGCG
GCTCCCTACGTCACTACACTAGAGACCCCTGGACAAATACAACGTGACTTCTGTGTTACGCGCAATGAC
ATCACCTGACTGTAGATGGCCGGACACCTATGAGGAAGTAAAGCAGGCTGGGAGTACAGAGAATGC
AAGCGCACGCAAGGGGTGTCCACCACAGACCTCGTGGCCGCATGCTGCTGGTAACCAAAGCCCATCAC
AGCAGCCAGGAGATGTCCCTGAGTACCGGGAGTATGCAGACAGTTTTGGCAAGCCCCCTCACCCGATA
CCCGCCGGGGACATACTTTCCTCAGAAGGCTGCTCCAGTGCCCTGGTGGGCGGAACCCCTGGACCGGG
GTATCCAGTTTCTGCAGACATCTCAGAAGATCATCCAGTTTGTCTTGGGAAGGAGCCCAAGCCAGGG
GAGACAGTCATCTATGTGGCTGGTGCTTCGACCTGTTCCACATCGGGCATGTGGACTTCTGGAGAAG
GTGCACAGGCTGGCAGAGAGGCCCTACATCATCGCGGGTTACACTTTGACCAGGAGGTCAATCACTAC
AAGGGGAAGAATAACCCATCATGAATCTGCATGAACGGACTCTGAGCGTGTGGCTGCCGTTACGTTG
TCAGAAGTGGTATTGGAGCCCGTACGCGGTACAGCAGAGCTCCTAAGTCACCTCAAGGTGGACCTG
GTGTGTCACGGCAAGACAGAAATTATCCCTGACAGGGATGGCTCCGACCCATACCAGGAGCCCAAGAGA
AGGGGCATCTTCCGTGAGATTGACAGTGGCAGCAACCTCACACAGACCTCATCGTCCAGCGGATATC
ACCAACAGGTTGGAGTATGAGGCGCAAACAGAAAGGAAGCAAGGAGCTGGCCTTCTGGAGGCT
GCCAGGCAGCAGGCGGCACAGCCCTGGGGAGCGCATGGTGACTTCTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-MluI



[View online >](#)

Plasmid Map:


ACCN: NM_001184917

Insert Size: 1224 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001184917.2](#)

RefSeq Size: 3295 bp

RefSeq ORF: 1224 bp

Locus ID: 5833

UniProt ID: [Q99447](#)

Cytogenetics: 17q25.3

Protein Pathways: Glycerophospholipid metabolism, Metabolic pathways

MW: 45.6 kDa

Gene Summary: This gene encodes an enzyme that catalyzes the formation of CDP-ethanolamine from CTP and phosphoethanolamine in the Kennedy pathway of phospholipid synthesis. Alternative splicing results in multiple transcript variants. [provided by RefSeq, May 2010]
Transcript Variant: This variant (1) encodes the longest isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.