

Product datasheet for SC328665

RSKR (NM_001174103) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	RSKR (NM_001174103) Human Untagged Clone
Tag:	Tag Free
Symbol:	RSKR
Synonyms:	SGK494
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC328665 representing NM_001174103. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGGAGCAGTAAGCTGTCGGCAGGGGCAGCACACCCAGCAGGGGGAACACACCCGGGTGGCTGTCCCT
CACAAGCAGGGTGGCAACATCCGGGGTCCCTGGGCCGAGGCTGGAAGAGCCTCTGGACAGGTTTGGGA
ACCATCAGGTCAGATCTGGAAGAACTCTGGAACTACGGGGCCACCACTATCTGCACCAGGAATCCCTA
AAGCCAGCCCCAGTACTGGTAGAGAAGCCTCTGCCAGAGTGGCCAGTGCCTCAGTTCATCAACCTCTTT
CTACCAGAGTTTCCATTAGGCCATTAGGGGGCAGCAGCAGCTGAAGATTTTAGCCTCGTGGCTAAA
GGCTCCTTTGAACTGTCTCAAGGTGCTAGATTGCACCCAGAAAGCTGTATTTGCAGTGAAGGTGGTG
CCCAAGGTAAAGTCTACAGAGGGATACCGTGAGGCAGTGCAAAGAGGAGGTTAGCATCCAGCGCAG
ATCAACCATCCCTTTGTACACAGCTTGGGGGACAGCTGGCAGGGAAAACGGCACCTTTTCATTATGTGT
AGCTACTGCAGCAGATCTGTACTCCCTTTGGTCCGGCTGTTGGCTGCTTTCTGAGGCTTCCATCCGT
CTCTTTGCTGCCGAGTTGGTGTGTTACTGTGTTATCTCCATGACTTGGGCATCATGCATCGAGATGTG
AAGATGGAGAATATTCTTAGATGAACGAGGCCATCTGAAACTGACAGACTTTGGTCTGTCCCGCCAC
GTGCCCGAGGAGCTCAAGCCTACACTATCTGTGGCACTTTCAGTACATGGCCCCAGAGGTCCTAAGT
GGAGGACCTTACAACCATGTGCTGATTGGTGGTCCCTGGGTGTCTTGTCTTCTCTGCGGACTGGA
AAGTTTCCAGTGGCTGCAGAGAGATCATGTGGCCATGTTGGCAAGTGTGACCCACAGTACTGTGAG
ATCCAGCTTCTCTTAACCGGGCCTCTCACTCCTGCTCCATGAGCTCTTATGCCAGAACCCTCCAT
CGTCTACGTTATCTGCATCACTCCAGGTCCACCTTTCTTTCCGGGTGTGGCCTTCGACCCAGAGCTC
CTACAGAAGCAGCCAGTGAACCTTTGTCACGGAGACACAAGCTACCCAGCCAGTTCAGCGGAGACCATG
CCCTTTGACGACTTTGACTGTGATCTGGAGTCTTCTTGTCTACCCTATCCCTGCTTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-Mlul



[View online >](#)

Plasmid Map:


ACCN: NM_001174103

Insert Size: 1233 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001174103.1](#)

RefSeq Size: 3779 bp

RefSeq ORF: 1233 bp

Locus ID: 124923

Cytogenetics: 17q11.2

Protein Families: Druggable Genome, Protein Kinase

MW: 46.2 kDa