

Product datasheet for **SC328682**

CELF3 (NM_001172648) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CELF3 (NM_001172648) Human Untagged Clone
Tag:	Tag Free
Symbol:	CELF3
Synonyms:	BRUNOL1; CAGH4; ERDA4; ETR-1; TNRC4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC328682 representing NM_001172648. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTGTAGTAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGGATGTGCCTTCTGACATACTGTGCCGGGATTAGCCCTGAAGGCCAGAGCGCCCTGCACGAA
CAGAAGACGCTTCCAGGGATGAACAGGCCGATCCAGGTCAAGCCAGCCGACAGCGAGAGCCGAGGAGAA
GACCGGAAGCTCTTTGTGGGGATGCTAGGGAAGCAGCAGACAGATGAGGACGTCCGGAAGATGTTTGAG
CCCTTCGGGACCATCGACGAGTGCCTGTCTCCGGGGCCAGATGGCACCAGCAAAGGCTGCGCCTTC
GTGAAGTTCAGACCCACGCTGAGGCCAGGCCGATCAACACCCTTACAGCAGCCGGACCCTGCCA
GGTGCCTCGTCCAGCCTGGTGGTGAAGTTTGTGACACTGAGAAGGAGCGAGGTCTCCGCCGATGCAG
CAGGTGGCCACCCAGTTGGGCATGTTAGCCCCATCGCCCTCCAGTTTGGAGCCTACAGCGCTACACC
CAGGCCCTGATGCAGCAGCAGGCCGCCCTGGTAGCGGCTCACAGTGCCTACCTCAGCCCCATGGCCACC
ATGGCTGCCGTGCAGATGCAGCACATGGCTGCCATCAATGCCAATGGCCTCATCGCCACCCCATCACC
CCATCCTCAGGAACCAGCACCCCTCCTGCCATCGCTGCCAGCCTGTCTCTGCCATTCCGGCTGCCCTG
GGCGTCAACGGCTACAGCCCGGTGCCACCCAGCCCACTGGGCAGCCTGCCCTGATGCTGTATCCC
AACGGGGTTCACCCCTACCCAGCCAGAGCCCCGCGGCCCGCTGGACCCCTGCAGCAGGCTACGCG
GGGATGCAGCACTACACAGCAGCCTACCCAGCAGCCTACAGCCTGGTTGCACCTGCGTTCGCCGAGCCT
CCAGCCCTGGTCCGCCAGCAGCCCCACCACCCTCAACAGCAGCAGCAGCAGCAGCAGCAACAG
CAGCAGCAGCAAAGAGAAGGCCCTGATGGCTGCAACATCTTCATCTACCACCTGCCCCAGGAGTTCACT
GACTCAGAGATCCTCCAGATGTTTGTCCCTTTGGCCACGTCATCTCAGCCAAAGTCTTTGTTGACCGA
GCCACCAATCAAAGCAAATGTTTTGGCTTTGTGAGTTTCGACAATCCGGCCAGTCCCAGGCTGCCATC
CAGGCCATGAATGGCTTCCAGATCGGCATGAAGCGCCTCAAAGTCCAGCTAAAGCGGCCTAAGGATGCC
AACCGGCCCTACTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



[View online »](#)

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001172648

Insert Size: 1257 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001172648.3](#)

RefSeq Size: 5368 bp

RefSeq ORF: 1257 bp

Locus ID: 11189

Cytogenetics: 1q21.3

Protein Families: Transcription Factors

MW: 45 kDa

Gene Summary: Members of the CELF/BRUNOL protein family contain two N-terminal RNA recognition motif (RRM) domains, one C-terminal RRM domain, and a divergent segment of 160-230 aa between the second and third RRM domains. Members of this protein family regulate pre-mRNA alternative splicing and may also be involved in mRNA editing, and translation. Multiple alternatively spliced transcript variants encoding different isoforms have been identified in this gene. [provided by RefSeq, Feb 2010]
Transcript Variant: This variant (2) has an alternate 5' exon, as compared to variant 1. The resulting isoform (2) is shorter at the N-terminus, as compared to isoform 1.