

Product datasheet for **SC328693**

p53R2 (RRM2B) (NM_001172477) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	p53R2 (RRM2B) (NM_001172477) Human Untagged Clone
Tag:	Tag Free
Symbol:	RRM2B
Synonyms:	MTDPS8A; MTDPS8B; P53R2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC328693 representing NM_001172477. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTGTAGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTTGCTGTTGCGTCTTCCCCCTCACCGCAGTCACGCCAGCCGTTAGATTGCAAGTTGCAGGACCGC
TGTAGGAAATGTTATTCGCCGCGGTGAGACAGGCCTGTCCGCCGCCCTCGCCGAGCCTGGTTCGT
CGTTGCGAGCGCCGGGAGGCGCTCCCCGGGAGGGCGGAGGAAGGAGCTGACTTTGGGTTTGCCTCC
GCTCGTCTGCTGCCCCGGGGCCAGCCAAGGACGACGCTTGGAGGCTCAGGCCGGGAGATCATCTTCA
GACACCAACGAAAGTAAATAAAGTCAAATGAAGAGCCACTCCTAAGAAAGATTCTCGCCGTTTGTGTC
ATCTTTCCAATCCAGTACCCTGATATTTGAAAATGTATAAACAGGCACAGGCTTCTTCTGGACAGCA
GAAGAGGTCGACTTATCAAAGGATCTCCCTCACTGGAACAAGCTTAAAGCAGATGAGAAGTACTTATC
TCTCACATCTTAGCCTTTTTTGCAGCCAGTGATGGAATTGTAATGAAAATTTGGTGGAGCGCTTGTAGT
CAGGAGGTGCAGGTTCCAGAGGCTCGCTGTTTCTATGGCTTTCAAATTCTCATCGAGAATGTTCACTCA
GAGATGTACAGTTTGTGATAGACACTTACATCAGAGATCCCAAGAAAAGGGAATTTTTATTTAATGCA
ATTGAAACCATGCCCTATGTTAAGAAAAAGCAGATTGGGCTTGGCATGGATAGCAGATAGAAAATCT
ACTTTTGGGAAAGAGTGGTGGCCTTGTCTGTGTAGAAGGAGTTTTCTTCTCAGGATCTTTTGTGCTGCT
ATATTCGGCTAAAGAAGAGAGGTTTATGCCAGGACTCACTTTTCCAATGAACTCATCAGCAGAGAT
GAAGGACTTCACTGTGACTTTGCTTGCCTGATGTTCCAATACTTAGTAAATAAGCCTTCAGAAGAAAGG
GTCAGGGAGATCATTGTTGATGCTGTCAAATGAGCAGGAGTTTTTAACAGAAGCCTTGCAGTTGGC
CTCATTGGAATGAATTGCATTTTGTGAAACAGTACATTGAGTTTGTAGCTGACAGATTACTTGTGGAA
CTTGGATTCTCAAAGTTTTTTCAGGCAGAAAATCCTTTTGTATTTATGAAAACATTTCTTTAGAAGGA
AAAACAAATTTCTTTGAGAAACGAGTTTCAGAGTATCAGCGTTTTGCAGTTATGGCAGAAACCACAGAT
AACGTCTTACCTTGGATGCAGATTTTAA
ACGCGTACGCGGCCGCTCGAGCAGAAAATCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



[View online »](#)

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001172477

Insert Size: 1272 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001172477.1](#)

RefSeq Size: 4913 bp

RefSeq ORF: 1272 bp

Locus ID: 50484

UniProt ID:	<u>Q7LG56</u>
Cytogenetics:	8q22.3
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Glutathione metabolism, Metabolic pathways, p53 signaling pathway, Purine metabolism, Pyrimidine metabolism
MW:	48.8 kDa
Gene Summary:	<p>This gene encodes the small subunit of a p53-inducible ribonucleotide reductase. This heterotetrameric enzyme catalyzes the conversion of ribonucleoside diphosphates to deoxyribonucleoside diphosphates. The product of this reaction is necessary for DNA synthesis. Mutations in this gene have been associated with autosomal recessive mitochondrial DNA depletion syndrome, autosomal dominant progressive external ophthalmoplegia-5, and mitochondrial neurogastrointestinal encephalopathy. Alternatively spliced transcript variants have been described.[provided by RefSeq, Feb 2010]</p> <p>Transcript Variant: This variant (2) differs in the 5' UTR, lacks a portion of the 5' coding region, and initiates translation at an alternate start codon, compared to variant 1. The encoded isoform (2) has a distinct N-terminus and is longer than isoform 1.</p>