

Product datasheet for **SC328694**

VSIG1 (NM_001170553) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	VSIG1 (NM_001170553) Human Untagged Clone
Tag:	Tag Free
Symbol:	VSIG1
Synonyms:	1700062D20Rik; dj889N15.1; GPA34
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC328694 representing NM_001170553. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGTGTTCGCATTTTGGAAAGTCTTTCTGATCCTAAGCTGCCTTGCAGGTCAGGTTAGTGTGGTCAA
GTGACCATCCCAGACGGTTTCGTGAACGTGACTGTTGGATCTAATGTCACCTCTCATCTGCATCTACACC
ACCACTGTGGCCTCCCGAGAACAGCTTTCCATCCAGTGGTCTTTCTTCCATAAGAAGGAGATGGAGCCA
ATTTCTCACAGCTCGTGCCTCAGTACTGAGGGTATGGAGGAAAAGGCAGTCAGTCAGTGTCTAAAAATG
ACGCACGCAAGAGACGCTCGGGGAAGATGTAGCTGGACCTCTGAGATTTACTTTTCTCAAGGTGGACAA
GCTGTAGCCATCGGGCAATTTAAAGATCGAATTACAGGGTCCAACGATCCAGGTAATGCATCTATCACT
ATCTCGCATATGCAGCCAGCAGACAGTGAATTTACATCTGCGATGTTAAACAACCCCCAGACTTTCTC
GGCCAAAACCAAGGCATCCTCAACGTCAGTGTGTTAGTGAACCTTCTAAGCCCTTTGTAGCGTTCAA
GGAAGACCAGAACTGGCCACACTATTTCCCTTTCCTGTCTCTCTGCGCTTGAACACCTTCCCCTGTG
TACTACTGGCATAAACTTGAGGGAAGAGACATCGTGCCAGTGAAGAAAACCTTCAACCCAACCCCGGG
ATTTTGGTCATTGGAAATCTGACAAATTTGAACAAGTTATTACCAGTGTACTGCCATCAACAGACTT
GGCAATAGTTCCTGCGAAATCGATCTCACTTCTCACATCCAGAAGTTGGAATCATTGTTGGGCGCTTG
ATTGGTAGCCTGGTAGGTGCCCCATCATCATCTCTGTTGTGTGCTTCGCAAGGAATAAGGCAAAAGCA
AAGGCAAAAGAAATTTCTAAGACCATCGCGAACTTGAGCCAATGACAAAGATAAACCCAAGGGGA
GAAAGCGAAGCAATGCCAAGAGAAGACGCTACCCAAGTGAAGTAACTCTACCATCTTCCATTCATGAG
ACTGGCCCTGATACCATCCAAGAACCAGACTATGAGCCAAGCCTACTCAGGAGCCTGCCCCAGAGCCT
GCCCCAGGATCAGAGCCTATGGCAGTGCCTGACCTTGACATCGAGCTGGAGCTGGAGCCAGAAACGCAG
TCGGAATTGGAGCCAGAGCCAGAGCCAGAGCCAGAGTCAGAGCCTGGGTTGTAGTTGAGCCCTTAAGT
GAAGATGAAAAGGGAGTGGTTAAGGCA TAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



[View online »](#)

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001170553

Insert Size: 1272 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001170553.1](#)

RefSeq Size: 3236 bp

RefSeq ORF: 1272 bp

Locus ID: 340547

UniProt ID: [Q86XK7](#)

Cytogenetics: Xq22.3

Protein Families: Transmembrane

MW: 45.8 kDa

Gene Summary: This gene encodes a member of the junctional adhesion molecule (JAM) family. The encoded protein contains multiple glycosylation sites at the N-terminal region, and multiple phosphorylation sites and glutamic acid/proline (EP) repeats at the C-terminal region. The gene is expressed in normal stomach and testis, as well as in gastric, esophageal and ovarian cancers. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2009]

Transcript Variant: This variant (1) encodes the longer isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.