

Product datasheet for **SC328717**

Repulsive Guidance Molecule A (RGMA) (NM_001166286) Human Untagged Clone

Product data:

| | |
|---------------------------|--|
| Product Type: | Expression Plasmids |
| Product Name: | Repulsive Guidance Molecule A (RGMA) (NM_001166286) Human Untagged Clone |
| Tag: | Tag Free |
| Symbol: | RGMA |
| Synonyms: | RGM |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |
| Fully Sequenced ORF: | >SC328717 representing NM_001166286. Blue=Insert sequence Red=Cloning site Green=Tag(s) |


```
GCTCGTTTGTAGTAACCGTCAGAATTTTGTAAATACGACTACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGGTATGGGGAGAGGGGCAGGACGTTTCAGCCCTGGGATTCTGGCCGACCCCTCGCCTTCTCTCTGC
AGCTTCCCGCAGCCACCTCCCGTGCAAGATCCTCAAGTGAACCTCTGAGTTCTGGAGCGCCACGTCG
GGCAGCCACGCCCCAGCCTCAGACGACACCCCCGAGTTCTGTGCAGCCTTGGCAGCTACGCCCTGTGC
ACGCGGGGACGGCCCGCACCTGCCGGGTGACCTGGCCTACCACTCGGCCGTCATGGCATAGAGGAC
CTCATGAGCCAGCACAACTGCTCAAGGATGGCCCCACCTCGCAGCCACGCTGCGCAGCTCCCACCG
GCCGGAGACAGCCAGGAGCGCTCGGACAGCCCCGAGATCTGCCATTACGAGAAGAGCTTTCACAAGCAC
TCGGCCACCCCAACTACACGACTGTGGCCTTTCGGGGACCCACACCTCAGGACTTTCACCGACCGC
TTCCAGACCTGCAAGGTGCAGGGCGCCTGGCCGCTCATCGACAATAATTACCTGAACGTGCAGGTCACC
AACACGCCTGTGCTGCCCGCTCAGCGGCCACTGCCACCAGCAAGCTCACCATCATCTTCAAGAACTTC
CAGGAGTGTGGACCAGAAGGTGTACCAGGCTGAGATGGACGAGCTCCCGGCCGCTTCGTGGATGGC
TCTAAGAACGGTGGGACAAGCACGGGGCAACAGCCTGAAGTCACTGAGAAGGTGTGAGGCCAGCAC
GTGGAGATCCAGGCCAAGTACATCGGCACCACCATCGTGGTGCGCCAGGTGGGCCGCTACCTGACCTTT
GCCGTCCGATGCCAGAGGAAGTGGTCAATGCTGTGGAGGACTGGGACAGCCAGGCTCTACCTCTGC
CTGCGGGGCTGCCCCCTCAACCAGCAGATCGACTTCCAGGCCCTCCACACCAATGCTGAGGGCACCGGT
GCCCGCAGGCTGGCAGCCGCGCAGCCTGCACCCACAGCCCCGAGACCTTCCCATACGAGACAGCCGTG
GCCAAGTGAAGGAGAAGCTGCCGGTGGAGGACCTGTAATACCAGGCCTGCGTCTTCGACCTCCTCACC
ACGGGCGACGTGAACCTCACACTGGCCGCTACTACGCGTTGGAGGATGTCAAGATGCTCCACTCCAAC
AAAGACAACTGCACCTGTATGAGAGGACTCGGGACCTGCCAGGCAGGGCGGCTGCGGGGCTGCCCTG
GCCCCCGGCCCTCCTGGGCGCCTCGTCCCGCTCCTGGCCCTGCTCCCTGTGTTCTGCTAG
AGGCGACCGACGCGTACGCGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGAT
ATCCTGGATTACAAGGATGACGACGATAAGGTTTAA
```



[View online »](#)

Restriction Sites: SgfI-RsrII

Plasmid Map:



ACCN: NM_001166286

Insert Size: 1305 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001166286.1](#)

RefSeq Size: 3060 bp

RefSeq ORF: 1305 bp

Locus ID: 56963

UniProt ID: [Q96B86](#)

Cytogenetics: 15q26.1

MW: 47.5 kDa

Gene Summary: This gene encodes a member of the repulsive guidance molecule family. The encoded protein is a glycosylphosphatidylinositol-anchored glycoprotein that functions as an axon guidance protein in the developing and adult central nervous system. This protein may also function as a tumor suppressor in some cancers. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Oct 2009]

Transcript Variant: This variant (2) differs in the 5' UTR, lacks a portion of the 5' coding region, and initiates translation at a downstream start codon, compared to variant 1. The encoded isoform (2) has a shorter N-terminus, compared to isoform 1. Variants 2, 3, 5 and 6 encode the same isoform (2). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.