

Product datasheet for **SC328723**

UBR2 (NM_001184801) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	UBR2 (NM_001184801) Human Untagged Clone
Tag:	Tag Free
Symbol:	UBR2
Synonyms:	bA49A4.1; C6orf133; dj242G1.1; dj392M17.3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC328723 representing NM_001184801. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCGTCGGAGCTAGAGCCAGAGGTGCAGGCCATCGACCGGAGCTTGCTGGAATGTTCCGGCCGAGGAG
ATTGCGGGGAAATGGCTGCAAGCAACTGACCTCACTAGAGAAGTGTACCAGCATTAGCCCACTATGTA
CCCAAAATCTACTGCAGGGTCCCAACCCTTTCCACAGAAAGAAGACATGCTGGCACAGCATGTTTTG
TTGGGACCAATGGAATGGTACCTTTGTGGTGAAGATCCTGCATTTGGATTTCCAAAACCTTGACAAAGCA
AACAAACCTTCTCATCTTTGTGGTGGTGTAAAAAGTAGGAGAGCCTACATATCTTGCAGAGACTGT
GCAGTTGATCCAACCTTGTGTTTTGTGCATGGAGTGCTTTTGGGAAGTATTCACAGAGATCATCGATAT
AGGATGACAACATCAGGAGGTGGAGTTTCTGTGACTGTGGTGACTGAAGCCTGGAAAGAGGGTCCCT
TACTGTCAAAAACATGAACTTAACACCTCTGAAATGAGGAAGAAGAGGATCCTCTGTTCATTTATCA
GAAGATGTGATAGCAAGAACTTATAACATTTTTGCTATTACGTTTCGGTATGCAGTAGAAATATTAACC
TGGGAAAAAGAAAGTGAATTGCCAGCAGATTTAGAGATGGTAGAGAAGAGTGACACCTACTATTGCATG
CTGTTTAAATGATGAGGTTACACCTATGAACAAGTTATTTATACTCTTCAGAAAGCTGTTAACTGTACA
CAAAAAGAAGCTATTGGTTTTGCAACTACAGTAGATCGAGATGGGCGTAGGTCTGTTTCGATATGGAGAT
TTTCAGTATTGTGAGCAAGCAAAATCAGTAATTGTGAGAAATACCAGTAGACAGACAAAGCCACTCAA
GTTCAAGTTATGCATTCGTCTATTGTGCACATCAGAATTTTGGTTTAACTTTTGTCTTGGCTGGGA
AGTATTATTGGATATTCAGATGGCCTTCGCCGGATTTTATGTCAAGTTGGTTTACAAGAAGGCCAGAT
GGTGA AAAACTCTTCTCTAGTGGACAGACTGATGCTTAGTGATTCCAAATTATGGAAGGTGCTAGGAGT
GTATATCATCAGTTGTTTCATGAGCAGTCTGCTTATGGATTTGAAATACAAGAACTATTTGCTGTTTCGA
TTTGCAAAAAACTATGAGCGTTTGCAGAGTGATTATGTGACAGATGACCACGACAGAGTTTTTCAGTC
GCAGACCTCTCGTTTCAGATATTCACGGTTCCTTCACTTTTCTCTATCTGCTGGGCGCAGTGGCTCA
CCTCTGTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```



[View online >](#)

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001184801

Insert Size: 1320 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001184801.1](#)

RefSeq Size: 2493 bp

RefSeq ORF: 1320 bp

Locus ID: 23304

UniProt ID: [Q8IWW8](#)

Cytogenetics: 6p21.1

Protein Families: Druggable Genome

MW: 50 kDa

Gene Summary: This gene encodes an E3 ubiquitin ligase of the N-end rule proteolytic pathway that targets proteins with destabilizing N-terminal residues for polyubiquitylation and proteasome-mediated degradation. Alternative splicing results in multiple transcript variants.[provided by RefSeq, May 2010]
Transcript Variant: This variant (2) differs in the 3' coding region and 3' UTR, compared to variant 1. The resulting protein (isoform 2) is shorter and has a distinct C-terminus, compared to isoform 1.