

## Product datasheet for **SC328731**

### LARP4 (NM\_001170804) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	LARP4 (NM_001170804) Human Untagged Clone
Tag:	Tag Free
Symbol:	LARP4
Synonyms:	PP13296
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001170804, the custom clone sequence may differ by one or more nucleotides

```

ATGTTGCTTTTCGTGGAGCAGGTAGCATCTAAAGGAACTGGTTTAAATCCTAATGCCAAA
GTATGGCAAGAAAATTGCTCCTGGAAATACTGATGCCACCCAGTAACTCATGGAAGTAA
AGCTCTTGGCATGAAATAGCAGCTACATCAGGTGCTCATCCTGAGGGTAAATGCAGAGCTC
TCAGAAGATATATGTAAGAATATGAAGTAATGTATTCTTCATCTTGTGAAACCACAAGA
AATACTACAGGCATTGAAGAATCAACTGATGGGATGATTTTAGGACCAGAAGATCTGAGT
TACCAAATATATGATGTTTCCGGAGAAAGCAATTCAGCAGTTTCTACAGAAGACCTAAAA
GAATGTCTGAAGAAACAATTAGAATTCTGTTTTTCAGGAGAAAATTTGTCAAAGGATCTT
TACTTGATATCTCAAATGGATAGTGATCAGTTCATCCCAATTTGGACAGTTGCCAACATG
GAAGAAAATAAAAAAGTTGACTACAGACCCTGATCTAATTCTTGAAGTGTAAAGATCTTCT
CCCATGGTACAAGTTGATGAGAAGGGTGAGAAAAGTGAGACCAAGTCATAAGCGTTGTATT
GTAATCTTAGAGAGATTCTGAAACAACCAATAGAGGAAAGTAAAAGTTTGTTCAAA
AGTGAAAAGTGCCTCAAGTATAAGCTGTGAGTTGCACACAATAGCAACTGGTATATC
ACTTTCCAGTCAGACACAGATGCACAACAGGCTTTTAAATACTTAAAGAGAAGAAGTTAAA
ACATTTCCAGGGCAAGCCAATTATGGCAAGGATAAAAGCCATCAATACATTTTTTGCTAAG
AATGGTTATCGATTAATGGATTCTAGTATCTATAGTACCCCCATTCAAACCAAGCACAG
TATGCCTCCCAAGTCTTTATGCAGCCTGTATATAATCCTCACCAACAGTACTCGGTCTAT
AGTATTTGCCTCAGTCTTGGTCTCCAAATCCTACACCTTACTTTGAAACACCAGTGGCT
CCCTTTCCCAATGGTAGTTTTGTGAATGGCTTTAATTCGCCAGGATCTTATAAAACAAAT
GCTGCTGCTATGAATATGGGTCGACCATTCCAAAAAATCGTGTGAAGCCTCAGTTTAGG
TCATCTGGTGGTTCAGAACACTCAACAGAGGGCTCTGTATCCTTGGGGGATGGACAGTTG
AACAGATATAGTTCAAGAAACTTTCCAGCTGAACGGCATAACCCACAGTAACTGGGCAT
CAGGAGCAAACCTTACCTTCAGAAGGAGACTTCCACTTTGCAGGTGGAACAGAATGGGGAC
TATGGTAGGGGCAGGTAA

```

Restriction Sites: Please inquire



[View online »](#)

<b>ACCN:</b>	NM_001170804
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u>NM_001170804.1, NP_001164275.1</u>
<b>RefSeq Size:</b>	1681 bp
<b>RefSeq ORF:</b>	1338 bp
<b>Locus ID:</b>	113251
<b>Cytogenetics:</b>	12q13.12
<b>Protein Families:</b>	Transcription Factors
<b>Gene Summary:</b>	<p>RNA binding protein that binds to the poly-A tract of mRNA molecules (PubMed:21098120). Associates with the 40S ribosomal subunit and with polysomes (PubMed:21098120). Plays a role in the regulation of mRNA translation (PubMed:21098120). Plays a role in the regulation of cell morphology and cytoskeletal organization (PubMed:21834987, PubMed:27615744). [UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (5) differs in the 3' UTR and lacks a portion of the 3' coding region, compared to variant 1. The encoded isoform (e) has a shorter C-terminus, compared to isoform a.</p>