

## Product datasheet for **SC328744**

### CP2c (TFCP2) (NM\_001173453) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	CP2c (TFCP2) (NM_001173453) Human Untagged Clone
Tag:	Tag Free
Symbol:	TFCP2
Synonyms:	LBP1C; LSF; LSF1D; SEF; TFCP2C
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC328744 representing NM_001173453. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCCTGGGCTCTGAAGCTGCCTCTGGCCGACGAAGTATTGAATCCGGGTTGGTGCAGGACTTTGAT
GCTAGCCTGTCCGGGATCGGCCAGGAAGTGGGTGCTGGTGCCTATAGCATGAGTATGTCCTTGCATTG
CCCATTTTAAAGCAAGAAGAGTCGAGTTTGCTCCTGATAATGAGAATAAAATCTGCCTTTTCAATAT
GTGCTTTGTGCTGCTACCTCTCCAGCAGTGAACTCCATGATGAAACCCTAACGTATCTCAATCAAGGA
CAGTCTTATGAAATTCGAATGCTAGACAATAGGAACTTGAGAACTCCAGAAATTAATGGCAAATTG
GTGAAGAGTATATCCGTGTGGTGTCCATGACAGAAGGCTTCAGTACACTGAGCATCAGCAGCTAGAG
GGCTGGAGGTGGAACCGACCTGGAGACAGAATTCTTGACATAGATATCCCGATGTCTGTGGGTATAATC
GATCCTAGGGCTAATCCAACCTCACTAAATACAGTGGAGTTCTGTGGGACCCTGCAAAGAGGACATCT
GTGTTTATTCAGCCCAAAGGTGCAGACAGAAAGCAAAAAACGGATAGGGAAAAATGGAGAAACGAACA
CCTCATGAAAAGGAGAAATATCAGCCTTCTATGAGACAACCACTACTCACAGAGTGTCTCCATGGCCC
GAGATCACGTATGTCAATAACTCCCCATCACCTGGCTTCAACAGTTCCCATAGCAGTTTTTCTCTGGG
GAAGGAAATGGTTCACCAAACCACAGCCAGAGCCACCCCTCCAGTCACAGATAACCTCTTACCAACA
ACCACACCTCAGGAAGCTCAGCAGTGGTTGCATCGAAATCGTTTTTCTACATTACAAGGCTTTTCACA
AACTTCTCAGGGCAGATTTATTGAAATTAAGTATGATGATGATGATGATGATGATGATGATGATGATGAT
GGAATCAGACTTTTTAATGCATTAAGGCGCGGATGGTGCCTCAAGGTTAACCATTTATGTTTGTGAG
GAATCACTGCAGTTGAGGAGCAGCAACAACAGCAGCAGCAACAGCAGCAGAAAGCATGAGGATGGAGAC
TCAAATGGTACTTTCTCGTTTACCATGCTATCTATCTAGAAGAACTAACAGCTGTTGAATTGACAGAA
AAAATTGCTCAGCTTTTCCAGCATTTCCCTTGGCAGATCAGCCAGATTTACAAGCAGGGGCCAACAGGA
ATTCATGTGCTCATCAGTATGAGATGATACAGAACTTTAGGAAGAAGCATGTTTTATTCTGGACACA
ATGAAAGAAACCAATGATAGCTATCATATCATACTGAAGTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```



[View online »](#)

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM\_001173453

Insert Size: 1353 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_001173453.1](#)

RefSeq Size: 3559 bp

RefSeq ORF: 1353 bp

Locus ID: 7024

UniProt ID: [Q12800](#)

**Cytogenetics:** 12q13.12-q13.13

**Protein Families:** Transcription Factors

**MW:** 51.3 kDa

**Gene Summary:** This gene encodes a transcription factor that binds the alpha-globin promoter and activates transcription of the alpha-globin gene. The encoded protein regulates erythroid gene expression, plays a role in the transcriptional switch of globin gene promoters, and it activates many other cellular and viral gene promoters. The gene product interacts with certain inflammatory response factors, and polymorphisms of this gene may be involved in the pathogenesis of Alzheimer's disease. [provided by RefSeq, Mar 2010]  
Transcript Variant: This variant (3) lacks an alternate in-frame exon in the central coding region, and uses an alternate in-frame splice site in the 3' coding region, compared to variant 1. This results in a shorter protein (isoform 3), compared to isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.