

Product datasheet for **SC328745**

SSH3BP1 (ABI1) (NM_001178120) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SSH3BP1 (ABI1) (NM_001178120) Human Untagged Clone
Tag:	Tag Free
Symbol:	ABI1
Synonyms:	ABI-1; ABLBP4; E3B1; NAP1BP; SSH3BP; SSH3BP1
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001178120, the custom clone sequence may differ by one or more nucleotides

```

ATGGCAGAGCTGCAGATGTTACTAGAGGAGGAGATCCCGTCTGGCAAGAGGGCGCTGATC
GAGAGTTACCAGAACCTGACTCGGGTGGCAGACTACTGTGAAAACAACTACATACAGGCT
ACAGACAAGAGAAAAGCTTTAGAGGAGACCAAAGCCTATACAACCAATCTCTAGCTAGT
GTTGCTTATCAAATAAATGCATTGGCCAACAATGTACTCCAGTTGCTGGATATCCAAGCC
TCTCAGCTTCGGAGAATGGAGTCTTCCATCAATCATATCTCACAGACTGTGGATATTCAT
AAGGAGAAAAGTGGCACGAAGAGAGATTGGTATTTTACAACAATAAGAATACATCAAGA
ACTCACAATAATAGCACCTGCGAATATGGAGCGCCCTGTAAGGTATATTCGGAAACCT
ATCGATTACACAGTTCTGGATGATGTGGCCATGGTGTCAAGTGGCTAAAAGCCAAGCAT
GGAAATAACCAGCCTGCAAGAAGTGGCACACTGTCGAGAACAATCCTCCTACTCAGAAA
CCGCCAAGTCTCCCATGTCAGGCCGGGAACACTGGGACGGAATACTCCTTATAAAACC
CTGGAACCTGTTAAACCCCAACAGTTCCTAATGACTATATGACCAAGTCTGCTAGGCTT
GGAAGTCAGCATAGTCCAGGCAGGACAGCATCTTAAATCAGAGACCAAGGACACACAGT
GGAAGTAGTGGAGGAAGTGGAAAGTCGAGAAAACAGTGGTAGCAGTAGTATTGGCATTCCC
ATTGCTGTGCCTACACCTTCGCCACCCACTATTGGACCAGCCCCGGGCTCAGTCTCTGGT
TCCCAGTATGGCACAATGACCAGGCAGATATCTCGACACAACCTACTACTTCTTCGACA
TCTTCTGGTGGATACAGACGAACCTCCCTCTGTGACTGCTCAATTTTCTGCTCAGCCTCAT
GTTAATGGAGGTCCACTTTATTCTCAAAATTCATTGCTGATAGTCCAACCTCCACCGCCA
CCACCTCCACGATGACATTCATGTTGATGACTCTCCACCTCCCCACCACCA
CCAGTGGATTATGAAGATGAGGAGGCTGCAGTAGTTCAGTATAATGATCCATATGCAGAT
GGGGATCTGCTTGGGCCCAAGAATTATATTGAGAAAGTTGTTGCAATATATGATTAT
ACAAAAGACAAGGATGATGAGCTGTCAATTTATGGAGGGTGAATCATTATGTTATAAAG
AAGAAATGATGATGGCTGGTATGAAGGAGTCTGCAATCGAGTACTGGTCTGTTCCCTGG
AACTATGTTGAATCAATCATGCACTATACTGATTA

```

Restriction Sites: Please inquire



[View online »](#)

ACCN:	NM_001178120
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001178120.1, NP_001171591.1</u>
RefSeq Size:	3554 bp
RefSeq ORF:	1356 bp
Locus ID:	10006
UniProt ID:	<u>Q8IZP0</u>
Cytogenetics:	10p12.1

Gene Summary:

This gene encodes a member of the Abelson-interactor family of adaptor proteins. These proteins facilitate signal transduction as components of several multiprotein complexes, and regulate actin polymerization and cytoskeletal remodeling through interactions with Abelson tyrosine kinases. The encoded protein plays a role in macropinocytosis as a component of the WAVE2 complex, and also forms a complex with EPS8 and SOS1 that mediates signal transduction from Ras to Rac. This gene may play a role in the progression of several malignancies including melanoma, colon cancer and breast cancer, and a t(10;11) chromosomal translocation involving this gene and the MLL gene has been associated with acute myeloid leukemia. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene, and a pseudogene of this gene is located on the long arm of chromosome 14. [provided by RefSeq, Sep 2011]

Transcript Variant: This variant (7) lacks two alternate in-frame coding exons and uses a different splice site, compared to variant 1. The resulting protein (isoform g) has the same N- and C-termini but is shorter when it is compared to isoform a. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.